




## 170K1=80 FreeWalk, 170K1=120 FreeWalk

<b>DE</b> Gebrauchsanweisung (Benutzer) .....	3
<b>EN</b> Instructions for use (user) .....	8
<b>FR</b> Instructions d'utilisation (Utilisateur) .....	13
<b>IT</b> Istruzioni per l'uso (Utilizzatore) .....	18
<b>ES</b> Instrucciones de uso (Usuario) .....	23
<b>PT</b> Manual de utilização (Usuário) .....	29
<b>NL</b> Gebruiksaanwijzing (Gebruiker) .....	34
<b>SV</b> Bruksanvisning (Användare) .....	39
<b>DA</b> Brugsanvisning (Bruger) .....	44
<b>NO</b> Bruksanvisning (Bruker) .....	49
<b>FI</b> Käyttöohje (Käyttäjä) .....	54
<b>HU</b> Használati utasítás (felhasználó) .....	59
<b>TR</b> Kullanma talimatı (Kullanıcı) .....	64

## INFORMATION

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### INFORMATION

Date of last update: 2022-10-27

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the use of the 170K1=80 and 170K1=120 Free Walk orthosis system.

## 2 Product description

### 2.1 Design

The Ottobock FreeWalk orthoses include an orthotic knee joint with lock. Three modes can be selected on the orthotic knee joint:

- **FreeWalk mode:** The orthotic knee joint locks in the stance phase and unlocks in the swing phase.
- **Permanent lock:** The orthotic knee joint is permanently locked.
- **Free motion mode:** The orthotic knee joint is permanently unlocked.

## 3 Intended use

### 3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

Reuse of the product on other persons is not permitted.

### 3.2 Indications

- Partial or total paralysis of the leg muscles
- Orthopaedic diseases of the lower limbs

Indications must be determined by the physician.

### 3.3 Contraindications

#### 3.3.1 Absolute Contraindications

- Ankle movement < 10°
- Knee flexion contractures > 10°

- Unstable knee joint with compensating movement in varus position
- Valgus deformity > 15°

### 3.3.2 Relative Contraindications

- Valgus deformity 0° to 15° (fitting only with 170D30 medial support)
- Spasticity resulting in involuntary movements in the ankle joint

### 3.4 Restrictions for use

The product is not suitable for use in sports.

### 3.5 Lifetime

The product is designed for a lifetime of **3 years** when used as intended and assembled professionally.

### 3.6 Qualification

Patients may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

## 4 Safety

### 4.1 Explanation of warning symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 4.2 General safety instructions

#### **WARNING**

##### **Operation of motor vehicles**

Risk of accident due to limited body function

- ▶ The product has not been approved by the manufacturer for the operation of motor vehicles.
- ▶ Observe the applicable legal regulations for the operation of motor vehicles and have your driving ability examined by an authorised test centre.

#### **CAUTION**

##### **Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

#### **CAUTION**

##### **Excessive strain on load-bearing components**

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

#### **CAUTION**

##### **Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.

- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

**⚠ CAUTION**

**Incorrect alignment, assembly or adjustment**

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

**⚠ CAUTION**

**Unallowable use of lubricants**

Risk of injury and loss of functionality as well as damage to the product

- ▶ The product is designed to be free of lubricants. Avoid contact with lubricants (e.g. grease, oil).
- ▶ After contact with lubricants, promptly clean the product or parts of the product with a degreasing cleaning agent (e.g. acetone or isopropyl alcohol).

**⚠ CAUTION**

**Use on more than one patient**

Skin irritation, formation of eczema or infections

- ▶ Only use the orthosis on one patient.

**⚠ CAUTION**

**Improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

**⚠ CAUTION**

**Improper handling**

Risk of falling due to accidentally opening the lock

- ▶ Prior to handing over the orthosis, ensure that the locks are engaged securely and in parallel.
- ▶ Avoid accidentally opening the locks.
- ▶ Instruct the patient in the handling of the orthotic joint.

**⚠ CAUTION**

**Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

**⚠ CAUTION**

**Improper moulding or application**

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.

- ▶ Always apply the orthosis according to the information in the instructions.

### **⚠ CAUTION**

#### **Using orthosis without padding**

Risk of local pressure points, superficial injuries or skin irritation due to skin contact with rigid or sharp orthosis components

- ▶ Only use the orthosis with undamaged padding.

### **⚠ CAUTION**

#### **Improper handling**

Risk of falling due to accidentally opening or engaging the lock

- ▶ Avoid accidentally opening or engaging the lock.
- ▶ Instruct the patient in the correct handling of the orthosis joints.

### **NOTICE**

#### **Contact with saltwater or water containing chlorine/soap, fresh water or abrasive substances (e.g. sand)**

Damage and premature product wear

- ▶ Following any contact with the substances identified above, promptly clean the product in accordance with the section "Cleaning".

### **NOTICE**

#### **Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

### **NOTICE**

#### **Exposure of the product to unsuitable environmental conditions**

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below  $-10\text{ °C}$  ( $14\text{ °F}$ ) or above  $+60\text{ °C}$  ( $140\text{ °F}$ ) (e.g. sauna, excessive sunlight, drying on a radiator).

### **NOTICE**

#### **Improper settings**

Loss of functionality due to damaged steel cable

- ▶ Do not use the set screw of the steel cable on the ankle joint bottom part to adjust the lock.

### **NOTICE**

#### **Thermal overloading of product components**

Loss of function due to improper thermal treatment

- ▶ Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

## 5 Use

### 5.1 Application

- > To apply the orthosis, the user must be sitting down.
- > The orthosis is in FreeWalk mode (see page 12).
- 1) Slide the foot component of the orthosis into the shoe (see fig. 1).
- 2) Put on the shoe (see fig. 2).
- 3) Slide the orthosis over the leg from the front.
- 4) If required, adjust the width of the closures. Adjust the hook-and-loop closures to do so.
- 5) Fasten the closures (see fig. 3).

### 5.2 Removal

- > To remove the orthosis, the user must be sitting down.
- > The orthosis is in Free Walk mode (see page 12).
- 1) Open the quick-release closures on the orthosis and pull it forward off the leg.
- 2) Remove the shoe.
- 3) Pull the foot component of the orthosis out of the shoe.

### 5.3 Selecting the function modes

#### Selecting FreeWalk mode

- ▶ Lock the orthosis in extension and slide the 3-phase switch to the upper position (see fig. 4).
  - The orthotic knee joint unlocks in the swing phase.

#### Selecting permanent locking

- ▶ Lock the orthosis in extension and slide the 3-phase switch to the lower position (see fig. 4).
  - The orthotic knee joint is permanently locked.

#### Selecting free motion mode

- > **Prerequisite:** The 3-phase switch is in the upper position.
- 1) Unlock and flex the orthotic knee joint.
- 2) Slide the 3-phase switch to the lower position (see fig. 5).
  - The orthotic knee joint can be flexed freely at all times.

## 6 Cleaning

Clean the orthosis regularly:

- 1) Remove the straps from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Hand-wash the straps and pads in warm water at **30 °C (86 °F)** using neutral detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

## 7 Maintenance

### INFORMATION

**The product may be exposed to increased loads by the patient.**

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

## 8 Disposal

Dispose of the product in accordance with national regulations.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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