
50S3 Malleo Sprint

DE	Gebrauchsanweisung	4
EN	Instructions for use	8
FR	Instructions d'utilisation	12
IT	Istruzioni per l'uso	17
ES	Instrucciones de uso	21
PT	Manual de utilização	25
NL	Gebruiksaanwijzing	30
SV	Bruksanvisning	34
DA	Brugsanvisning	38
NO	Bruksanvisning	42
FI	Käyttöohje	46
PL	Instrukcja użytkowania	50
HU	Használati utasítás	54
CS	Návod k použití	59
HR	Upute za uporabu	63
SK	Návod na používanie	67
BG	Инструкция за употреба	71
TR	Kullanma talimatı	75
EL	Οδηγίες χρήσης	79
ET	Kasutusjuhhis	84
RU	Руководство по применению	88
ZH	使用说明书	93

	Kennzeichen Reference number	Farbe Colour	Größe Size	Fersenumfang Heel circumference	
				cm	inch
	50S3=XS-1	grau grey	XS	25,5– 28	10–11
	50S3=S-1	grau grey	S	28– 30,5	11–12
	50S3=M-1	grau grey	M	30,5– 33	12–13
	50S3=L-1	grau grey	L	33– 35,5	13–14
	50S3=XL-1	grau grey	XL	35,5– 38	14–15

Material	TPU, SBR, PES, PA, PU, EA
-----------------	---------------------------

INFORMATION

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INFORMATION

Date of last update: 2022-11-15

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50S3 Malleo Sensa ankle orthosis.

2 Intended use

2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

2.2 Indications

- Functional treatment of ligament instability
- Post-operative treatment of ligament ruptures
- Ankle sprains
- Functional, conservative therapy of lateral ligament ruptures

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

Not known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; pronounced varicose veins, especially with return flow impairment, lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support will be applied; sensory and circulatory disorders in the foot area, e.g. in case of diabetic neuropathy.

2.4 Lifetime

The product is designed for a maximum lifetime of **1 year**.

2.5 Effects

This orthosis stabilises the ankle with lateral elements and straps. Eversion and inversion of the ankle joint are restricted without negatively affecting plantar and dorsal flexion.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General safety instructions

 **CAUTION**

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

 **CAUTION**

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Size Selection

- ▶ Select the orthosis size based on the circumference of the ankle (see size chart).

4.2 Fitting and Application

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.

- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

NOTICE

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

The orthosis is delivered ready for trial fitting. Adjustment is not required.

> The strap on the orthosis is open and the fasteners are loosened.

1) Pull the orthosis over your foot (see fig. 1) and tie it like a normal shoe (see fig. 2).

Optional: Use the quick fasteners (see fig. 3).

2) Pass the outer end of the strap over the instep and pull it under the foot from the inner (medial) side.

3) Fasten the outer end of the strap over the ankle on the outer (lateral) side (see fig. 4).

4) Pass the inner end of the strap over the instep and pull it under the foot from the outer (lateral) side.

5) Fasten the inner end of the strap over the ankle on the inner (medial) side (see fig. 5).

6) If the physician has allowed weight-bearing on the foot:
Loosen the ends of the strap and re-tighten (see fig. 6).

4.3 Cleaning

1) Close all hook-and-loop closures.

2) Hand wash the orthosis in warm water at 40°C with mild laundry detergent.

3) Rinse well.

4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

6.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

6.3 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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