
50R50 Lumbo Direxa 50R51 Lumbo Direxa Women 50R52 Lumbo Direxa High

DE	Gebrauchsanweisung	3
EN	Instructions for use	8
FR	Instructions d'utilisation	13
IT	Istruzioni per l'uso	19
ES	Instrucciones de uso	24
PT	Manual de utilização	30
NL	Gebruiksaanwijzing	35
SV	Bruksanvisning	41
DA	Brugsanvisning	46
NO	Bruksanvisning	51
FI	Käyttöohje	57
IS	Notkunarleiðbeiningar	62
PL	Instrukcja użytkowania	67
HU	Használati utasítás	73
CS	Návod k použití	78
SK	Návod na používanie	83
TR	Kullanma talimatı	88
ET	Kasutusjuhhis	93
RU	Руководство по применению	98
JA	取扱説明書	104
ZH	使用说明书	109

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INFORMATION

Date of last update: 2023-01-11

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50R50 Lumbo Direxa, 50R51 Lumbo Direxa Women and 50R52 Lumbo Direxa High lumbar braces.

2 Intended use

2.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the back and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

2.2 Indications

50R50 Lumbo Direxa, 50R51 Lumbo Direxa Women

- Moderate low back pain and sciatica
- Myostatic insufficiency/muscular imbalance of the lumbar trunk muscles
- Moderate degenerative changes of the lumbar spine (e.g., osteochondrosis, facet joint syndrome)
- Stabilization after discectomy L4-S1
- Intervertebral disc protrusion L4-S1
- Moderate to severe irritation of the sacroiliac joint

The indication must be determined by the physician.

50R52 Lumbo Direxa High

- Slight to moderate low back pain, thoracic spine/lumbar spine syndrome, sciatica, thoracic and lumbar ligament sprain
- Myostatic insufficiency/muscular imbalance of the lumbar and lower thoracic trunk muscles
- Minor to moderate degenerative changes of the lumbar spine and lower thoracic spine (e.g., osteochondrosis, facet joint syndrome)
- Intervertebral disc protrusion L1-S1
- Stabilization after discectomy L1-S1
- Moderate irritation of the sacroiliac joint

The indication must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

The following conditions require consultation with a physician: skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the device is fitted, tumours, lymphatic congestion – including unclear soft tissue swelling distal to the device, and dysaesthesia in the trunk.

2.4 Lifetime

The product is designed for a maximum lifetime of **1 year**.

2.5 Mechanism of Action

The brace promotes sensorimotor function and activates muscular stabilization of the lumbar spine. It helps alleviate pain and relieves the lumbar spine.

3 Safety

3.1 Explanation of warning symbols

 **WARNING** Warning regarding possible serious risks of accident or injury.

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

CAUTION

Brace is put on too tightly

Risk of local pressure points and restriction of adjacent blood vessels or nerves

- ▶ Check that the brace is correctly positioned and fitted.

NOTICE**Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the Size

- 1) Measure the waist circumference approx. **2 cm** above the iliac crest.
- 2) Determine the brace size (see size chart).

4.2 Fitting and Application

⚠ CAUTION**Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

⚠ CAUTION**Direct skin contact with the product**

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.

NOTICE**Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

Optional – inserting the pad:

A pad is also available in addition to the brace and can be ordered as a service part if necessary.

- ▶ **Optional – inserting the pad:** Use the hook and loop closure to attach the pad to the inside of the brace so that it is at least 1 cm above the seat surface when seated.

Fitting and application

- 1) **Optional – inserting the pad:** Use the hook and loop closure to attach the pad to the inside of the brace. In the sitting position, the pad has to be at least 1 cm above the seat surface.
- 2) **Applying** (see fig. 1): Place the brace around the trunk. Place the fingers in the hand straps of the abdominal closure. Apply equal force to both sides to pull the brace forward.
- 3) **Hook and loop material** (see fig. 2): Attach the right part to the left part of the abdominal closure. Remove the fingers from the left and then from the right hand strap. Make sure that the hook and loop surfaces are closed evenly against each other and not in contact with clothing or skin.
- 4) **Increasing compression:** Pull the ends of the straps forward (see fig. 3) and attach them to the closure (see fig. 4).

Final steps

- 1) **Final inspection:** Verify the correct fit and proper positioning of the brace during handover to the patient. The brace is positioned correctly when the lower edge of the brace extends just below the groin.
- 2) **Instruct the patient:** Explain to the patient how to apply the brace and ensure that it is positioned correctly. Practice the application and removal of the brace with the patient.

4.3 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the brace regularly:

- 1) Fasten all hook-and-loop closures.

- 2) Recommendation: Use a laundry bag or net.
- 3) Wash the support in warm water at **40 °C** (104 °F) with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

6.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

6.2 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

6.3 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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