

## 17LK3

<b>DE</b> Gebrauchsanweisung .....	5
<b>EN</b> Instructions for use .....	14
<b>FR</b> Instructions d'utilisation .....	24
<b>IT</b> Istruzioni per l'uso .....	35
<b>ES</b> Instrucciones de uso .....	45
<b>PT</b> Manual de utilização .....	56
<b>NL</b> Gebruiksaanwijzing .....	66
<b>SV</b> Bruksanvisning .....	76
<b>DA</b> Brugsanvisning .....	85
<b>NO</b> Bruksanvisning .....	95
<b>FI</b> Käyttöohje .....	104
<b>HU</b> Használati utasítás .....	114
<b>CS</b> Návod k použití .....	124
<b>SK</b> Návod na používanie .....	134
<b>RU</b> Руководство по применению .....	143
<b>ZH</b> 使用说明书 .....	154


1






2

3

## INFORMATION

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### INFORMATION

Date of last update: 2023-02-17

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.

- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the processing of the 17LK3 unilateral knee joint.

## 2 Product description

### 2.1 Available sizes

The orthosis joint is available in five sizes and can be used unilaterally or bilaterally. The choice of size is based on the patient's body weight and the conditions of use.

#### INFORMATION

- ▶ When using the product unilaterally, and in case of flexion contracture in the knee or hip  $> 10^\circ$  or distinct torsion or valgus/varus instabilities or valgus/varus malpositions  $> 10^\circ$  or increased bodily activity, the next higher size must be used.
- ▶ For tubercle seat the product must be fitted bilaterally.

	Maximum patient weight [kg]				
	17LK3=10	17LK3=12	17LK3=14	17LK3=16	17LK3=20
Unilateral use (lateral or medial)	15	20	50	85	110
Bilateral use (lateral and medial)	25	40	80	120	160

## 2.2 Components

Scope of delivery (see fig. 1)				
Item	Pcs.	Description	Joint size 17LK3=L* , 17LK3=R* *=10,=12,=14,= 16,=20	Reference number
1	1	Upper joint section		30X94=*
2	1	Joint screw		501F9=M6x8
3	1	Axle	*=10,=12,=14 *=16 *=20	4A101=5x10.6 4A101=6x12.7 4A101=6x13.7
4	1	Lower joint section		30U132=*
5	2	Bushing	*=16 *=20	30Y87=2 30Y87=1
Not illustrated	2	Bushing	*=10,=12,=14	30Y87=3
6	2	Bearing washer	*=16,=20	170Z93=1
8	4	Screw	*=10,=12 *=14 *=16,=20	501S137=M4x8 501S137=M4x1-0-2 501S137=M5x1-0-1
9	1	Protective plug	*=10,=12,=14 *=16,=20	30Y91=3 30Y91=1
10	1	Lock washer	*=10,=12,=14 *=16,=20	507S96=3.2 507S96=4
11	1	Bearing nut	*=10,=12,=14 *=16 *=20	17Y93=9x9.2xM6 17Y93=9x10.75x-M6 17Y93=9x11.75x-M6

<b>Scope of delivery (see fig. 2)</b>				
<b>Item</b>	<b>Pcs.</b>	<b>Description</b>	<b>Joint size 17LK3=L* , 17LK3=R* *=10,=12,=14,=- 16,=20</b>	<b>Reference number</b>
1	1	Lamination dummy		30Y268=10,=12,- =14,=16,=20

<b>Scope of delivery (see fig. 3)</b>				
<b>Item</b>	<b>Pcs.</b>	<b>Description</b>	<b>Joint size 17LK3=L* , 17LK3=R* *=10,=12,=14,=- 16,=20</b>	<b>Reference number</b>
1	1	Knee lever		30Y265=*-A
2	2	Washer		30Y267=3.3
3	1	Axle	17LK3=10,=12,- =14 17LK3=16 17LK3=20	4A101=3.2x15.1 4A101=3.2x17.1 4A101=3.2x18.1
4	1	Threaded sleeve		21A12
5	1	Axle	17LK3=10,=12,- =14 17LK3=16 17LK3=20	4A101=4x15.1 4A101=4x17.1 4A101=4x18.1
6	2	Washer		30Y267=4.2
7	1	Lock wedge	17LK3=10,=12,- =14 17LK3=16 17LK3=20	17Y156=* 17Y157=* 17Y158=*
8	2	Lock washer		507S96=2.3
9	1	Compression spring	17LK3=10,=12 17LK3=14,=16,- =20	513D83=0.75x4.- 3x19 513D83=0.75x4.- 6x20
10	1	Perlon cable		21A18=2

17Y162* optional temporary switch scope of delivery (see fig. 4)		
Joint size	Qty.	Article number
17LK3=L10/17LK3=R10	1	17Y162=3
17LK3=L10-T/17LK3=R10-T		
17LK3=L12/17LK3=R12		
17LK3=L12-T/17LK3=R12-T		
17LK3=L14/17LK3=R14		
17LK3=L14-T/17LK3=R14-T		
17LK3=L16/17LK3=R16	1	17Y162=2
17LK3=L16-T/17LK3=R16-T		
17LK3=L20/17LK3=R20	1	17Y162=1
17LK3=L20-T/17LK3=R20-T		

Not included in scope of delivery			
Item	Pcs.	Description	Reference number
(see fig. 1), item 7	2	Lamination bar	17LS3=*
Not illustrated	1	Side bar	17LV3=*
Not illustrated	1	Side bar with contoured calf	17LV3=L-*, =R-*
(see fig. 2), item 2	1	Shoulder screw (part of 743R6 adjustment set)	30Y89

### 3 Intended use

#### 3.1 Indications for use

The manufacturer recommends that the orthosis joint be processed using lamination technology or carbon prepreg technology if used on one side of the body.

The orthosis joint is intended **exclusively** for orthotic fittings of the lower limbs on one patient.

#### 3.2 Indications

Partial or total paralysis of the leg muscles.

Indications must be determined by the physician.

### 3.3 Lifetime

The product is designed for the following lifetime, given intended use and proper installation:

Joint sizes 10, 12: **1 year**

Joint sizes 14, 16, 20: **3 years**

### 3.4 Qualification

Patients may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

## 4 Safety

### 4.1 Explanation of warning symbols

 **WARNING** Warning regarding possible serious risks of accident or injury.

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 4.2 General safety instructions

#### **WARNING**

##### **Operation of motor vehicles**

Risk of accidents due to restricted body function

- ▶ Observe the applicable legal and insurance regulations for the operation of motor vehicles and have your driving ability examined and certified by an authorised agency.

#### **CAUTION**

##### **Excessive strain on load-bearing components**

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

#### **CAUTION**

##### **Incorrect alignment, assembly or adjustment**

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

**⚠ CAUTION**

**Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

**⚠ CAUTION**

**Reaching into the area of the joint mechanism**

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

**⚠ CAUTION**

**Improper handling**

Risk of falling due to accidentally opening or engaging the lock

- ▶ Before placing weight on the orthosis, verify that the lock is in the desired position.
- ▶ Avoid accidentally opening or engaging the lock.
- ▶ Instruct the patient in the correct handling of the orthotic joints.

**⚠ CAUTION**

**Excessive strain due to use on more than one patient**

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

**NOTICE**

**Thermal overloading of the product**

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.

- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

#### NOTICE

##### **Thermal overloading of product components**

Loss of function due to improper thermal treatment

- ▶ Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

#### NOTICE

##### **Exposure of the product to unsuitable environmental conditions**

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

## 5 Preparing the product for use

#### INFORMATION

The instructions for fabricating an orthosis with unilateral joints are included in the 646T5441 technical information. The laminate structure described here is a recommendation from Ottobock and must be adjusted to the needs of the patient.

### **Laminating the orthosis**

- > **Required tools and materials:** 616G15 woven carbon stockinette, 17Y106 PVC profile material, 616G12 carbon fibre cloth, 99B81=\* PVA bag, 623T3=\* Perlon stockinette, 636K8 Plastaband, 633W8 wax, 17LS3=\* lamination bars, 617H55 C-Orthocryl lamination resin, 743R6 adjustment set for orthotic joints.
- > **Optional:** 30Y268=\* lamination dummy
  - 1) Place cut pads or separating layers on plaster model.
  - 2) Pull a layer of Perlon stockinette over the plaster model.
  - 3) Align the bars on the plaster model.
- 4) **If using a lamination dummy:** Fix the lamination bars to the 743Y56=4 alignment insert using the lamination dummy and the shoulder screw.

- 5) **If not using a lamination dummy:** Fix the lamination bars to the 743Y56=4 alignment insert using the base body (not illustrated).
- 6) Pull a PVA bag over the plaster model to isolate it.
- 7) Pull a layer of Perlon stockinette over the plaster model.
- 8) Secure the reinforced bars on the plaster model. Carry out the reinforcement in accordance with the patient's needs.
- 9) Seal the lamination dummy or base body of the joint with Plastaband. Fill in hollow areas.
- 10) Pull another layer of Perlon stockinette over the plaster model.
- 11) Affix further reinforcement layers for the support surfaces as described under point 8.
- 12) Complete reinforcement with a layer of Perlon stockinette and designer fabric.
- 13) Pull on a PVA bag and apply vacuum.
- 14) Then laminate with C-Orthocryl lamination resin.

### Assembly

	Torque [Nm]				
	17LK3=10	17LK3=12	17LK3=14	17LK3=16	17LK3=20
Joint screw (item 15)	8	8	8	8	8
Screw (item 22)	3.5	3.5	3.5	7	7

- 1) Demould the model and cut free the lamination dummy or the orthotic joint.
- 2) **17LK3=10,=12,=14:** Insert the lower joint section with the bushings into the upper joint section.
- 3) **17LK3=16,=20:** Insert the lower joint section with the bushings and washers into the upper joint section.
- 4) Fasten the joint screw and bearing nut; secure the screw connection with Loctite 241 and tighten to the specified tightening torque (see table).
- 5) Fasten the joint on the lamination bars using the screws; secure the screw connections with Loctite 241 and tighten to the specified tightening torque (see table).

## Disassembling the lock

- > Disassemble the lock to use the orthotic knee joint as a free-moving knee joint.
- 1) Unload the orthotic knee joint using the knee lever.
  - 2) Remove the lock rings and lock washers.
  - 3) Remove the axles.
  - 4) Pull off the knee lever.
  - 5) Push the lock wedge with the compression spring sideways out of the upper joint section.

## Optional: Installing the 17Y162=\* temporary switch

A temporary switch that can be subsequently installed in the orthotic joint is included with the orthotic joint. The temporary switch prevents the lock from engaging, e.g. while riding a therapy bike (see fig. 8).

- 1) Join the button to the tilt lever and insert it into the orthotic joint (see fig. 4, see fig. 5).
- 2) Insert the dowel pin into the orthotic joint and slide it through the hole in the tilt lever in order to secure the temporary switch (see fig. 6, see fig. 7).

## 6 Cleaning

Clean the orthosis joint promptly after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Rinse the orthosis joint with clean, fresh water.
- 2) Dry the orthosis joint with a cloth or allow it to air-dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

## 7 Maintenance

### INFORMATION

**The product may be exposed to increased loads by the patient.**

- Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

Do not use lubricants.

Spare parts are listed in the section "Components/design".

If you have any doubts, please contact the manufacturer.

## 8 Disposal

Dispose of the product in accordance with national regulations.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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