1C11 Terion K2

DE	Gebrauchsanweisung	3
EN	Instructions for use	10
FR	Instructions d'utilisation	17
IT	Istruzioni per l'uso	24
ES	Instrucciones de uso	31
PT	Manual de utilização	39
NL	Gebruiksaanwijzing	46
SV	Bruksanvisning	53
DA	Brugsanvisning	59
NO	Bruksanvisning	66
FI	Käyttöohje	73
PL	Instrukcja użytkowania	80
HU	Használati utasítás	87
CS	Návod k použití	94
RO	Instrucțiuni de utilizare	101
HR	Upute za uporabu	109
SL	Navodila za uporabo	116
SK	Návod na používanie	122
BG	Инструкция за употреба	129
TR	Kullanma talimatı	137
EL	Οδηγίες χρήσης	143
RU	Руководство по применению	151
JA	取扱説明書	159
ZH	使用说明书	165
KO	사용 설명서	171

1

INFORMATION

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INFORMATION

Date of last update: 2023-02-05

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 1C11 Terion K2 prosthetic foot was developed for the needs of users with mobility grades 1 and 2. The foot is suitable for walking on various surfaces.

The functional properties of the prosthetic foot are achieved through the spring element made of fibre-reinforced polymer and a functional foam base. The prosthetic foot permits perceptible plantar flexion at heel strike and a natural rollover movement.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

The MOBIS classification describes the mobility grade and body weight, and makes it easy to identify compatible components.

The product is recommended for mobility grade 1 (indoor walker) and mobility grade 2 (restricted outdoor walker).

- The maximum approved body weight is specified in the technical data (see page 16).

2.3 Environmental conditions

Storage and transport
Temperature range -20 °C to +60 °C, relative humidity 20 % to 90 %, no mechanical vibrations or impacts
Allowable environmental conditions
Temperature range: -10 °C to +45 °C
Chemicals/liquids: Fresh water, salt water, perspiration, urine, soapsuds, chlorine water
Moisture: Submersion: max. 1 h in 3 m depth, relative humidity: no restrictions
Solids: Dust, occasional contact with sand
Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage (see page 15).
Prohibited environmental conditions
Solids: Highly liquid-binding particles (e.g. talcum), continuous contact with sand
Chemicals/liquids: Acids, continuous use in liquid media

2.4 Lifetime

Prosthetic foot

The product was tested by the manufacturer with 2 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 3 years.

Footshell, protective sock

The product is a wear part, which means it is subject to normal wear and tear.

3 General safety instructions

Risk of injury and risk of product damage

- To avoid the risk of injury and product damage, do not use the product beyond the tested lifetime.

- ▶ To avoid the risk of injury and product damage, only use the product for a single patient.
- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 10).
- ▶ Note the combination possibilities/combo exclusions in the instructions for use of the products.

NOTICE!

Risk of product damage and limited functionality

- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

Signs of changes in or loss of functionality during use

Decreased forefoot resistance or changes in roll-over behaviour are noticeable indications of loss of functionality.

4 Scope of delivery

1C11 Terion K2		
Quantity	Designation	Reference number
1	Instructions for use	–
1	Prosthetic foot	–
1	Protective sock (size dependent)	22 to 28: SL=SPECTRA-SOCK2-7 29 to 30: SL=Spectra-Sock-7

5 Preparing the product for use

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

NOTICE

Altering the prosthetic foot or footshell

Premature wear resulting from damage to the product

- ▶ Do not alter the prosthetic foot or footshell.

INFORMATION

There is a plastic protector on the product's pyramid. It protects the connecting section from scratches during the alignment and fitting of the prosthesis.

- ▶ Remove the protector before the patient leaves the workshop/fitting area.

5.1 Applying/removing the footshell

INFORMATION

- ▶ Pull the protective sock over the prosthetic foot to prevent noises in the footshell.

▶ **CAUTION! Always use the prosthetic foot with the footshell.**

Apply or remove the footshell as described in the footshell instructions for use.

5.2 Bench Alignment

TT bench alignment

Bench alignment process

Required materials: 662M4 goniometer, 743S12 heel height measuring apparatus, 743A80 50:50 gauge, alignment apparatus (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly)

Assemble and align the prosthetic components in the alignment apparatus according to the following specifications:

①

Heel height: **Effective heel height** (shoe heel height - sole thickness of forefoot) + 5 mm

Bench alignment process	
2	Sagittal alignment: Alignment reference line on a-p markings of the connection cap/footshell (##06611)
3	Frontal alignment: Alignment reference line on m-l markings of the connection cap/footshell (##06611)
4	Connect the prosthetic foot and prosthetic socket using the chosen adapters. Be sure to follow the adapter instructions for use.
5	Sagittal alignment: Determine the centre of the prosthetic socket with the 50:50 gauge. Align the prosthetic socket centrally to the alignment reference line. Socket flexion: Individual residual limb flexion + 5°
6	Observe the abduction or adduction position.

TF bench alignment

- Observe the information in the prosthetic knee joint instructions for use.

5.3 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.
- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) may be requested from Ottobock.

5.4 Dynamic Trial Fitting

- Adapt the alignment of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) to ensure an optimum gait pattern.
- **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal plane is achieved when the leg begins to bear weight after the heel strike. Avoid medial movement of the knee joint. If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction. If the medial movement occurs in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.
- Remove the plastic protector from the pyramid after completing the dynamic fitting and the walking exercises.

5.5 Optional: Installing the foam cover

The foam cover sits between the prosthetic socket and prosthetic foot. It is cut longer in order to compensate for the movements of the prosthetic foot

and prosthetic knee joint. During flexion of the prosthetic knee joint, the foam cover undergoes posterior compression and anterior elongation. The foam cover should be stretched as little as possible in order to increase its service life. There is a connecting element (such as a connection plate, connection cap or connection cover) on the prosthetic foot.

> **Required materials:** degreasing cleaner (e.g. 634A58 isopropyl alcohol), 636N9 contact adhesive or 636W17 plastic adhesive

1) Measure the length of the foam cover on the prosthesis and add the length allowance.

TT prostheses: Distal allowance for movement of the prosthetic foot.

TF prostheses: Allowance proximal of the knee rotation point for flexion of the prosthetic knee joint and distal allowance for movement of the prosthetic foot.

2) Cut the pre-shaped foam cover to length and fit it in the proximal area on the prosthetic socket.

3) Pull the foam cover over the prosthesis.

4) Set the connecting element onto the footshell or prosthetic foot. Depending on the version, the connecting element engages in the edge or rests on the foot adapter.

5) Install the prosthetic foot on the prosthesis.

6) Mark the outer contour of the connecting element on the distal face of the foam cover.

7) Disassemble the prosthetic foot and remove the connecting element.

8) Clean the connecting element using a degreasing cleaner.

9) Glue the connecting element onto the distal face of the foam cover according to the marked outer contour.

10) Let the glue dry (approx. **10 minutes**).

11) Install the prosthetic foot and adapt the exterior cosmetic shape. Take into account compression caused by cosmetic stockings or SuperSkin.

6 Cleaning

1) Clean the product with a damp, soft cloth.

2) Dry the product with a soft cloth.

3) Allow to air dry in order to remove residual moisture.

7 Maintenance

▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.

▶ During regular inspection: Inspect the prosthesis for signs of wear and check its function.

- ▶ Conduct annual safety inspections.

8 Disposal

Do not dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

10 Technical data

Sizes [cm]	22	23	24	25	26	27	28	29	30
Heel height [mm]	10 ±5								
Footshell	narrow			standard					
System height with footshell [mm]	40		41		46		48	51	
Structural height [mm]	58		59		65		67	70	
Product weight without footshell [g]	237		283		325			423	
Product weight with footshell [g]	410	445	491	509	573	580	635	709	744
Max. body weight [kg]	100		150		175				
Mobility grade	1 + 2								

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