










28U11, 28U22, 28U23, 28U33

DE Gebrauchsanweisung	6
EN Instructions for use	12
FR Instructions d'utilisation	18
IT Istruzioni per l'uso	24
ES Instrucciones de uso	30
PT Manual de utilização	36
NL Gebruiksaanwijzing	42
SV Bruksanvisning	48
DA Brugsanvisning	54
NO Bruksanvisning	60
FI Käyttöohje	65
PL Instrukcja użytkowania	71
HU Használati utasítás	77
CS Návod k použití	83
RO Instrucțiuni de utilizare	89
HR Upute za uporabu	95
SK Návod na používanie	101
BG Инструкция за употреба	107
TR Kullanma talimatı	113
RU Руководство по применению	119
JA 取扱説明書	125
ZH 使用说明书	131

28U11* WalkOn	Artikelnummer / Article number	Seite / Side	Schuhgröße / Shoe size EU	Wadenhöhe / Calf height	
	28U11=L36-39	links/left	36-39	35 cm	13.8 in
	28U11=R36-39	rechts/right			
	28U11=L39-42	links/left	39-42	36 cm	14.2 in
	28U11=R39-42	rechts/right			
	28U11=L42-45	links/left	42-45	37,5 cm	14.8 in
	28U11=R42-45	rechts/right			
	28U11=L45-48	links/left	45-48	39 cm	15.4 in
	28U11=R45-48	rechts/right			

28U22* WalkOn Flex	Artikelnummer / Article number	Seite / Side	Schuhgröße / Shoe size EU	Wadenhöhe / Calf height	
	28U22=L36-39	links/left	36-39	32,5 cm	12.8 in
	28U22=R36-39	rechts/right			
	28U22=L39-42	links/left	39-42	35 cm	13.8 in
	28U22=R39-42	rechts/right			
	28U22=L42-45	links/left	42-45	37 cm	14.6 in
	28U22=R42-45	rechts/right			
	28U22=L45-48	links/left	45-48	39 cm	15.4 in
	28U22=R45-48	rechts/right			


28U22* WalkOn Flex Junior	Artikelnummer / Article number	Seite / Side	Schuhgröße / Shoe size EU	Wadenhöhe / Calf height	
	28U22=L24-27	links/left	24-27	20,5 cm	8.0 in
	28U22=R24-27	rechts/right			
	28U22=L27-30	links/left	27-30	24 cm	9.5 in
	28U22=R27-30	rechts/right			
	28U22=L30-33	links/left	30-33	27 cm	10.6 in
	28U22=R30-33	rechts/right			
	28U22=L33-36	links/left	33-36	29,5 cm	11.6 in
	28U22=R33-36	rechts/right			

28U23* WalkOn Trimable	Artikelnummer / Article number	Seite / Side	Schuhgröße / Shoe size EU	Wadenhöhe / Calf height	
	28U23=L36-39	links/left	36-39	35 cm	13.8 in
	28U23=R36-39	rechts/right			
	28U23=L39-42	links/left	39-42	36 cm	14.2 in
	28U23=R39-42	rechts/right			
	28U23=L42-45	links/left	42-45	37,5 cm	14.8 in
	28U23=R42-45	rechts/right			
	28U23=L45-48	links/left	45-48	39 cm	15.4 in
	28U23=R45-48	rechts/right			

28U33* WalkOn Lateral	Artikelnummer / Article number	Seite / Side	Schuhgröße / Shoe size EU	Wadenhöhe / Calf height	
	28U33=L36-39	links/left	36-39	32,5 cm	12.8 in
	28U33=R36-39	rechts/right			
	28U33=L39-42	links/left	39-42	35 cm	13.8 in
	28U33=R39-42	rechts/right			
	28U33=L42-45	links/left	42-45	37 cm	14.6 in
	28U33=R42-45	rechts/right			
	28U33=L45-48	links/left	45-48	39 cm	15.4 in
28U33=R45-48	rechts/right				

28U33* WalkOn Lateral Junior	Artikelnummer / Article number	Seite / Side	Schuhgröße / Shoe size EU	Wadenhöhe / Calf height	
	28U33=L24-27	links/left	24-27	20,5 cm	8.0 in
	28U33=R24-27	rechts/right			
	28U33=L27-30	links/left	27-30	24 cm	9.5 in
	28U33=R27-30	rechts/right			
	28U33=L30-33	links/left	30-33	27 cm	10.6 in
	28U33=R30-33	rechts/right			
	28U33=L33-36	links/left	33-36	29,5 cm	11.6 in
28U33=R33-36	rechts/right				

Material	28U11, 28U23, 28U33 Strukturteil/structural part: Carbonfaserverbundwerkstoff/carbon fibre composite; Wadengurt und Polster/calf strap and padding: Polyurethan/polyurethane
	28U22 Strukturteil/structural part: Glasfaserverbundwerkstoff/glas fibre composite; Wadengurt und Polster/calf strap and padding: Polyurethan/pPolyurethane
	Latex-frei: Dieses Produkt enthält kein Naturgummi./ Latex-free: This product does not contain natural rubber.

INFORMATION

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INFORMATION

Date of last update: 2023-05-09

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 28U11 WalkOn, 28U22 WalkOn Flex/WalkOn Flex Junior, 28U23 WalkOn Trimable and 28U33 WalkOn Lateral/WalkOn Lateral Junior ankle-foot orthoses.

2 Product overview

(see fig. 1)

Item	Designation	Item	Designation
1	Lower leg shell	4	Calf strap with padding
2	Orthosis	5	Sole
3	Y-hook-and-loop		

3 Intended use

3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

Ottobock recommends wearing the product with a stocking and not directly on the skin. A thin stocking made of a washable, moisture-absorbent material is suitable to wear under the product. To avoid pressure points, make sure that the stocking is free of wrinkles.

The orthosis must be used in accordance with the indications.

3.2 Indications

- Drop foot, without or with mild spasticity

Caused for example by:

- Peroneal paralysis
- Stroke

- Traumatic brain injury
- Multiple sclerosis
- Neuromuscular atrophy

Indications must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

- Moderate to severe spasticity in the lower leg
- Moderate to severe oedema
- Moderate to severe foot deformities

3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area to which the medical device will be applied; sensory and circulatory disorders in the lower leg and the foot area.

- Severe deviations of the leg axis (genu varum/valgum/recurvatum)
- Contractions at the knee and/or ankle joint

3.4 Lifetime

When used as intended, the product is designed for:

Junior **1 year** (shoe sizes EU 24-27, 27-30, 30-33, 33-36)

Adults **2 years** (shoe sizes EU 36-39, 39-42, 42-45, 45-48)

When the end of the lifetime is reached, have the product inspected for functionality by qualified personnel. Use of the product beyond the specified lifetime is the responsibility of the user.

The doctor decides the daily duration of use for the product.

3.5 Restrictions for use

The orthosis must not be used to prevent contractures.

The orthosis must not be used for sports activities that involve jumping, sudden movements or rapid step sequences (e.g. basketball, badminton, riding sports). Sports activities should generally be discussed with the patient.

3.6 Effects

- Supports dorsiflexion
- Promotes dynamic rollover
- Promotes a physiological gait pattern
- Reduces the risk of falling

3.7 Accessories

28U11, 28U22, 28U23

- 28Z10 Lateral pronation strap

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 General safety instructions

CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

5 Use

WARNING

Operation of motor vehicles

Risk of accident due to limited body function

- ▶ The product has not been approved by the manufacturer for the operation of motor vehicles.
- ▶ Observe the applicable legal regulations for the operation of motor vehicles and have your driving ability examined by an authorised test centre.

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

5.1 Size Selection

- ▶ Select the orthosis size based on the shoe size (see sizing table).

5.2 Adaptation

CAUTION

Improper moulding or application

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

NOTICE

Damage due to improper handling (heat treatment, drilling holes)

Loss of orthosis functionality

- ▶ Do not carry out any heat treatment on the orthosis. The orthosis is not thermoformable.
- ▶ Do not drill any holes in the bottom of the sole or the connecting elements as holes break the fibres and weaken the component.

NOTICE

Use without suitable footwear

Restriction of functionality due to insufficient stabilisation

- ▶ Only use the orthosis with closed footwear, taking the allowable effective heel height into account.

Attaching the calf strap to the lower leg shell

28U11, 28U23, 28U33

- 1) Open the strap and loosen it from the padding (see fig. 2, see fig. 3).
- 2) Pull the padding over the wings of the lower leg shell, positioning the padding on the inside and the strap on the outside of the shell (see fig. 4, see fig. 5, see fig. 6).
- 3) Pull the strap tight through the loop (see fig. 7, 4082282).

28U22

- 1) Open the strap and detach it from the hook-and-loop (see fig. 2).
- 2) Pull the padding over the long wing of the lower leg shell, positioning the padding on the inside and the strap on the outside of the shell (see fig. 4).
- 3) Pull the padding over the short wing of the lower leg shell (see fig. 5, see fig. 6).
- 4) Fasten the strap to the hook-and-loop on the short wing of the lower leg shell (see fig. 7, 4082282).

Selecting footwear

- ▶ Choose a sturdy lace-up shoe with solid heel cup to ensure the optimum effectiveness of the orthosis. The effective heel height should be between 0.5 cm and 1.5 cm.

Adapting the orthosis

- 1) Put the orthosis on the patient's foot.
- 2) If the insole in the patient's shoe can be removed, use it as a template to transfer the correct size to the foot plate or mark the contours of the foot directly on the plate with a pen (see fig. 9).
- 3) If using the insole, pay attention to the position of the insole on the foot plate. The orthosis should not come into contact with the ankle bone, and the movement of the patient's Achilles tendon should not be restricted.

- 4) **28U11:** Sand the foot plate up to the marking (see fig. 10).
28U22, 28U23, 28U33: Use scissors to trim the foot plate along the marking (see fig. 10).
NOTICE! When adjusting the width of the foot plate, do not remove more material than necessary on the lateral/medial side to prevent the spring from pressing against the ankle bone. If the plate needs to be narrower, remove material on both sides.
- 5) **28U11:** To adjust the length of the foot plate, focus on sanding the posterior section of the plate so that the spring is positioned behind the ankle and pressure on the ankle joint or ankle is prevented (see fig. 11).
28U22, 28U23, 28U33: To adjust the length of the foot plate, trim the posterior section of the plate so that the spring is positioned behind the ankle and pressure on the ankle joint or ankle is prevented (see fig. 11).
- 6) **Optional:** If the lower leg shell needs to be adapted to the patient, sand the shell.
- 7) If the orthosis has sharp edges, smooth down the sanded edges.
- 8) **For foot deformities:** Correct the foot position with an insole or a custom moulded device.
INFORMATION: If the foot cannot be corrected by an insole, the orthosis and a sturdy shoe, the WalkOn should not be used.

Shortening the calf strap

INFORMATION

If the strap is too long, unfavourable loading of the pad can loosen the orthosis and impair the function; therefore always adjust the strap.

- > Put the orthosis on the patient.
- 1) Determine the length of the strap. **INFORMATION: After shortening the strap, ensure that the Y-hook-and-loop closure is flush on the side of the strap and that the full surface of the Y-hook-and-loop adheres to the loop area of the pad (see fig. 8).**
- 2) Remove the Y hook-and-loop and shorten the strap using scissors.
- 3) Fasten the Y hook-and-loop to the strap.

28U22 WalkOn Flex Junior, 28U33 WalkOn Lateral Junior

- ▶ **CAUTION! The Y hook-and-loop of the child-size calf strap can be swallowed by young children.** Sew the Y-hook-and-loop to the calf strap after fitting.

Checking the fit of the orthosis

- ▶ Check that the orthosis is fitted properly before handing it over to the patient.
- ▶ Conduct trial walking with the patient. Testing on inclines, ramps and stairs is mandatory. Fine tune as required.

5.3 Application

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- 1) Remove the insole from the shoe on the side to be fitted.

- 2) Open the calf strap.
- 3) Insert the orthosis into the shoe.
- 4) Put the orthosis on together with the shoe.
- 5) Fasten the calf strap.

5.4 Removal

- 1) Open the calf strap.
- 2) Remove the shoe together with the orthosis.
- 3) **Optional:** Remove the orthosis from the shoe and re-insert the insole of the shoe.

6 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.
- ▶ Ottobock recommends cleaning the orthosis twice a week.

Textile component

▶ **NOTICE! Detergent residues may cause skin irritation and wear of material.**

- 1) Remove the textile component from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) **Optional:** Use a mesh laundry bag.
- 4) Wash the textile component in warm water at **30 °C** with a standard mild detergent. Rinse out detergent residues thoroughly.
- 5) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

Orthosis frame

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

7 Packaging

⚠ WARNING

Packaging in the hands of children

Risk of suffocation if the carrier bag is pulled over the mouth and nose, risk of strangulation by the cord

- ▶ Store the packaging out of reach of children.

The product is delivered in a carrier bag. The carrier bag is intended **exclusively** for transporting and storing the product.

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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