

## 28U70 GoOn

<b>DE</b>	Gebrauchsanweisung .....	4
<b>EN</b>	Instructions for use .....	9
<b>FR</b>	Instructions d'utilisation .....	14
<b>IT</b>	Istruzioni per l'uso .....	20
<b>ES</b>	Instrucciones de uso .....	25
<b>PT</b>	Manual de utilização .....	30
<b>NL</b>	Gebruiksaanwijzing .....	36
<b>SV</b>	Bruksanvisning .....	41
<b>DA</b>	Brugsanvisning .....	46
<b>NO</b>	Bruksanvisning .....	52
<b>FI</b>	Käyttöohje .....	56
<b>PL</b>	Instrukcja użytkowania .....	62
<b>HU</b>	Használati utasítás .....	67
<b>CS</b>	Návod k použití .....	73
<b>RO</b>	Instrucțiuni de utilizare .....	78
<b>HR</b>	Upute za uporabu .....	83
<b>SK</b>	Návod na používanie .....	88
<b>BG</b>	Инструкция за употреба .....	94
<b>TR</b>	Kullanma talimatı .....	99
<b>RU</b>	Руководство по применению .....	104
<b>JA</b>	取扱説明書 .....	110
<b>ZH</b>	使用说明书 .....	115

<b>Material</b>	PA, PES, PU, TPU, Federstahl / Spring steel
-----------------	---

## INFORMATION

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### INFORMATION

Date of last update: 2023-07-17

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use contain important information for fitting and setting up your 28U70 GoOn ankle foot orthosis.

## 2 Components (see fig. 1)

Item:	Designation	Item:	Designation
1	Upper calf strap	2	Lower calf strap

Item:	Designation	Item:	Designation
3	Foot strap	4	Orthosis
5	Hook-and-loop attachment point for foot strap	6	Y-hook-and-loop for calf straps

### 3 Intended use

#### 3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

#### 3.2 Indications

- Dorsiflexor weakness (e.g. through fatigue, balance problems, dementia, stroke, traumatic brain injury, multiple sclerosis, neuromuscular atrophy or peroneal paralysis)

Indications must be determined by the physician.

#### 3.3 Contraindications

##### 3.3.1 Absolute Contraindications

- Spasticity in the lower leg
- Moderate to severe oedema
- Leg ulcers
- Foot deformities
- Unstable ankle joint in the ML direction
- Unstable ankle joint in the AP and ML direction

##### 3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area to which the medical device will be applied; sensory and circulatory disorders in the lower leg and the foot area.

#### 3.4 Restrictions for use

The orthosis must not be used for the prevention of contractures (shortening) etc.

#### 3.5 Effects

The ankle foot orthosis provides lifting of the foot and limits plantar flexion.

#### 3.6 Lifetime

The product is designed for a maximum lifetime of **1 year**.

## 4 Safety

### 4.1 Explanation of warning symbols

**⚠ WARNING** Warning regarding possible serious risks of accident or injury.

**⚠ CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 4.2 General safety instructions

**⚠ CAUTION**

**Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

**⚠ CAUTION**

**Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

**⚠ CAUTION**

**Reuse on other persons and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

**NOTICE**

**Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

## 5 Handling

### **WARNING**

#### **Operation of motor vehicles**

Risk of accidents due to restricted body function

- ▶ Observe the applicable legal and insurance regulations for the operation of motor vehicles and have your driving ability examined and certified by an authorised agency.

### **INFORMATION**

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 5.1 Size Selection

The product is available in a universal size.

### 5.2 Application

#### **CAUTION**

#### **Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

#### **CAUTION**

#### **Improper moulding or application**

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

## NOTICE

### **Mechanical overload**

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

### **Applying the orthosis with a shoe**

> Closed shoes are used.

- 1) Put on the orthosis together with the shoe. Ensure that the short side of the orthosis is located at the instep (see fig. 2).
- 2) Move the foot slightly towards the back. This will cause the orthosis to move slightly in the direction of the tibia.
- 3) Fasten the calf straps to the medial or lateral side (see fig. 3). Fasten the upper calf strap to the proximal side.
- 4) Pass the calf strap around the calf and fasten it (see fig. 4).

### **Applying the orthosis without a shoe**

- 1) Place the orthosis against the foot. Ensure that the short side is located at the instep (see fig. 5).
- 2) Fasten the foot strap to the orthosis.
- 3) Pass the foot strap around the foot.
- 4) Affix the foot strap to the attachment point and then close the foot strap (see fig. 6).
- 5) Move the foot slightly towards the back. This will cause the orthosis to move slightly in the direction of the tibia.
- 6) Fasten the calf straps to the medial or lateral side (see fig. 7). Fasten the upper calf strap to the proximal side.
- 7) Pass the calf strap around the calf and fasten it (see fig. 8).

## **5.3 Cleaning**

Clean the orthosis regularly:

- 1) Remove the plastic splint.
- 2) Fasten all hook-and-loop closures.
- 3) Recommendation: Use a laundry bag.
- 4) Wash the product in 30 °C warm water with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 5) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

6) Reapply the plastic splint.

## 6 Disposal

Dispose of the product in accordance with national regulations.

## 7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 7.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

### 7.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 7.3 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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