## 17B203=\* E-MAG Active

<b>DE</b> Gebrauchsanweisung (Benutzer) .....	3
<b>EN</b> Instructions for use (user) .....	10
<b>FR</b> Instructions d'utilisation (Utilisateur) .....	18
<b>IT</b> Istruzioni per l'uso (Utilizzatore) .....	25
<b>ES</b> Instrucciones de uso (Usuario) .....	33
<b>PT</b> Manual de utilização (Usuário) .....	40
<b>NL</b> Gebruiksaanwijzing (Gebruiker) .....	48
<b>SV</b> Bruksanvisning (Användare) .....	55
<b>DA</b> Brugsanvisning (Bruger) .....	62
<b>NO</b> Bruksanvisning (Bruker) .....	70
<b>FI</b> Käyttöohje (Käyttäjä) .....	77
<b>HU</b> Használati utasítás (felhasználó) .....	84
<b>CS</b> Návod k použití (Uživatel) .....	92
<b>HR</b> Upute za uporabu (Korisnik) .....	99
<b>SK</b> Návod na používanie (Uživateľ) .....	107
<b>RU</b> Руководство по применению (Пользователь) .....	114
<b>ZH</b> 使用说明书 (用户) .....	123

# 1 Vorwort

Deutsch

## INFORMATION

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### INFORMATION

Date of last update: 2024-01-10

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the use of the 17B203 E-MAG Active orthotic knee joint.

## 2 Product description

### 2.1 Available sizes

Article number	Material	Pre-flexion	17B206 medial support	Max. body weight
17B203=L/R, 17B203=L/R-7.5	Steel	5° and 7.5°	Without medial support	Up to 85 kg
17B203=L/R, 17B203=L/R-7.5	Steel	5° and 7.5°	With medial support	Up to 100 kg

### 2.2 Construction and Function


The E-MAG Active is an orthotic knee joint with stance phase control which offers you a free swing phase and locks the orthotic knee joint for the stance phase prior to heel strike. The orthotic knee joint determines the position of your leg with the help of the electronics and enables the swing phase. Read the 646D1043 "Quick start guide" for more information.

The E-MAG Active is an orthotic knee joint with a concealed locking mechanism, which means that the joint is less susceptible to dirt and also prevents clothing from becoming trapped. Not only is this gentle on your clothing but, above all, offers you enhanced safety.

The E-MAG Active has an additional locking function (PreLock) that is activated at 15° knee flexion. This secures the orthotic knee joint before the swing phase ends. The orthotic knee joint is then fully locked when it reaches the extension stop.

Two buttons allow you to electronically unlock the orthotic knee joint (once) or engage the lock (permanently). In addition, the (temporary) mechanical unlocking mechanism supports a wider area of application, for example when cycling. These instructions for use describe the functions and operation of the E-MAG Active as well as how to charge it and care for it.

### 2.3 Product overview

(see fig. 1), item	Component
1	Orthotic knee joint
2	Unlocking mechanism
3	Battery
4	Electronics

## 3 Intended use

### 3.1 Indications for use

The product is intended **solely** for orthotic fittings of the lower limbs.

### 3.2 Conditions of use

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.), sports activities that involve jumping, sudden movements or fast sequences of steps (e.g. basketball, badminton, sporty riding).

Possible sports activities must be discussed with the O&P professional.

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

### 3.3 Environmental conditions

<b>Storage temperature range</b>	-20 °C to +70 °C (-4 °F to 176 °F)
<b>Operating temperature range</b>	-15 °C to +50 °C (-5 °F to 122 °F)
<b>Relative humidity for both ranges</b>	15% to 93%
<b>Power supply for joint</b>	NiMh battery, 4.8 V nominal voltage
<b>Range of a battery pack</b>	Approx. 5,000 steps, corresponds to about 5 km

### 3.4 Indications

- Partial or total paralysis of the leg muscles

Indications must be determined by the physician.

### 3.5 Contraindications

**Contraindications:**

- Knee flexion contractures in excess of 15° that prevent the knee joint lock from being relieved
- Hip flexion contractures
- Uncontrollable spasticity
- Ischial support
- Non-physiological deviation in the frontal plane

### 3.6 Lifetime

The product is designed for a lifetime of **3 years** when used as intended and assembled professionally.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 4.2 General safety instructions

<p> <b>CAUTION</b></p> <p><b>Mechanical damage to the product</b> Injuries due to changes in or loss of functionality</p> <ul style="list-style-type: none"> <li>▶ Use caution when working with the product.</li> <li>▶ Check the product for proper function and readiness for use.</li> <li>▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.</li> </ul>
---

<p> <b>CAUTION</b></p> <p><b>Exposure of the product to unsuitable environmental conditions</b> Patient injury, damage, brittleness or destruction due to improper handling</p> <ul style="list-style-type: none"> <li>▶ Do not expose the product to condensing ambient humidity or liquids.</li> <li>▶ Do not expose the product to abrasive substances (e.g. sand, dust).</li> <li>▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).</li> </ul>
---

**⚠ CAUTION****Excessive strain on load-bearing components**

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

**⚠ CAUTION****Incorrect alignment, assembly or adjustment**

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

**⚠ CAUTION****Improper handling**

Risk of falling due to accidentally opening or engaging the lock

- ▶ Before placing weight on the orthosis, verify that the lock is in the desired position.
- ▶ Avoid accidentally opening or engaging the lock.
- ▶ Instruct the patient in the correct handling of the orthotic joints.

**⚠ CAUTION****Knee not fully extended at heel strike**

Risk of injury due to unlocked orthosis joint

- ▶ Fully extend the knee prior to heel strike.

**⚠ CAUTION****Risk of injury due to electromagnetic interference**

Malfunctions may occur in the vicinity of high-voltage power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores). This can cause the patient to fall.

Avoid proximity to strong sources of magnetic and electric interference (e.g. transformer stations, transmitters).

**⚠ CAUTION****Risk of falling due to improper operation of the service switches**

The A and T service switches may only be adjusted by qualified personnel.

**⚠ CAUTION****Risk of falling while using the control unit**

To avoid the risk of falling, you may only activate the locking and sitting function while standing and not while walking.

**⚠ CAUTION****Risk of falling due to use of grease (or other lubricants)**

This could contaminate the joint and jam the lock. The bearing washers and joint locking mechanism should never be lubricated.

**NOTICE****Risk of damage due to incorrect environmental conditions**

Solid particles and fluid in the joint, lock or electronics of the E-MAG Active may lead to malfunctioning. Do not expose the components to excessive smoke or dust in the environment, mechanical vibrations or shocks.

## 5 Preparing the product for use

### INFORMATION

The E-MAG Active can only be operated with the intended 317B20 battery. It was developed especially for the E-MAG Active. Both the receptacle and lock unit are designed to fit with the form of this battery.

### INFORMATION

Before putting the E-MAG Active into operation, read the instructions for the 317B20 battery and the 317L20 battery charger.

This section will describe how the E-MAG Active can be activated/deactivated and how it can be reactivated from the energy-saving mode. The section also explains how the electronic system starts and which signals the system emits during activation.

### 5.1 Charging the battery

### INFORMATION

After charging, the E-MAG Active can be used for around 5,000 steps with one battery charge. Nevertheless, charging the battery daily (ideally at night) is recommended.

### INFORMATION

The battery should always be removed when the orthosis is not in use.

### NOTICE

#### Continuous charging of the battery

Damage to the battery

- ▶ Avoid uninterrupted charging for several days

### NOTICE

#### Charging the battery with soiled or damaged contacts

Loss of battery functionality due to insufficient charging functionality

- ▶ Ensure that the contacts are always clean and free of grease
- ▶ Clean all electrical contacts of the battery regularly with a lint-free cloth
- ▶ Take care to avoid damaging the contacts with pointed or sharp objects

- 1) Insert one or two 317B20 batteries into the charger (see fig. 2). **INFORMATION: The corresponding LED lights up for each battery. The charging time for a discharged battery is 8-10 hours.**
- 2) Take the battery out of the battery charger after charging.

## 5.2 Activating/deactivating the E-MAG Active

### ⚠ CAUTION

#### Risk of falling due to defective control device/rechargeable battery

Injuries due to changes in or loss of functionality

- ▶ Pay attention to the signals when activating.
- ▶ If signals are absent, do not continue using the product and have it inspected by authorised, qualified personnel.

### INFORMATION

Failure to unlock the orthotic knee joint for more than one hour automatically causes the control electronics to switch to standby mode. The system is reactivated by simply pressing the top button on the control electronics.

Activation is performed by inserting the rechargeable battery into the receptacle. Every time the E-MAG Active is activated, the control electronics unit runs a short self-test that emits the following signals:

Signal	Visual feedback (LEDs)	Audible feedback	Tone/length
System ready > knee joint locked	All LEDs flash briefly (white)		Rising
System not ready > knee joint flexed (open)	Both LEDs flash continuously (white and green)		Tone-pause-tone series

If signals are absent, remove and reinsert the rechargeable battery. If the signals fail, there is a fault.

When the E-MAG Active is flexed, an audible signal-pause-audible signal sequence sounds as a warning signal until the E-MAG Active is extended/locked. Deactivation is performed by removing the battery.

## 6 Use

This section explains how to operate the knee joint control unit by button and tells you about the warnings used by the feedback system of the E-MAG Active to inform you of malfunctions.

### 6.1 Control and use of 17B203 electronics


### ⚠ CAUTION

#### Risk of falling due to unlocking/locking

Injuries due to changes in or loss of functionality

- ▶ Only operate the single unlocking mechanism or permanent locking mechanism when standing securely.

The electronic control unit measures the position of the orthosis during walking. It unlocks the orthotic knee joint before toe-off and locks it prior to heel strike at the end of the swing phase so the user can step down safely. Your O&P professional adapted the electronic control unit of the E-MAG Active to your needs during the initial fitting. Two manual switching possibilities are available to you to give you full control over the E-MAG Active (see fig. 3):


### 6.1.1 Electronic locking (permanent)

If you feel unsafe, you can lock the joint at any time by pressing the top button for locking the orthotic knee joint. This may be the case when walking on a slope or if you want to stand for an extended period of time without focusing on the electronic control device. In order to do so, briefly press the top button twice in quick succession (like double-clicking a computer mouse). This function is confirmed by an audible signal. The locking function is deactivated by pressing the top button twice in quick succession until you hear an audible signal.

### 6.1.2 Electronic unlocking (once)

The patient can use the lower button to unlock the orthotic knee joint once (e.g. in order to sit down). In order to do so, briefly press the lower button twice in quick succession (like double-clicking a computer mouse). This function is **not** confirmed by an audible signal.

### 6.1.3 Mechanical unlocking (temporary)

For certain activities (such as cycling), it may be beneficial to temporarily unlock the orthotic knee joint. Mechanical unlocking of the orthotic knee joint is indicated by the symbol of a cyclist (switch to "🚲").

The unlocking feature needs to be cancelled in order to restore the functionality of the orthotic knee joint. In order to secure the orthotic knee joint again, the switch is operated once more in the opposite direction. The normal state of the joint is indicated by the symbol of a walking person (switch to "🚶").

## 6.2 Warning signals, use of 17B203 electronics

### 6.2.1 Battery charge level

**INFORMATION**  
The electronics will turn off if the battery is discharged completely. The lock can be released manually.

Should the battery charge level fall below the preset charge capacity, the electronics will emit the warning signals listed here:

Signal	Visual feedback (LEDs)	Acoustic feedback	Tone/length
Battery low	Alternate flashing (white and red)		Series of high-low tones/individual signals
Battery empty			Falling/long

The warning signal is repeated every 15 minutes. Replacing the battery immediately is highly recommended. How quickly the battery is discharged depends on the ambient temperature: The colder it is, the faster the battery will be drained. The electronics will turn off if the battery power is insufficient. If this occurs, the orthotic knee joint can only be unlocked manually (see fig. 4).

## 6.2.2 Failure of the electronics

The E-MAG Active is not vulnerable to failure. However, should disruptions occur, the electronics recognise the failure and send a warning signal.

Signal	Visual feedback (LEDs)	Audible feedback	Tone/length
Electronics malfunction	Continuously lit (all LEDs lit, except green)		High/continuous tone

## 7 Cleaning

Clean the orthosis regularly:

### Orthosis frame

- 1) Clean the orthosis or E-MAG Active components from the outside with a slightly damp cloth.
- 2) Dry the orthosis or E-MAG Active components with a cloth or allow to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

### Pads

- 1) Remove the pad from the orthosis.
- 2) Wash by hand in warm water at 30 °C with pH-neutral soap and rinse thoroughly.
- 3) Allow the pad to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).
- 4) Refit the pad in the orthosis.
- 5) Once it has been refitted, wipe the pad with a damp cloth and dry.

## 8 Maintenance

Have your O&P professional check your E-MAG Active for wear and tear as well as proper functioning of the electronic components every six months.

## 9 Disposal

Dispose of the product in accordance with national regulations.

## 10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

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