1C56 Taleo Adjust


DE	Gebrauchsanweisung	7
EN	Instructions for use	16
FR	Instructions d'utilisation	25
IT	Istruzioni per l'uso	35
ES	Instrucciones de uso	45
PT	Manual de utilização	55
NL	Gebruiksaanwijzing	65
SV	Bruksanvisning	75
DA	Brugsanvisning	84
NO	Bruksanvisning	93
FI	Käyttöohje	102
PL	Instrukcja użytkowania	111
HU	Használati utasítás	121
CS	Návod k použití	130
RO	Instrucțiuni de utilizare	139
HR	Upute za uporabu	149
SL	Navodila za uporabo	158
SK	Návod na používanie	167
BG	Инструкция за употреба	177
TR	Kullanma talimatı	187
EL	Οδηγίες χρήσης	196
RU	Руководство по применению	207
JA	取扱説明書	217
ZH	使用说明书	226
KO	사용 설명서	233

1

2

Bench alignment, transtibial

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INFORMATION

Date of last update: 2024-03-19

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The heel height of the Taleo Adjust prosthetic foot 1C56 is infinitely adjustable (**0 mm** to **70 mm**) to enable shoes with different heel heights to be worn.

Size 27 and 28: Set the prosthetic foot to a maximum heel height of 25 mm.

Carbon and polymer spring elements permit perceptible plantar flexion at heel strike, a natural rollover and high energy return. A removable heel wedge absorbs the heel strike.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

The MOBIS classification describes the mobility grade and body weight, and makes it easy to identify compatible components.

The product is recommended for mobility grade 2 (restricted outdoor walker) and mobility grade 3 (unrestricted outdoor walker).

The following table shows the appropriate spring stiffness for the level of load exerted by the user on the prosthetic foot.

TIP: Users requiring extra stability can test the stiffness for a higher load level.

Spring stiffness relative to body weight and load level		
Body weight [kg]	Low load level	Medium load level
Up to 51	1	2
52 to 58	2	3
59 to 67	3	4
68 to 77	4	5
78 to 88	5	6
89 to 100	6	7
101 to 115	7	

2.3 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2.4 Environmental conditions

Allowable environmental conditions
Temperature range: -10 °C to +45 °C

Allowable environmental conditions
Chemicals/liquids: fresh water, soap solution, chlorine water, salt water ≤ 3.5 % salt content
Moisture: Submersion: max. 1 h in 2 m depth, relative humidity: no restrictions
Solids: Dust, occasional contact with sand
Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage (see page 24).
Prohibited environmental conditions
Solids: highly liquid-binding particles (e.g. talcum), dust in high concentrations (e.g. construction site), intensive contact with sand
Chemicals/liquids: perspiration, urine, acids, continuous use in liquid media
Storage and transport
Storage temperature: -20 °C to $+60$ °C, relative humidity: 20 % to 90 %, no mechanical vibrations or impacts

2.5 Lifetime

Prosthetic foot

Depending on the patient's activity level, the maximum lifetime of the product is 3 years.

Footshell, protective sock

The product is a wear part, which means it is subject to normal wear and tear.

3 General safety instructions

CAUTION! Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ Note the combination possibilities/comboination exclusions in the instructions for use of the products.
- ▶ Observe the maximum lifetime of the product.
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not expose the product to excessive strain (section: "Area of application" see page 17) or to prohibited environmental conditions (see page 17).
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ The product may only be used by one person; reuse on another person is prohibited.

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover) or delamination of the spring are indications of loss of functionality. Play in the hydraulic ankle joint, oil leakage and unusual noises may indicate a loss of functionality.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	-
1	Prosthetic foot	-
1	Protective sock	SL=Spectra-Sock-7
1	Heel wedge set	2F50=*

Spare parts/accessories

Designation	Reference number
Footshell	2C8=*
Connection cover	2C21=*
Footshell, low-cut	2C15=*-L

5 Preparing the product for use

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- Observe the alignment, assembly and adjustment instructions.

CAUTION

Reaching into the ankle joint

Risk of injury due to pinching of the fingers

- Do not reach into the ankle joint area.

NOTICE

Altering the prosthetic foot or footshell

Premature wear resulting from damage to the product

- Do not alter the prosthetic foot or footshell.

INFORMATION

A plastic part protects the connecting section of the product from scratches during the alignment and fitting.

- ▶ Remove the plastic part before the patient leaves the fitting area.

5.1 Applying/removing the footshell

INFORMATION

- ▶ Pull the protective sock over the prosthetic foot to prevent noises in the footshell.

▶ CAUTION! Always use the prosthetic foot with the footshell.

Apply or remove the footshell as described in the footshell instructions for use.

5.2 Bench Alignment

INFORMATION

Use the prosthetic foot only with an installed heel wedge.

INFORMATION

If the user needs more stability: Move the prosthetic foot approx. 1 cm further in the anterior direction. This lengthens the forefoot and increases the knee securing effect.

Bench alignment, transtibial prosthesis

Bench alignment process

Required materials: 662M4 goniometer, 743S12 heel height measuring device, 743A80 50:50 gauge, alignment apparatus (e.g. 743A200 PROS.A. Assembly)

Assemble and align the prosthetic components in the alignment tool:

- 1 Heel height: **0 mm**
Bring the prosthetic foot into maximum dorsiflexion (see page 23) so the complete adjustment range can be accessed.
- 2 **Sagittal plane**
Footshell with markings: Place alignment reference line on markings (see fig. 2)
Footshell without markings: Place the alignment reference line in a posterior position in relation to the middle of the foot (**25 mm** for size 22 to 24, **30 mm** for size 25 to 28)

Bench alignment process	
3	<p>Frontal plane Footshell with markings: Place alignment reference line on markings (see fig. 2) Footshell without markings: Set exterior rotation to approx. 5° (alignment reference line in the middle of the second toe)</p>
4	Install the structural components between the prosthetic foot and prosthetic socket according to their instructions for use.
5	<p>Sagittal plane Determine the centre of the prosthetic socket proximally and distally with the 50:50 gauge and mark the centre line. Mark the socket reference point on the centre line: at the height of the centre of the patella. The alignment reference line runs vertically through the socket reference point. Adjust the socket flexion by rotating around the socket reference point: individual residual limb flexion + 5°</p>
6	Observe the abduction or adduction position.

Bench alignment, thigh/knee disarticulation prosthesis

1) Heel height: **0 mm**

Bring the prosthetic foot into maximum dorsiflexion (see page 23) so the complete adjustment range can be accessed.

2) **Sagittal plane**

→ Footshell with markings: Place alignment reference line on markings.

→ Footshell without markings: Place the alignment reference line in a posterior position in relation to the middle of the foot (**25 mm for size 22 to 24, 30 mm for size 25 to 28**)

3) **Frontal plane**

→ Footshell with markings: Place alignment reference line on markings

→ Footshell without markings: Set exterior rotation to approx. 5° (alignment reference line in the middle of the second toe)

► Observe the instructions for use for the prosthetic knee joint.

5.3 Static Alignment

The static alignment of the Taleo Adjust consists of two stages. You have to perform both stages in order to set the minimum and maximum heel height.

Ottobock recommends checking the alignment of the prosthesis during both stages using the L.A.S.A.R. Posture.

1) **Without footwear (see fig. 3):**

- To perform the static alignment without footwear, bring the prosthetic foot into maximum dorsiflexion.
- Use the set screws to make any necessary adjustments.

2) **With everyday footwear (see fig. 4):**

- Carry out the static alignment with the everyday footwear.
- Only adjust the angle of the sagittal plane by changing the heel height (see page 23), not using the set screws. Using the set screws changes the maximum/minimum adjustable heel height.

The alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) can be requested from Ottobock as needed.

5.4 Dynamic Trial Fitting

The optimal gait pattern is developed during dynamic fitting. The prosthesis is aligned optimally in the frontal plane and the sagittal plane for this purpose.

- ▶ Use everyday footwear for dynamic fitting.
- ▶ **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal planes is achieved after the leg begins to bear weight following the heel strike. Avoid medial movement of the knee joint.
 - If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction.
 - If the knee joint moves in the medial direction in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.
- ▶ Once the gait pattern has been optimised with everyday footwear, try additional heel heights. The heel height affects rollover as it changes the angle at which the prosthetic foot is loaded.
- ▶ Remove the plastic part from the connecting section after dynamic fitting and the walking exercises are complete.

5.4.1 Replacing the heel wedge

The behaviour of the prosthetic foot at heel strike and during heel contact in the mid-stance phase can be adapted by replacing the heel wedge. Heel wedges of various hardness are included in the scope of delivery.

Heel wedge hardness: The colour of the heel wedge indicates the hardness (3648805). Ottobock recommends starting with the preinstalled heel wedge.

- 1) Bend the prosthetic foot slightly and remove the installed heel wedge.

- 2) Align the other heel wedge so the Ottobock lettering is right side up and the tip points in the anterior direction.
- 3) Insert the heel wedge into the prosthetic foot (Fersenkeil einsetzen Taleo). To do so, first insert the tip into the front part of the frame ❶, then push the heel wedge into the rear of the frame ❷.

6 Use

6.1 Setting the heel height

INFORMATION

First-time use or after an extended break: there may be stiffness when operating the lock axle and/or adjusting the ankle angle.

- To loosen: slide the lock axle back and forth several times and move the ankle joint through its entire scope of movement.

The heel height of the prosthetic foot can be smoothly adjusted from **0 mm** to **70 mm**.

Size 27 and 28: Set the prosthetic foot to a maximum heel height of 25 mm.

Adjustment method

Set the heel height while seated. Then hold on to a suitable object, stand up and test whether the set heel height is appropriate or needs to be corrected again.

Setting the heel height

- 1) Note how your prosthesis feels when you walk barefoot.
- 2) Sit down.
- 3) Note the angle of your prosthetic leg when you sit barefoot and the prosthetic foot is flat on the floor.
- 4) Put on the shoes you want to wear now.
- 5) **CAUTION! Do not put any weight on the prosthesis while setting the heel height. There is a risk of injury given the freely movable nature of the prosthetic foot.**
Unlock the ankle joint: Push the lock axis in from the medial direction as far as it will go (see fig. 9). The ankle joint can now be moved.
- 6) Push the tip of the prosthetic foot towards the sole of the foot as far as it will go (plantar flexion).
- 7) Place your prosthetic leg on the floor. The shoe should be flat on the floor and the prosthetic leg should have the same angle as it has when you are sitting barefoot.

- 8) **Lock the ankle joint:** Push in the lateral lock axis as far as it will go (see fig. 10). The ankle joint is now locked and can bear weight once more.
- 9) Hold on to a suitable object, stand up and test whether the heel height setting is appropriate.

7 Cleaning

- > **Permissible cleaning agent:** pH-neutral soap (e.g.453H10 Derma Clean)
- 1) **NOTICE! To avoid product damage, only use the permissible cleaning agents.**
Clean the prosthetic foot and footshell with clear fresh water and pH-neutral soap.
- 2) **If present:** Use a toothpick to remove dirt from the contours that drain off the adapters or footshell and rinse.
- 3) Rinse the soap away with clear fresh water. In doing so, rinse the footshell until all dirt has been removed.
- 4) **In case of heavy soiling and occurring sounds:** Pull the footshell and protective sock off the prosthetic foot. Clean all products individually with water.
- 5) Dry the product with a soft cloth.
- 6) Allow to air dry in order to remove residual moisture.

8 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ During regular inspection: Inspect the prosthesis for signs of wear and check its function.
- ▶ Conduct annual safety inspections.

9 Disposal

Do not dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

11 Technical data

Sizes [cm]	22	23	24	25	26	27	28
Heel height [mm]	0 to 70					0 to 25	
System height with 2C8* [mm]	92		95		100		
System height with 2C15*-L [mm]	87				100	-	
Build height with 2C8* [mm]	110		113		118		
Build height with 2C15*-L [mm]	105				118	-	
Weight with 2C8* [g]	785	825	845	890	930	955	1000
Weight with 2C15*-L [g]	740	770	810	850	900	-	
Max. body weight at low load level [kg]	88	100		115			
Max. body weight at medium load level [kg]	77	88		100			

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