21Y2=M10

DE Gebrauchsanweisung	10
EN Instructions for use	13

**Order your
free printed copy**

order-ifu@ottobock.com

Document: 647G2026 Version: 01

<https://product-documents.ottobock.com/IFU/INT/21Y2=M10/647G2026/01/O/S/F>

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
Basic UDI-DI: 403276700000000021Y2V2

1

Drill and deburr a hole with a diameter of 8.5 mm.
Cut an M10 thread.

2


Apply adhesive to the thread.
Check socket: 617H46
Definitive socket: 636W28 or 617H46

3

Carefully screw in until the O-ring seals against the outside of the prosthetic socket.

4

Line up the markings on the cover and base.
Place the cover on the base.

5

Turned clockwise to the stop: valve locked, air can only come out

Turned anticlockwise to the stop: bypass is open, air can flow in and out

Middle position between locked and open: cover can be removed or attached

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INFORMATION

Date of last update: 2024-04-04

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These are the instructions for use for the 21Y2=M10 outlet valve.

The outlet valve is used in combination with a sealing sleeve or a sealing liner to generate a vacuum in the prosthetic socket. This reduces friction and shear forces during the swing phase. During the stance phase, the body weight pushes the air in the prosthetic socket out of the socket via a check valve.

The outlet valve can be opened to remove the prosthetic socket or increase comfort during extended sitting.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

- The product can be used for vacuum-formed and laminated prosthetic sockets.
- Only use the product in combination with a liner.

2.3 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2.4 Environmental conditions

Storage and transport

Storage temperature: -20 °C to +60 °C, relative humidity: 20 % to 90 %, no mechanical vibrations or impacts

Allowable environmental conditions

Temperature during use: -10 °C to +45 °C

Chemicals/liquids: Fresh water, salt water, perspiration, urine, soapsuds, chlorine water

Moisture: Submersion: max. 1 h in 3 m depth, relative humidity: no restrictions

Solids: Dust, occasional contact with sand

Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and tear and damage (see page 15).

Prohibited environmental conditions

Solids: Highly liquid-binding particles (e.g. talcum), continuous contact with sand

Chemicals/liquids: Acids, continuous use in liquid media

2.5 Lifetime

The manufacturer has specified a maximum lifetime of 2 years for the product.

3 General safety instructions

CAUTION! Risk of injury and risk of product damage

- ▶ To avoid damage to the product, observe the alignment and assembly instructions.
- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ Note the combination possibilities/combination exclusions in the instructions for use of the products.
- ▶ Observe the maximum lifetime of the product.
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not expose the product to excessive strain (section: "Area of application" see page 13) or to prohibited environmental conditions (see page 13).
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

CAUTION! Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may only be used by one person; reuse on another person is prohibited.
- ▶ Observe the cleaning instructions.

NOTICE! Danger of product damage and restricted functionality

- ▶ Do not allow the product to come into contact with oils, ointments or lotions containing oils or acids.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	-
1	Valve cover	21Z42=1
1	Valve base	-

5 Preparing the product for use

The outlet valve is assembled in the distal area of the prosthetic socket, offset to the posterior or lateral. The base of the outlet valve is screwed into the prosthetic socket, and the cover is positioned on the base.

5.1 Assembling the valve

Preparing for assembly

- > **Required materials:** drill bit 8.5 mm, M10 screw tap set
- 1) Drill a hole (**8.5 mm**) in the prosthetic socket with the posterior or lateral offset in the distal area. (see fig. 1).
- 2) Carefully deburr the bore hole on the inside and the outside of the prosthetic socket. Remove chips and dust.
- 3) Cut an M10 thread into the bore (see fig. 1).
- 4) **Test assembly:** Carefully screw the valve base into the prosthetic socket until the O-ring touches the prosthetic socket.
 - **The O-ring does not lie flush:** Sand the prosthetic socket slightly flat in this area.
 - **The thread of the valve base projects into the prosthetic socket.** Carefully sand the thread a bit shorter.

Assembling the valve base

To seal the valve base, treat the thread with silicone bonding agent or special adhesive (**only for the definitive socket**) prior to assembly.

- > **Required materials:** 617H46 silicone bonding agent or 636W28 or 636W18 and 636W19 Ottobock special adhesive
- > **Prerequisite:** The thread of the valve base is clean and free of grease.
- 1) Apply a bit of silicone bonding agent or special adhesive to the thread of the valve base (see fig. 2).
- 2) Carefully screw the valve base into the prosthetic socket (see fig. 3). The O-ring has to seal on the prosthetic socket without being pushed away.
- 3) Allow the material to cure.

Assembling the valve cover

Assembly or disassembly of the valve cover is only possible in one position. Markings on the valve base and valve cover indicate the position.

- 1) Align the markings one above the other and set the valve cover onto the valve base (see fig. 4, see fig. 5).

TIP: The valve cover covers the marking on the base. Transfer the marking of the base to the prosthetic socket with a pen so it can be easily found again.

- 2) Turn (lock) the valve cover clockwise to the stop.

6 Use

NOTICE! If you turn the valve cover to the middle between “open” and “locked”, you can remove the valve cover from the base.

Open

- ▶ Turn the valve cover counter-clockwise to the stop.
- The valve is open, and air can flow into the prosthetic socket.

Locking

- ▶ Turn the valve cover clockwise to the stop.
- The valve is locked; air can only come out of the prosthetic socket.

7 Cleaning

Clean the valve if it is dirty or no longer sealed reliably.

- 1) Turn the valve cover to the middle between “open” and “locked”, and pull it off the base.
- 2) Flush the valve with compressed air or water from the underside of the valve cover.
- 3) Rinse the valve base with water or compressed air.
- 4) Align the markings on the valve base and valve cover on top of each other, and set the valve cover back on the valve base.
- 5) Turn the valve cover clockwise to the stop to lock the valve.

8 Disposal

Do not dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ During regular inspection: Inspect the prosthesis for signs of wear and check its function.
- ▶ During regular inspection: Clean the outlet valve.
- ▶ Conduct annual safety inspections.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

Ottobock SE & Co. KGaA
Max-Näder-Straße 15 · 37115 Duderstadt · Germany
T +49 5527 848-0 · F +49 5527 848-3360
healthcare@ottobock.de · www.ottobock.com