RS2-00 Runway, RS2-K1 Runway HX

DE	Gebrauchsanweisung	5
EN	Instructions for use	15
FR	Instructions d'utilisation	24
IT	Istruzioni per l'uso	34
ES	Instrucciones de uso	45
PT	Manual de utilização	55
NL	Gebruiksaanwijzing	65
SV	Bruksanvisning	75
DA	Brugsanvisning	85
NO	Bruksanvisning	94
FI	Käyttöohje	104
PL	Instrukcja użytkowania	113
HU	Használati utasítás	123
CS	Návod k použití	133
RO	Instrucțiuni de utilizare	143
HR	Upute za uporabu	153
SL	Navodila za uporabo	163
SK	Návod na používanie	173
BG	Инструкция за употреба	183
TR	Kullanma talimatı	193
EL	Οδηγίες χρήσης	203
RU	Руководство по применению	214
JA	取扱説明書	224
ZH	使用说明书	234
KO	사용 설명서	241

ISO 10328 - P n - m kg*)

*) Body mass limit not to be exceeded!

Stiffness	Weight (kg)	Label text
1	59	ISO 10328-P3-59-kg
2	68	ISO 10328-P3-68-kg
3	77	ISO 10328-P3-77-kg
4	88	ISO 10328-P4-88-kg
5	100	ISO 10328-P5-100-kg
6	116	ISO 10328-P5-116-kg


1

2

3

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INFORMATION

Date of last update: 2024-04-16

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Construction and Function

This document applies to the following products: RS2-00 Runway, RS2-K1 Runway HX

The prosthetic foot has carbon spring elements. The heel stiffness can be increased with a heel wedge.

The heel height of the prosthetic foot can be adjusted from **0 cm** to **5 cm**. Thanks to the design of the mechanism, the ratio between the forefoot lever and hindfoot lever remains constant. The user may change the settings independently.

- **RS2-00 Runway:** Adjustment in ten stages by pressing a release button.
- **RS2-K1 Runway HX:** Continuous adjustment using an Allen key.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Not compatible Ottobock components

- **Up to 68 kg:** 4R98, 4R120, 4R121=30
- **From 69 kg:** 2R2, 2R45=S, 4R21*, 4R23, 4R41, 4R43, 4R52*, 4R68, 4R85, 4R86, 4R89, 4R95, 4R100, 4R104=*, 4R120, 4R121=30, 4R121=34, 4R147*, 5R6=1, 5R6=2, 6A53, 6A54

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 2 (restricted outdoor walker) and mobility grade 3 (unrestricted outdoor walker).

The table that follows shows the suitable spring stiffness of the prosthetic foot, matching the patient's body weight and activity.

Spring stiffness relative to body weight and activity level			
Body weight [kg]	Low	Normal	High
44 to 52	1	1	2
53 to 59	1	2	3
60 to 68	2	3	4
69 to 77	3	4	5
78 to 88	4	5	6
89 to 100	5	6	–
101 to 116	6	–	–

2.3 Environmental conditions

Storage and transport
Temperature range $-20\text{ }^{\circ}\text{C}$ to $+60\text{ }^{\circ}\text{C}$ ($-4\text{ }^{\circ}\text{F}$ to $+140\text{ }^{\circ}\text{F}$), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

Allowable environmental conditions
Temperature range: -10 °C to +45 °C (14 °F to 113 °F)
Moisture: relative humidity: 20% to 90%, non-condensing
Chemicals/liquids: fresh water as dripping water, occasional contact with salty air (e.g. near the ocean)
Solids: dust
Prohibited environmental conditions
Chemicals/moisture: salt water, perspiration, urine, acids, soapsuds, chlorine water
Solids: dust in high concentrations (e.g. construction site), sand, highly hygroscopic particles (e.g. talcum)

2.4 Lifetime

Prosthetic foot

Depending on the patient's activity level, the maximum lifetime of the product is 3 years.

Footshell, protective sock

The product is a wear part, which means it is subject to normal wear and tear.

3 Safety

3.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION! Risk of injury and risk of product damage

- ▶ Note the combination possibilities/comboination exclusions in the instructions for use of the products.
- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 16).
- ▶ To avoid the risk of injury and product damage, do not use the product beyond the tested lifetime.
- ▶ To avoid the risk of injury and product damage, only use the product for a single patient.
- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.

- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

NOTICE! Danger of product damage and restricted functionality

- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover behaviour) or delamination of the spring are indications of loss of functionality. Unusual noises can indicate a loss of functionality.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Prosthetic foot	–
1	Protective sock (black)	S0-NPS-200*
1	Heel wedge set	KIT-00-1131*
1	RS2-K1 Runway HX only: Adjustment wrench	ACC-00-10400-4M
1	Shim wrench (1/8")	ACC-00-10400-00

Additional accessories/spare parts (not included in the scope of delivery)	
Designation	Reference number
Footshell (without cap)	Runway: FTC-3D-1*, Runway HX: FTC-3E-1*
Footshell with sandal toe (without cap)	Runway: FTC-2D-1*, Runway HX: FTC-2E-1*
Tool to change the footshell	ACC-00-10300-00
Protective sock (white)	SL=Spectra-Sock

5 Preparing the product for use

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

NOTICE

Grinding the prosthetic foot or footshell

Premature wear resulting from damage to the product

- ▶ Do not grind the prosthetic foot or footshell.

5.1 Alignment

INFORMATION

There is a plastic protector on the product's pyramid. It protects the connecting section from scratches during the alignment and fitting of the prosthesis.

- ▶ Remove the protector before the patient leaves the workshop/fitting area.

5.1.1 Applying/removing the footshell

INFORMATION

- ▶ Pull the protective sock over the prosthetic foot to prevent noises in the footshell.
- ▶ Always use the prosthetic foot with the footshell.

- ▶ Apply or remove the footshell as described in the footshell instructions for use.

5.1.2 Checking the heel height setting

INFORMATION

Make sure that the release button does not touch the footshell at any heel height.

- ▶ **If the release button touches the footshell:** Sand the footshell in the relevant spot until there is no longer any contact.

The entire range of motion for the heel height adjustment is tested before the heel height is set for the alignment.

- ▶ Set the heel height (see page 22).

5.1.3 Bench Alignment

TT bench alignment

Bench alignment process	
Required materials: 662M4 goniometer, 743S12 heel height measuring device, 743A80 50:50 gauge, alignment tool (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PRO.S.A. Assembly)	
Assemble and align the prosthetic components in the alignment apparatus according to the following specifications:	
1	Dorsiflexion/plantar flexion: Adjust the prosthetic foot to the heel height of the shoe.
2	Adduction/abduction: Adjust the angle of the prosthetic socket in the frontal plane.
3	Flexion/extension: Adjust the angle of the prosthetic socket in the sagittal plane.
4	Linear shifting of the prosthetic socket: Shift the prosthetic socket so that the plumb line runs along the front edge of the adapter on the prosthetic foot (see fig. 3).

TF bench alignment

- ▶ Observe the information in the prosthetic knee joint instructions for use.

5.1.4 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.

5.1.5 Dynamic Trial Fitting

- Adapt the alignment of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) to ensure an optimum gait pattern.
- **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal plane is achieved when the leg begins to bear weight after the heel strike. Avoid medial movement of the knee joint. If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction. If the medial movement occurs in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.
- Remove the plastic protector from the pyramid after completing the dynamic fitting and the walking exercises.

Heel too soft	
Symptoms	Possible solutions

Heel too soft	
<ul style="list-style-type: none"> • Foot makes full ground contact too quickly • Forefoot feels too stiff • Knee goes into hyperextension 	<ul style="list-style-type: none"> • Shift prosthetic socket forward relative to the foot • Use a heel wedge

Heel too hard	
Symptoms	Possible solutions
<ul style="list-style-type: none"> • Fast knee flexion, low stability • Transition from heel strike to toe-off too fast • Energy return feels low 	<ul style="list-style-type: none"> • Shift prosthetic socket back relative to the foot • Reduce the heel stiffness (shift or remove heel wedge)

Prosthetic foot too stiff	
Symptoms	Possible solutions
<ul style="list-style-type: none"> • Minimal rollover movement of the prosthetic foot at low walking speed (long full ground contact) 	<ul style="list-style-type: none"> • Select prosthetic foot in lower stiffness

Prosthetic foot too soft	
Symptoms	Possible solutions
<ul style="list-style-type: none"> • Clicking noise at the start of ground contact. • Very pronounced deformation of the forefoot at high activity levels 	<ul style="list-style-type: none"> • Select prosthetic foot in higher stiffness

5.1.5.1 Optimising the heel characteristics

The behaviour of the prosthetic foot at heel strike and during heel contact in the mid-stance phase can be adapted with the placement of a heel wedge. The heel wedge is attached with adhesive tape for testing. For definitive assembly, it is glued to the prosthetic foot.

Fitting

- 1) Use the adhesive tape on the underside of the heel wedge.
- 2) Place the heel wedge on the base spring in the recommended position (Dummybild_klein2).
- 3) Set the desired stiffness by shifting the heel wedge (anterior = harder, posterior = softer).

Definitive assembly

- 1) Remove the adhesive tape from the heel wedge. This can be done using a cloth with acetone.
- 2) Roughen the contact surface on the prosthetic foot slightly using sandpaper. Remove the dust after sanding.
- 3) Glue the heel wedge to the prosthetic foot with cyanoacrylate adhesive.

6 Use

6.1 Setting the heel height

The heel height may only be set with no load on the prosthetic foot.

CAUTION! Sit down to set the heel height. Otherwise, there is a risk of falling.

RS2-00 Runway

- 1) Push and hold the release button on the side of the prosthetic foot.
- 2) Shift the prosthetic foot in the ankle module to set the heel height.
 - Shifting the prosthetic foot forward in relation to the prosthetic socket = increases the heel height
 - Shifting the prosthetic foot backward in relation to the prosthetic socket = decreases the heel height
- 3) Let go of the release button.
- 4) Push the tip of the prosthetic foot up and down to verify that the ankle module is locked. If the release button sticks, this may release it.

RS2-K1 Runway HX

- 1) Insert an Allen key (**4 mm**) on the rear of the ankle module and turn it.
INFORMATION: 14 turns change the heel height by 2.5 cm
 - Turning clockwise = increases the heel height
 - Turning anticlockwise = decreases the heel height
- 2) Remove the Allen key.

7 Cleaning

- > **Allowable cleaning agent:** pH neutral soap (e.g. 453H10 Derma Clean)
- 1) **NOTICE! To avoid product damage, only use the allowable cleaning agents.**
Clean the product with clear fresh water and a pH neutral soap.
 - 2) **If present:** Remove dirt from water drainage contours using a toothpick and rinse.
 - 3) Rinse the soap away with clear fresh water. In doing so, rinse the foot-shell until all dirt has been removed.
 - 4) Dry the product with a soft cloth.
 - 5) Allow to air dry in order to remove residual moisture.

8 Maintenance

- A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.

- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Inspect the product for visible damage every six months. Make additional inspection appointments as needed (e.g. for highly active or heavy users).

8.1 Maintaining the ankle joint

Cleaning and greasing the mechanism

The mechanism of the ankle joint may become stiff, especially if the prosthetic foot is used in damp or dirty conditions.

- > **Required materials:** Lithium grease
 - 1) Remove the footshell from the prosthetic foot.
 - 2) Rinse the prosthetic foot thoroughly while moving the ankle joint through its full range of motion. Clean the footshell.
 - 3) Remove the old grease from the ankle joint mechanism.
 - 4) **O&P professionals only:** In case of corrosion, clean the mechanism with compressed air and acetone or WD40.
 - 5) Lubricate all sliding surfaces with lithium grease and move the mechanism through its full range of motion.
 - 6) Wipe off excess grease.

O&P professionals only: Shim adjustment

If the ankle joint mechanism makes a clicking noise while walking (load on the forefoot or heel), the shims have to be adjusted. To do so, first tighten and then loosen the screws. Setting the heel height is difficult if the screws are too tight. A clicking noise is produced if the screws are too loose.

- > **Prerequisite:** The ankle joint mechanism has been cleaned and greased.
- > **Prerequisite:** The footshell has been removed.
- > **Required materials:** Shim wrench (1/8")
 - 1) Tighten the screws on the underside of the sliding carriage until there is firm resistance.
 - 2) Gradually loosen the screws until:
 - No clicking noise occurs.
 - The heel height can be adjusted without significant resistance.

9 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

11 Technical data

RS2-00 Runway, RS2-K1 Runway HX							
Sizes [cm]	22	23	24	25	26	27	28
Heel height [mm]	0 to 50						
System height [mm]	80			82			
Build height [mm]	98			100			
Average weight with footshell [g]	690			792			
Max. body weight [kg]	116						
Mobility grade	2, 3						

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