


17AT1=1, 17AT2=*^{*}

DE Gebrauchsanweisung	4
EN Instructions for use	10
FR Instructions d'utilisation	16
IT Istruzioni per l'uso	22
ES Instrucciones de uso	29
PT Manual de utilização	35
NL Gebruiksaanwijzing	42
SV Bruksanvisning	48
DA Brugsanvisning	54
NO Bruksanvisning	60
FI Käyttöohje	66
PL Instrukcja użytkowania	72
HU Használati utasítás	78
CS Návod k použití	85
HR Upute za uporabu	91
SK Návod na používanie	96
TR Kullanma talimatı	102
RU Руководство по применению	109
JA 取扱説明書	115
ZH 使用说明书	121

Material	Strukturteil/structural part: Carbonfaserverbundwerkstoff/carbon fibre composite; Wadengurt und Polster/calf strap and padding: Polyurethan/polyurethane
	Latex-frei: Dieses Produkt enthält kein Naturgummi./ Latex-free: This product does not contain natural rubber.

INFORMATION

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INFORMATION

Date of last update: 2024-04-17

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the use of the 17AT1=1 and 17AT2=* test orthosis.

2 Product description

2.1 Available sizes

The 17AT1=1 test orthosis is available in a set with 4 17AT2=* orthoses and 2 additional 29F209=* foot components. Orthoses are chosen based on shoe size and weight.

17AT1=1 test orthosis set consisting of:

Article number	Shoe size	Allowable body weight
17AT2=L/R-14-1	36–41	up to max. 50 kg
29F209=L/R-20-2	36–41	up to max. 85 kg
17AT2=L/R-20-1	42–45	up to max. 110 kg

17AT2=* test orthoses can be delivered separately:

Article number	Shoe size	Allowable body weight
17AT2=L/R-14-1	36–41	up to max. 50 kg
17AT2=L/R-20-2	36–41	up to max. 85 kg
17AT2=L/R-20-1	42–45	up to max. 110 kg

2.2 Components/design

Scope of delivery (see fig. 1)			
Item	Quantity [Piece(s)] 17AT1	Quantity [Piece(s)] 17AT2=* 17AT2=*	Designation
1	4	1	Lower leg shell
2	4	1	17AD300=* orthotic ankle joint
3	6	1	Foot component
4	4	1	Calf strap with pad
5	8	2	Y-hook-and-loop
not illustrated	2	-	Axial washers set
not illustrated	4	2	Screws

Not included in scope of delivery			
Item	Quantity [Piece(s)]	Designation	Article number
not illustrated	1	Fit-Kit for 17AD300	17AT3=1

17AT2=* spare parts

Not included in scope of delivery (see fig. 1)			
Item	Quantity [Piece(s)]	Designation	Article number
3	1	Foot component	29F209=L/R-14-1 29F209=L/R-20-1 29F209=L/R-20-2
4	1	Calf strap with pad and Y-hook-and-loop	29U91=L/R-14 29U91=L/R-20
5	1	Y-hook-and-loop	29U25=1 (50 mm belt width) 29U25=2 (40 mm belt width)

3 Intended use

3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

3.2 Area of application

- Maximum permissible body weight up to **110 kg**
- The test orthosis is to be used **only** indoors and under the constant supervision of qualified personnel.
- Use the test orthosis unilaterally on one foot or in pairs on two feet.

3.3 Indications

- Partial or total paralysis of the leg muscles.
- Orthopaedic diseases of the lower limbs

Indications must be determined by the physician.

3.4 Lifetime

The product is designed for a lifetime of **5 years** when used as intended and assembled professionally.

3.5 Qualification

Users may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

4 Safety

4.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

4.2 Safety instructions

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

Improper moulding or application

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.

- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

5 Preparing the product for use

⚠ CAUTION

Incorrect alignment, assembly or adjustment

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

Preparing the orthotic joint

Use the 17AT1=1/17AT2=* test orthosis with the 17AT3=1 Fit Kit for the 17AD300=* orthotic ankle joint. Process the modules of the 17AD300=* orthotic ankle joint according to the 647G1788=* instructions for use.

5.1 Size Selection

- ▶ Select the orthosis size based on the shoe size and body weight (see page 10).

Replacing the foot component for size 20

The 29F209=L/R-20-1 foot component is pre-installed on delivery of size 20. The lower leg shell of size 20 can be optionally combined with the 29F209=L/R-20-2 foot component.

- 1) Loosen the joint screw and remove the joint nut from the orthotic joint.
- 2) Remove the foot component.
- 3) Install the other foot component as described in the **647G1788=*** instructions for use.

5.2 Adaptation

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Damage due to improper handling (heat treatment, drilling holes)

Loss of orthosis functionality

- ▶ Do not carry out any heat treatment on the orthosis. The orthosis is not thermoformable.

- ▶ Do not drill any holes in the bottom of the sole or the connecting elements as holes break the fibres and weaken the component.

NOTICE

Use without suitable footwear

Restriction of functionality due to insufficient stabilisation

- ▶ Only use the orthosis with closed footwear, taking the allowable effective heel height into account.

Selecting footwear

- ▶ Choose a sturdy shoe with laces or hook-and-loop fasteners and a solid toe cap to ensure the optimal effectiveness of the orthosis. The effective heel height should be between 0.5 cm and 1.5 cm.

Attaching the pad to the lower leg shell

- 1) Open the straps on the pad.
- 2) Open the triangular hook-and-loop closure on the pad (see fig. 2).
- 3) Slide the pad onto the lower leg shell starting from the inside (see fig. 3, see fig. 4).
- 4) Slide on the pad on the lower leg shell and fasten the triangular hook-and-loop closure on the pad (see fig. 5, see fig. 6).

Adapting the orthosis

The test orthosis can initially be adapted without assembling the modules of the 17AD300=* orthotic ankle joint (see fig. 7). For static and dynamic fitting of the test orthosis, install the modules in the orthotic joint, use the **17AT3=1** Fit Kit for this purpose and observe the 647G1788=* instructions for use. **INFORMATION: In order to assemble the modules, the lower leg component can be removed from the orthotic ankle joint by disassembling the lamination bar.**

- 1) If the user wears footwear with removable insoles, take the insole out of the shoe on the fitting side and store it.
- 2) Sit the user on a chair with shoe(s). Position the foot at an angle of **90°** to the lower leg. **INFORMATION: Put on the orthosis in such a way that the ankle bone has sufficient freedom of movement and does not come in contact with the orthotic ankle joint. Note the alignment of the orthotic joint relative to the compromise pivot point of the upper ankle joint.**
- 3) Mark the position of the ankle joint on the shoe with a small piece of adhesive tape (see fig. 8).
- 4) **Optionally** fabricate an insole that matches the removable insole and temporarily attach it to the foot component of the orthosis with double-sided adhesive tape (see fig. 9).
- 5) **Optionally**, foot malpositions can be corrected with an insole or a custom moulded device.
- 6) Put on the orthosis together with the shoe and check the fit of the orthosis in the shoe (see fig. 10).

Checking the fit of the orthosis

- ▶ Check the fit of the test orthosis prior to fitting on the patient.

5.3 Dynamic fitting

- > **Prerequisite:** The alignment and static fitting of the test orthosis have been checked and carried out in accordance with the 647G1788=* instructions for use.
- > The test orthosis has been applied to the user.
 - 1) Have the user stand and walk on level surfaces, on slopes and stairs.
 - 2) Adjust the modules of the orthotic ankle joint to the desired result.

Dynamic fitting process	
Standing	The orthosis generates resistance against dorsal extension in the ankle joint and extends the knee.
Walking on level surfaces	The orthosis generates resistance during rollover in the mid to terminal stance phase. The orthosis lifts the tip of the foot in the swing phase.
Walking on slopes and stairs	The orthosis controls the plantar strike downhill and dorsal extension uphill.

5.4 Removal

- 1) Open the calf straps.
- 2) Remove the shoe together with the orthosis.
- 3) Remove the orthosis from the shoe and re-insert the insole of the shoe.

6 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Textile component

▶ **NOTICE! Detergent residues may cause skin irritation and wear of material.**

- 1) Remove the textile component from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) **Optional:** Use a mesh laundry bag.
- 4) Wash the textile component in warm water at **30 °C** with a standard mild detergent. Rinse out detergent residues thoroughly.
- 5) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

Orthosis frame

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

7 Disinfection

⚠ CAUTION

Reuse on other persons and improper cleaning, disinfection

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may only be used on another person after thorough cleaning and disinfection.
- ▶ Clean the product regularly.

- 1) Prior to use on another user, clean the product with a conventional disinfectant.
- 2) Wipe all parts of the product with a disinfectant.

8 Maintenance

INFORMATION

The product may be exposed to increased loads by the patient.

- ▶ Shorten the maintenance intervals according to the expected loads.

Spare parts are listed in the section "Components/design".

9 Disposal

Dispose of the product in accordance with national regulations.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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