## Genium X4 3B5-4=P, 3B5-4=ST

<b>DE</b> Gebrauchsanweisung [kurz] (Fachpersonal) .....	3
<b>EN</b> Instructions for use [short] (qualified personnel) .....	6

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**INFORMATION**

Date of last update: 2024-05-15

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report any serious incident related to the product, in particular any deterioration of health, to the manufacturer and the competent authority in your country.

**Note the instructions for use**

This document contains the sections "Intended use" and "Safety", and does not replace the instructions for use. The instructions for use are available in digital form (PDF format) (see page 2 of this document). To view PDF files, you need an Internet browser or other program capable of opening this file format.

You may request a printed copy of the instructions for use at the e-mail address on page 2. This will be delivered to the specified address within 7 days after the order is received.

The product "Genium X4 3B5-4=" is called the product/prosthesis/knee joint/component in the following.

The product "USB charging adapter 757L47=1" is referred to as the charging adapter below.

Basic UDI-DI: 4064411000000003B5-4HJ

**2 Intended use****2.1 Indications for use**

The product is to be used **solely** for lower limb exoprosthetic fittings.


**2.2 Conditions of use**

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.).

Permissible ambient conditions are described in the technical data (Technical data).


The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

The MOBIS classification describes the mobility grade and body weight, and makes it easy to identify compatible components.

**Knee joint with attached 2R68=280 Axon tube adapter**

The product is recommended for mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands). Approved for a body weight of up to **150 kg (330 lbs)**.

## Knee joint with attached 2R69=280 Axon tube adapter with torsion

The product is recommended for mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands). Approved for a body weight of up to **125 kg (275 lbs)**.

### 2.3 Indications

- For patients with knee disarticulation, transfemoral amputation or hip disarticulation
- For unilateral or bilateral amputation
- Those affected by dysmelia, with the characteristics of the affected body part corresponding to a knee disarticulation, transfemoral amputation or hip disarticulation
- Osseointegration
- The patient must fulfil the physical and mental requirements for perceiving visual/acoustic signals and/or mechanical vibrations.
- The patient must be able to understand usage and safety messages and put them into practice.

### 2.4 Contraindications

#### 2.4.1 Absolute Contraindications

- Body weight over 150 kg

#### 2.4.2 Relative Contraindications

- Body weight less than 35 kg

### 2.5 Qualification

The product may be fitted only by qualified personnel authorised by Ottobock after completing the corresponding training.

If the product is to be connected to an osseointegrated implant system, the qualified personnel must also be authorised for the connection to the osseointegrated implant system.

### 3 Safety

Ottobock developed this product according to applicable standards and rules and tested it multiple times. In order for you to enjoy the product, we need your help. Only by complying with the following instructions can we guarantee safe treatment and operation.

#### 3.1 Meaning of warning levels

**WARNING!** Failure to follow the instructions can lead to serious accidents and injuries.

**CAUTION!** Failure to follow the instructions can lead to accidents and injuries.

#### 3.2 Before treatment

**WARNING! Possibility of influences on the human body**

- ▶ Observe the areas of application and operating conditions of your product in combination with osseointegrated implant systems according to the information of the manufacturer.
- ▶ Note the instructions of the clinical personnel that indicated the use of the osseointegrated implant system.

#### 3.3 During treatment

**CAUTION! Possibility of falling**

- ▶ If the product is to be used in water, salt water or chlorinated water, the entire prosthesis has to be suitable for such use. Resistance to these liquids has to be checked for each prosthetic component.
- ▶ The listed prosthetic feet may only be combined with the respective foot sizes [cm] described, depending on the patient's body weight. Please contact Ottobock customer service if you would like a combination outside the approved ranges.
- ▶ Observe the alignment and assembly instructions.
- ▶ At maximum flexion under full load, it is essential to maintain a minimum distance of 3 mm (5 mm when using a functional cosmesis) between the frame of the knee joint and the socket.
- ▶ If the minimum distance of 3 mm (5 mm when using a functional cosmesis) is not met, a flexion stop has to be installed. If a flexion stop has already been installed, it has to be replaced with the next larger flexion stop. If the largest flexion stop is already installed, please contact Ottobock Customer Service.
- ▶ There should be no contact between the knee joint and the socket at maximum extension (obtained through manual extension).
- ▶ If there is contact, either install an additional, appropriate adapter between the knee joint and the socket, or change the socket position by optimising alignment.

A series of 21 horizontal lines spanning the width of the page, providing space for writing.

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

Otto Bock Healthcare Products GmbH  
Brehmstraße 16 · 1110 Wien · Austria  
T +43-1 523 37 86 · F +43-1 523 22 64  
info.austria@ottobock.com · www.ottobock.com