
1C70 Evanto

DE Gebrauchsanweisung	6
EN Instructions for use	12
FR Instructions d'utilisation	17
IT Istruzioni per l'uso	24
ES Instrucciones de uso	30
PT Manual de utilização	36
NL Gebruiksaanwijzing	42
SV Bruksanvisning	48
DA Brugsanvisning	53
NO Bruksanvisning	59
FI Käyttöohje	65
PL Instrukcja użytkowania	70
HU Használati utasítás	77
CS Návod k použití	82
RO Instrucțiuni de utilizare	88
HR Upute za uporabu	94
SL Navodila za uporabo	100
SK Návod na použitie	106
BG Инструкция за употреба	112
TR Kullanma talimatı	118
EL Οδηγίες χρήσης	124
RU Руководство по применению	130
JA 取扱説明書	136
ZH 使用说明书	142
KO 사용 설명서	147

1

Bench Alignment, transtibial

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INFORMATION

Date of last update: 2024-05-28

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The 1C70 Evanto prosthetic foot has two exchangeable heel wedges. The heel wedges influence the behaviour of the prosthetic foot in different phases of the gait cycle. This allows the prosthetic foot to be individually adapted.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

The MOBIS classification describes the mobility grade and body weight, and makes it easy to identify compatible components.

The product is recommended for mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands).

- **The prosthetic foot is not suitable for sports prostheses.**
Certain movement patterns in sports (e.g. soccer, handball or tennis) can damage the prosthetic foot: pronounced stopping when running and spontaneous changes in direction

Body weight [kg]	Spring stiffness
Up to 58	1
59 to 72	2
73 to 86	3
87 to 100	4
101 to 125	5

TIP: Also try the next higher stiffness so the user can experience the difference.

2.3 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Limited combination options for Ottobock components

Combination with prosthetic knee joints	
From stiffness 4 and foot size 28	Do not combine with 3R85 Dynion
From stiffness 5 and foot size 26	Do not combine with 3R55

The prosthetic foot generates high moments in the ankle area. Use structural components with higher weight limits:

Body weight [kg]	73 to 86	87 to 100		101 to 125
Foot size [cm]	From 28	From 26	From 28	From 24
Structural component weight limit [kg]	≥125	≥125	≥150	≥150

2.4 Environmental conditions

Storage and transport
Storage temperature: -20 °C to +60 °C, relative humidity: 20 % to 90 %, no mechanical vibrations or impacts

Allowable environmental conditions
Temperature range: -10 °C to +45 °C
Chemicals/liquids: fresh water, soap solution, chlorine water, salt water ≤3.5 % salt content
Moisture: Submersion: max. 1 h in 2 m depth, relative humidity: no restrictions
Solids: Dust, occasional contact with sand
Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage (see page 16).

Prohibited environmental conditions
Solids: highly liquid-binding particles (e.g. talcum), dust in high concentrations (e.g. construction site), intensive contact with sand
Chemicals/liquids: perspiration, urine, acids, continuous use in liquid media

2.5 Lifetime

Depending on the patient's activity level, the maximum lifetime of the product is 3 years.

Footshell, protective sock

The products are wear and tear parts that are subject to normal degeneration.

3 General safety instructions

CAUTION! Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ Note the combination possibilities/combination exclusions in the instructions for use of the products.
- ▶ Observe the maximum lifetime of the product.
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not expose the product to excessive strain (section: "Area of application" see page 12) or to prohibited environmental conditions (see page 13).
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ The product may only be used by one person; reuse on another person is prohibited.

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover behaviour) or delamination of the spring are indications of loss of functionality. Unusual noises can indicate a loss of functionality.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	-
1	Quick start guide	647G1809=*
1	Prosthetic foot	-
1	Protective sock	Size 22 to 25: 2U3=1-7 Size 26 to 30: 2U3=2-7
Additional heel wedges		
1	Upper heel wedge	2F70=*
1	Lower heel wedge	2F71=*

5 Preparing the product for use

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

NOTICE

Altering the prosthetic foot or footshell

Premature wear resulting from damage to the product

- ▶ Do not alter the prosthetic foot or footshell.

5.1 Applying/removing the footshell

INFORMATION

- ▶ Pull the protective sock over the prosthetic foot to prevent noises in the footshell.

▶ **CAUTION! Always use the prosthetic foot with the footshell.**

Apply or remove the footshell as described in the footshell instructions for use.

5.2 Bench Alignment

Bench alignment, transtibial prosthesis

> **Required materials:** 662M4 goniometer, 743S12 heel height measuring device, 743A80 50:50 gauge, alignment apparatus (e.g. 743A200 PROS.A. Assembly)

- 1) Position the prosthetic foot in the alignment apparatus.
 - **Heel height: effective heel height of everyday shoe + 5 mm**
- 2) **Sagittal plane:** Align the prosthetic foot with the alignment reference line.
 - Position the alignment reference line on the markings of the footshell.
- 3) **Frontal plane:** Align the prosthetic foot with the alignment reference line.
 - Position the alignment reference line on the markings of the footshell (**external rotation approx. 5°**).
- 4) Install the structural components between the prosthetic foot and prosthetic socket.
- 5) **Sagittal plane:** Align the prosthetic socket with the alignment reference line.
 - Determine the centre of the prosthetic socket proximally and distally with the 50:50 gauge and mark the centre line.
 - Mark the socket reference point on the centre line: **For transtibial prostheses, at the height of the centre of the patella.** The alignment reference line runs vertically through the socket reference point.

→ Adjust the socket flexion by rotating around the socket reference point: **individual residual limb flexion + 5°**

- 6) **Frontal plane:** Align the prosthetic socket with the alignment reference line.
 → Take the abduction or adduction of the residual limb into account.

Bench alignment, thigh/knee disarticulation prosthesis

- ▶ Observe the instructions for use for the prosthetic knee joint.

5.3 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.
- **Check the length of the prosthesis:** Compression of the prosthetic foot and other components (e.g. a vacuum pump) may require correction of the length.

The alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) can be requested from Ottobock as needed.

5.4 Dynamic Trial Fitting

The optimal gait pattern is developed during dynamic fitting. The prosthesis is aligned optimally in the frontal plane and the sagittal plane for this purpose.

Walking and standing with Evanto looks and feels different. The user typically immediately feels increased flexibility when walking and standing.

Recommended action

- Let the user walk for a few minutes to get used to the new feel and flexibility of the Evanto. In doing so, observe the process of heel strike, rollover and toe-off.
- Have the user walk up ramps. Compared to the previously used prosthetic foot, this should be perceived as equally heavy or lighter.
- ▶ **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal planes is achieved after the leg begins to bear weight following the heel strike. Avoid medial movement of the knee joint.
 - If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction.
 - If the knee joint moves in the medial direction in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.

5.4.1 Adjusting stiffness

The stiffness of the prosthetic foot can be adjusted by replacing the two heel wedges:

- **Upper heel wedge:** Determines the heel stiffness.
 - **Lower heel wedge:** Determines the forefoot stiffness and heel stiffness.
- The heel stiffness influences the knee flexion effect (soft heel wedge = less knee flexion effect. Hard heel wedge = more knee flexion effect). Forefoot stiffness influences the rollover of the prosthetic foot.

The table shows some examples:

Situation	Measures (several measures individually or in combination)
The user feels a stepped rollover (second half of the stance phase).	First insert a harder lower heel wedge. To make the heel softer again, insert a softer upper heel wedge.
The user feels that the resistance is too high when walking up a slope.	Move the prosthetic foot gradually in the posterior direction.
The user feels a double rocking during load transfer.	Insert a harder upper heel wedge.

Situation	Measures (several measures individually or in combination)
The user feels a dead point during the rollover (first half of the stance phase).	Insert a harder upper heel wedge.
Some users (mobility grade 3 and especially mobility grade 2) prefer a foot position slightly further anterior. This reduces the knee flexing effect, gives a secure feeling at heel strike and a stable rollover.	

Replacing heel wedges

Heel wedge degrees of hardness: The heel wedges are marked with numbers. The smaller the number, the softer the heel wedge.

To replace the heel wedges, the strap looped around the heel must be relaxed. The seat of the strap is then checked and the strap is re-tightened.

- 1) Lift off the cover of the prosthetic foot towards the heel (see fig. 3).
- 2) Fully unscrew the strap tensioner (see fig. 5).
- 3) **NOTICE! Make sure the strap stays in position. This prevents damage to the base spring of the prosthetic foot.**
Remove the lower heel wedge (see fig. 6).
- 4) **Optional:** Replace the upper heel wedge, aligning the recess in the heel wedge upwards (see fig. 7, see fig. 8).

INFORMATION: As soon as the upper heel wedge has been removed, the seat of the strap tensioner can be easily checked from the bottom.

- 5) Insert the lower heel wedge (see fig. 9).
- 6) Check whether the strap lies correctly in the guide:
 - The strap is positioned all around the guide (see fig. 10).
 - The Ottobock logo is centred on the side of the foot (see fig. 11).
 - The strap tensioner is not at an angle (see fig. 12).
- 7) Tighten the strap tensioner (tightening torque: **4 Nm**) (see fig. 13).
- 8) Put the cover of the prosthetic foot back on (see fig. 14).

6 Cleaning

> **Allowable cleaning agent:** pH-neutral soap (e.g. 453H10 Derma Clean)


- 1) **NOTICE! To avoid product damage, only use the allowable cleaning agents.**
Clean the prosthetic foot and footshell with clear fresh water and pH-neutral soap.
- 2) **If present:** Use a toothpick to remove dirt from the contours that drain off the adapters or footshell and rinse.
- 3) Rinse the soap away with clear fresh water. In doing so, rinse the footshell until all dirt has been removed.
- 4) **In case of heavy soiling and noises:** Pull the footshell and protective sock off the prosthetic foot. Clean all products individually with water or compressed air.
- 5) Dry the product with a soft cloth.
- 6) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ During regular inspection: Inspect the prosthesis for signs of wear and check its function.
- ▶ Direct the patient to observe the following maintenance instructions:
Check the entire prosthesis for damage after each use.
In case of changes in or loss of functionality and damage, do not continue using the prosthesis and have it checked by qualified personnel.

7.1 Strap tensioner – replacing the screw

The strap tensioner consists of a thrust piece and a screw. The screw presses the thrust piece against the strap and thus generates the tension. To compensate for fabrication tolerances of the strap, the screw is available in 6 different lengths. During the fabrication of the prosthetic foot, the appropriate length is selected.

Point indentations on the screw head indicate the length of the screw. They are arranged around the Allen head. The more indentations (1 to 6), the longer the screw.

- ▶ Replace a defective screw with a new screw of the same length.
- ▶ If you are not sure of the correct length of the screw: Contact Ottobock.

8 Disposal

Do not dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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