## 5A60=\* Varos

<b>DE</b> Gebrauchsanweisung (Fachpersonal) .....	3
<b>EN</b> Instructions for use (qualified personnel) .....	31
<b>FR</b> Instructions d'utilisation (Personnel spécialisé) .....	59
<b>IT</b> Istruzioni per l'uso (Personale tecnico specializzato) .....	87
<b>ES</b> Instrucciones de uso (Personal técnico especializado) .....	115
<b>PT</b> Manual de utilização (Pessoal técnico) .....	143
<b>NL</b> Gebruiksaanwijzing (Vakmensen) .....	171
<b>SV</b> Bruksanvisning (Fackpersonal) .....	197
<b>DA</b> Brugsanvisning (Faguddannet personale) .....	223
<b>NO</b> Bruksanvisning (Fagpersonell) .....	249
<b>FI</b> Käyttöohje (Ammattihenkilöstö) .....	275
<b>PL</b> Instrukcja użytkowania (Personel fachowy) .....	301
<b>HU</b> Használati utasítás (szakszemélyzet) .....	329
<b>BG</b> Инструкция за употреба (Специалисти) .....	355
<b>JA</b> 取扱説明書 (有資格担当者) .....	383

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<b>1</b>	<b>Foreword .....</b>	<b>33</b>
<b>2</b>	<b>Product description .....</b>	<b>33</b>
2.1	Design .....	33
2.2	Function .....	33
2.3	Combination possibilities.....	33
<b>3</b>	<b>Intended use .....</b>	<b>34</b>
3.1	Indications for use .....	34
3.2	Area of application .....	34
3.3	Conditions of use .....	34
3.4	Indications.....	35
3.5	Contraindications .....	35
3.6	Qualification .....	35
<b>4</b>	<b>Safety.....</b>	<b>35</b>
4.1	Explanation of warning symbols.....	35
4.2	Structure of the safety instructions.....	35
4.3	General safety instructions.....	35
<b>5</b>	<b>Scope of Delivery and Accessories.....</b>	<b>38</b>
<b>6</b>	<b>Preparing the product for use.....</b>	<b>39</b>
6.1	Required machines, tools and materials .....	39
6.2	Preparing for socket fabrication.....	39
6.2.1	Determining residual limb circumferences, residual limb length, flexion contracture ....	39
6.2.2	Determining the socket size and socket length .....	42
6.2.3	Selecting the lateral side rail .....	42
6.3	Fabricating the socket.....	43
6.3.1	Shortening socket components .....	43
6.3.2	Cutting off and gluing in the fabric adhesive strips.....	45
6.3.3	Gluing in rubber nubs .....	45
6.3.4	Routing the cables.....	45
6.3.5	Installing the axle.....	48
6.3.6	Screwing together the side rails and socket shells .....	48
6.3.7	Installing the lace crossover parts.....	49
6.3.8	Installing the socket brim pad.....	50
6.3.9	Screwing the bottom plate to the connection adapter .....	50
6.4	Finishing the prosthetic socket .....	51
6.5	Verification of function and fit .....	51
<b>7</b>	<b>Use .....</b>	<b>53</b>
7.1	Putting on the liner and socket .....	53
7.2	Taking off the socket and liner .....	53
<b>8</b>	<b>Cleaning .....</b>	<b>53</b>
<b>9</b>	<b>Maintenance .....</b>	<b>54</b>

**10 Legal information .....54**  
10.1 Liability .....54  
10.2 Trademarks .....54  
10.3 CE conformity .....54

**11 Technical data .....54**

**12 Appendices .....55**  
12.1 Symbols Used .....55  
12.2 Spare parts .....56  
12.3 Accessories .....57

# 1 Foreword

## INFORMATION

Date of last update: 2024-07-01

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The product "5A60=\* Varos" is referred to simply as the product/socket below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.



Instruct the patient in the proper use and care of the product. The product may not be transferred to the patient without prior instruction.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Design

The product consists of the following components:

1. Lateral socket shell
2. Upper Boa lace
3. Medial socket shell
4. Lace crossover part
5. Lateral side rail
6. Connection adapter
7. Lower Boa lace
8. Lower Boa dial
9. Upper Boa dial
10. Socket brim pad
11. Medial side rail
12. Magnetic bottom shell
13. Fabric adhesive strips

### 2.2 Function

The product is a transfemoral socket with a modular structure. The product connects the residual limb to the knee joint. The product can be adapted to the volume of the residual limb while it is being worn.

The length of the product is adjusted by qualified personnel. The circumference can be individually adjusted by the user in the course of day-to-day use.

### 2.3 Combination possibilities

This product can be combined with the following Ottobock components:

- All electronic, mechatronic and mechanical knee joints

## Adapters

- 4R54 Socket adapter with pyramid adapter
- 4R74, 4R74=AL Socket adapter with pyramid adapter
- 4R23 Socket adapter with pyramid adapter
- 4R77 Socket adapter with pyramid adapter, rotating
- 4R73=A, 4R73=D Socket adapter with pyramid adapter, eccentric
- 4R55 Socket adapter with pyramid receiver
- 4R95 Socket adapter with pyramid receiver
- 4R22 Socket adapter with pyramid receiver
- 4R51 Socket adapter with pyramid receiver, rotating
- 4R37 Socket adapter with pyramid receiver, rotating

## 3 Intended use

### 3.1 Indications for use

The product is to be used solely for the prosthetic treatment of the lower limbs. The product is to be used for the treatment of transfemoral amputations.

### 3.2 Area of application

The MOBIS classification describes the mobility grade and body weight, and makes it easy to identify compatible components.

The product is recommended for mobility grade 1 (indoor walker), mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands).

Size	S	SC	M	MC	L
Max. body weight [kg]	100	100	100	100	125

### 3.3 Conditions of use

The product is intended **exclusively** for use on **one** patient.

#### **⚠ CAUTION**

#### **Reuse on another patient and improper cleaning**

Skin irritation, development of eczema or infections due to contamination with germs.

- ▶ The product should be used only on the patient for whom it was fitted.
- ▶ If the product is to be used on another patient, it must first be cleaned with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean) and disinfected.
- ▶ If a surface disinfection process using an alcohol-based disinfectant without additional additives, for example "Terallin liquid", is used for disinfection, the application instructions and safety notices of the disinfectant manufacturer must be observed.

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports that require continuous high walking speeds (running ...) or extreme sports (free climbing, parachuting, paragliding, etc.).

#### **Reusing the product**

With additional fall protection, the product can be used multiple times for internal testing.

If the product is reused by more than one patient (check socket), proper disinfection has to be ensured.

## INFORMATION

A product that has been used as a check socket may no longer be used for the long-term treatment of the patient.

### 3.4 Indications

- For users with unilateral or bilateral transfemoral amputation
- For users with a normal weight bearing residual limb of conical or cylindrical shape
- For users with or without daily volume fluctuations of the residual limb
- For users whose residual limb length and residual limb circumference are in the specified range (see the instructions for use for the 6Y200=\* Varos liner).

### 3.5 Contraindications

- Skin changes in the fitting area (e.g. open sores, skin diseases)
- Sensory disorders on the residual limb
- Pronounced atrophic scars or very bulgy residual limb

### 3.6 Qualification

Only qualified personnel (such as an O&P professional) may prepare the product for use and fit patients with the product.

## 4 Safety

### 4.1 Explanation of warning symbols

 **WARNING** Warning regarding possible serious risks of accident or injury.

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 4.2 Structure of the safety instructions

#### **WARNING**

##### **The heading describes the source and/or the type of hazard**

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

### 4.3 General safety instructions

#### **WARNING**

##### **Non-observance of safety notices**

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

#### **WARNING**

##### **Use of the product near active implanted systems**

The product's magnetic field can interfere with active implanted systems (e.g. pacemakers, defibrillators etc.).

- ▶ When using the product in the immediate vicinity of active implanted systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

**⚠ CAUTION**

**Independent manipulation of system components**

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired or replaced by qualified personnel.

**⚠ CAUTION**

**Mechanical overload during transportation**

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Only use the transport packaging for transportation.

**⚠ CAUTION**

**Mechanical overloading of the product**

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Do not subject the product to excessively strong mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

**⚠ CAUTION**

**Signs of wear on system components**

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Annual safety inspections have to be carried out in the interest of the patient's safety and in order to maintain operational reliability.

**⚠ CAUTION**

**Reuse on another patient and improper cleaning**

Skin irritation, development of eczema or infections due to contamination with germs.

- ▶ The product should be used only on the patient for whom it was fitted.
- ▶ If the product is to be used on another patient, it must first be cleaned with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean) and disinfected.
- ▶ If a surface disinfection process using an alcohol-based disinfectant without additional additives, for example "Terallin liquid", is used for disinfection, the application instructions and safety notices of the disinfectant manufacturer must be observed.

**⚠ CAUTION**

**Risk of pinching when inserting the liner into the socket**

Injury due to powerful magnetic attraction between the magnetic cap of the matching 6Y200=\* Varos liner and the magnetic bottom shell of the product.

- ▶ When inserting the liner into the socket, ensure that no fingers/body parts are between the liner and socket.

**⚠ CAUTION**

**Storage of the product in direct sunlight**

- > Falling due to the breakage of load-bearing components caused by plastic parts of the product becoming brittle due to UV exposure.
- > Discolouration of plastic parts.
- ▶ Avoid storing the product in direct sunlight.

**⚠ CAUTION**

**Contact of the product with embers or fire**

Injuries due to contact with molten material.

- ▶ Keep the product away from embers and open flames.

**⚠ CAUTION**

**Incorrectly secured screws**

Falling due to breakage of load-bearing components caused by screw connections coming loose.

- ▶ Observe the specified torque values.
- ▶ Observe the instructions for securing the screw connections.

**⚠ CAUTION**

**Use of a damaged product**

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Regularly check the socket, Boa cable, Boa dials, socket brim pad and fabric or plastic adhesive strips for wear and tear or damage.
- ▶ If components are damaged or a Boa cable is torn, the socket has to be repaired by qualified personnel.

**NOTICE**

**Improper product care**

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).

## 5 Scope of Delivery and Accessories

### Scope of Delivery

- 1 pc. medial side rail (item 1)
- 1 pc. lateral side rail - straight (item 2)
- 1 pc. lateral side rail - curved (item 3)
- 3 pc. clamping pieces (item 4)
- 7 pc. slot nut (item 5)
- 7 pc. M4x8 countersunk head screw (item 6)
- 4 pc. Boa cable (item 7)
- 6 pc. lace crossover part 3 pc. "L", 3 pc. "R" (item 8)
- 1 pc. Boa tool (item 9)
- 2 pc. Boa reel (item 10)
- 2 pc. Boa dial (item 11)
- 1 pc. axle (item 12)
- 1 pc. connection adapter (item 13)
- 6 pc. cap screw (item 14)
- 4 pc. M6x12 countersunk screw (item 15)
- 4 pc. M6x16 countersunk screw (item 16)
- 1 pc. magnetic bottom shell (item 17)
- 14 pc. rubber nub (item 18)
- 1 pc. socket brim pad (item 19)
- 1 pc. fabric adhesive strip set (item 20)
- 1 pc. lateral socket shell, includes 2 pc. Boa housings (item 21)
- 1 pc. medial socket shell (item 22)
- 1 pc. Instructions for use (qualified personnel)
- 1 pc. instructions for use (user)

### Accessories

The following components are not included in the scope of delivery and must be supplied for operation:

#### INFORMATION

The product always has to be used with the 6Y200=\* Varos liner in the appropriate size.

- 6Y200=\* Varos liner

The following components are not included in the scope of delivery and may be ordered separately:

- Magnetic cap, threaded version, art. no. 5E10=\*  
Contains: magnetic cap and M10 mounting screw
- Set of plastic adhesive strips, art. no. 5E9=\*  
Contains: 10 pc. plastic adhesive strips

## 6 Preparing the product for use

### 6.1 Required machines, tools and materials

- 743S10 Callipers
- Retractable measuring tape (e.g.: 743B4)
- 636K13 adhesive (Loctite 241)
- 714U4 Pin punch set (ø 3 mm, 4 mm, 5 mm pin punches)
- 710D20 Torque wrench (up to 25 Nm)
- 709Y1 Torque screwdriver
- 710D4 torque wrench (up to 30 Nm)
- 710Y24 Torx bit set
- 719S7 Sizing scissors
- 756B2 Electric jigsaw (1x 230 V N/PE, 50 Hz)
- 756Y19=2 Jigsaw blades for plastic/non-ferrous metals
- 706Z6=160 High-leverage side cutting pliers
- 627B6=30 Crepe adhesive tape
- 634A58 Isopropyl alcohol
- 705B2=300 Ball-peen hammer
- 756E6 Hot air gun (1x 230 V N/PE, 50/60 Hz)
- 701P250=1 Belt sander (1x 230 V, 50/60 Hz)
- 701F39 Grinding and milling machine
- 701S3=G Band saw (3x 400 V N/PE, 50 Hz)
- 749Z12=150X25 Sanding brush
- 749F16=\* Silicone sanding cone set
- 749F6 Sanding drum
- Tape measure (e.g.: 743B2)
- 710H30 Kraffform screwdriver set

### 6.2 Preparing for socket fabrication

#### 6.2.1 Determining residual limb circumferences, residual limb length, flexion contracture



##### Determining the residual limb circumference

Complete the following steps to determine the circumference of the residual limb:

- 1) Have the patient sit on a chair.
- 2) Measure the proximal residual limb circumference as the top circular circumference (see fig. 1).
- 3) Measure the distal residual limb (see fig. 2).

**INFORMATION: The measurement is taken at the transition from the hemispherical residual limb end to the conical section of the residual limb.**

**INFORMATION: If there is flaccid soft tissue at the end of the residual limb, the residual limb end has to be formed into a hemisphere while measuring.**



### Determining the residual limb length


Complete the following steps to determine the residual limb length:

- ▶ Have the patient sit on a chair.

The upper body and lower leg must be at a 90° angle to the thigh while measuring (see fig. 3).


- 1) Have the patient spread their legs slightly (see fig. 4).


- 2) Position the callipers on the patient's pubic bone (pubis) (see fig. 5).  
**CAUTION! Apply sufficient pressure so the callipers do not slip, but not so much as to cause pain.**


- 3) Have the patient adduct their legs (see fig. 6).

- 4) If there is flaccid, hanging soft tissue, shape the end of the residual limb with your hands (see fig. 7).  
Apply only a bit of pressure to push the manual slide of the callipers (743S10) against the end of the patient's residual limb.


- 5) Align both thighs parallel with the callipers and read the measurement (see fig. 8).

**INFORMATION: Measuring with the legs abducted would lead to a deviation in the measurement so the prosthetic socket would be cut to length incorrectly.**

### Determination of flexion contracture

Complete the following steps to determine the flexion contracture:

- 1) Have the patient lie on his or her back.
- 2) Ask the patient to bend the sound leg.  
→ The lumbar spine area has to lie flat on the supporting surface.
- 3) Determine the flexion contracture (angle) (see fig. 9).

### 6.2.2 Determining the socket size and socket length

The reference value for the required socket length is:

Required socket length = residual limb length

#### INFORMATION

For this product, the socket length is defined as the distance between the bottom of the socket and the medial socket brim, measured in the direction of the longitudinal axis of the socket.

The exact socket length can be determined by trying several socket lengths, preferably in the course of using a check socket.

#### INFORMATION

Pressure on the medial socket brim can be reduced by reducing the socket length.

Pressure on the residual limb end can be reduced by increasing the socket length.

The socket size can be determined based on the measured residual limb circumferences and the residual limb length. Use the chart below for this purpose.

	<b>VAROS S</b>	<b>VAROS SC</b>	<b>VAROS M</b>	<b>VAROS MC</b>	<b>VAROS L</b>
<b>Proximal residual limb circumference</b>	410–510 mm	460–580 mm	460–580 mm	520–650 mm	520–650 mm
<b>Distal residual limb circumference</b>	320–390 mm	320–390 mm	350–440 mm	350–440 mm	400–490 mm
<b>Residual limb length</b>	200–320 mm	200–320 mm	200–320 mm	200–320 mm	200–320 mm
<b>Adjustable socket length</b>	200–320 mm	200–320 mm	200–320 mm	200–320 mm	200–320 mm

### 6.2.3 Selecting the lateral side rail

The distal residual limb circumference determines which lateral side rail is appropriate. Use the chart below to select the appropriate lateral side rail.

		<b>Straight lateral side rail</b> (see the figure “Scope of delivery”, item 2)	<b>Curved lateral side rail</b> (see the figure “Scope of delivery”, item 3)
<b>VAROS S and SC</b>	Distal residual limb circumference < 370 mm	x	
	Distal residual limb circumference > 370 mm		x
<b>VAROS M and MC</b>	Distal residual limb circumference < 410 mm	x	
	Distal residual limb circumference > 410 mm		x
<b>VAROS L</b>	Distal residual limb circumference < 460 mm	x	
	Distal residual limb circumference > 460 mm		x


## 6.3 Fabricating the socket

### 6.3.1 Shortening socket components

#### INFORMATION

The medial side rail and the medial and lateral socket shells must be shortened to adjust the socket length to less than 300 mm.

The lateral side rail must be shortened to adjust the socket length to less than 320 mm.

#### INFORMATION

The parts to be modified have corresponding cutting marks and measurements on the inside/reverse side. All parts always must be shortened to the same measurement. Parts must always be shortened along the cutting marks provided.

Perform the following steps to shorten the socket components (medial side rail, lateral side rail, medial socket shell, lateral socket shell):

- 1) Apply adhesive tape around the circumference of the socket shells at the level of the measured residual limb length.

- 2) Transfer the cutting marks from the inside of the socket shells to the outside (onto the adhesive tape) (see fig. 10).
- 3) Extend the cutting marks along the entire circumference of the parts.
- 4) Use a hand saw, jigsaw or band saw to trim the socket components (see fig. 11).

**INFORMATION: When using a jigsaw, select a low cutting speed to avoid melting the plastic at the cut edge.**

**NOTICE! When shortening the side rails with a hand saw, the clamping fixture (e.g. a vice) is not permitted to put strain on the thinnest section (row of bores) on the side rails, otherwise the side rails may break.**

- 5) Round the edges and corners with a suitable tool.
- 6) Smooth the cut edges with the silicone sanding cone (see fig. 12).

**⚠ CAUTION**

**Damage to the cable ducts and Boa housing**

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Re-shaping additional product components other than the parts identified in this section is not permitted (in particular, note the area of the cable guides and Boa housing).

The socket shells can be re-shaped with a hot air gun if needed.

The following areas can be re-shaped:

- Lateral socket shell in the area between the Boa dials (to bring the lateral socket brim into contact with the body) (see fig. 13)
- Anterior section of the socket brim (increase the outward flare) (see fig. 14)
- Medial and lateral socket shell with distal area directly next to the side rail guides (reduce socket cross-section in the distal area)
- Area of the socket brim pad (to optimise pressure distribution in this area)

**INFORMATION**

When heating the socket shells with the hot air gun, make sure to heat them evenly (from the inside and outside of the socket). Excessive heating causes glossy areas on the matt surface of the socket shells.

### 6.3.2 Cutting off and gluing in the fabric adhesive strips

Complete the following steps to attach the fabric adhesive strips:

- 1) Clean the recesses on the inside of the socket shells with isopropanol.
- 2) Position the fabric adhesive strips in the corresponding recesses and make a mark 2 mm away from the cut edge (see fig. 15).
- 3) Trim the fabric adhesive strips at the markings, rounding the edges.
- 4) Pull the fabric adhesive strips off the backing and glue them into the recesses (see fig. 16).

**INFORMATION: Note the alignment of the fabric adhesive strips: The hairs on the fabric adhesive strips point down, which makes donning easier and improves the hold while wearing.**

### 6.3.3 Gluing in rubber nubs

Rubber nubs are glued in to reduce sliding movements between the medial and lateral socket shells that result in wear and tear.

Complete the following steps to glue in the rubber nubs:

- 1) Clean the recesses on the inside of the medial socket shell with isopropanol.
- 2) Pull the rubber nubs off the backing material and glue them into the recesses (see fig. 17).

**INFORMATION: In particular, the rubber nubs should be glued in where a cable is expected to come out.**

### 6.3.4 Routing the cables

#### Selecting the cable length


Four different lace lengths are included with the product. Two of them are selected according to the respective socket length, as follows:

Socket length [in mm]	200/210	220/230	240/250	260/270	280/290	300/310	320
Proximal lace length [in mm]	900	900	900	900	900	1000	1000
Distal lace length [in mm]	1000	1000	1300	1300	1600	1600	1600
Number of lace crossovers in the distal lace	1	1	2	2	3	3	3

### INFORMATION

To keep the load on the cable as low as possible, thereby increasing the lifetime, make sure that the distal cable has the maximum number of cable crossings according to the table.





### Threading the Boa cables

Complete the following steps to thread the Boa cables:

- 1) Thread the distal Boa cable into the cable guides on the medial socket shell (see fig. 18, orange cable guide), producing the number of cable crossings specified in the table.
- 2) Thread the proximal Boa cable into the medial socket shell, initially on one side.
- 3) Insert the lateral socket shell into the medial socket shell.
- 4) Thread the proximal Boa cable into the cable guide above the dial housing (see fig. 18, red cable guide).
- 5) Thread the proximal Boa cable into the medial socket shell on the second side.
- 6) Thread the two ends of the two Boa cables through the corresponding openings into the dial housing.

## Installing the Boa reels

Complete the following steps to install the Boa reel:

- 1) Thread one end of the proximal Boa cable through the bore marked in red (see fig. 19).
- 2) Fold over the cable end and thread it back through a parallel bore (see fig. 20).
- 3) Fold the cable end over again and insert it to the stop into a blind hole next to the bore marked in red (see fig. 21).

**INFORMATION: The Boa cable is secured against slipping out by means of self-locking due to the two sharp turns.**

- 4) Tighten the two resulting cable loops (see fig. 22, see fig. 23).
- 5) Perform the described steps for the second cable end on the opposite bore marked in red.
- 6) Carefully pull on the Boa cable on both sides of the Boa housing.
- 7) Lightly push the reel into the Boa housing with your finger (see fig. 24).

**INFORMATION: The spur gearing on the reel still has to be visible after insertion, which means it has to face out.**

- 8) Repeat steps 1-7 for the second reel.

## Installing the Boa dials

Complete the following steps to install the Boa dials:

- Install the Boa dials using the Boa tool.

**INFORMATION: Make sure that the Boa dials are parallel to the socket shell.**

**INFORMATION: Take care to fully tighten the screws to the stop.**

### 6.3.5 Installing the axle

Complete the following steps to install the axle:



- 1) Insert the connecting lugs of the lateral side rail into the medial side rail.
- 2) With a hammer, tap the axle into the bore in the medial and lateral side rails.
- 3) Use a cotter pin drive to tap the axle all the way in (see fig. 26).

**INFORMATION: The axle is thicker at the end in the form of a knurling. This knurling must only be tapped into the bore as the final step. This side is also marked with a face end bore so the required removal direction can be determined if the axle has to be removed later.**

- 4) Tap the axle in until the end of the axle only just disappears into the bore.

If the axle has to be removed, the knurled end has to be pushed out of the bore first. To use the medial side rail and axle again, the axle has to be installed into the medial side rail from the other side the second time.

### 6.3.6 Screwing together the side rails and socket shells

Complete the following steps to install the side rails:

- 1) Slide the two side rails simultaneously into the side rail guides on the socket shell.
- 2) The medial side rail has to extend at least 2 cm into the pocket of the medial side rail guide.
- 3) Screw the side rails to the side rail guide on the socket shell with slot nuts, clamping pieces and countersunk screws (M4x8).
- 4) Tighten the countersunk screws to a tightening torque of 1.5 Nm and 2 Nm (see fig. 28 and see fig. 27).
- 5) Re-tighten all screws after 30 minutes and secure with 636K13 adhesive.

**INFORMATION: When securing the screws, loosen only one screw at a time, apply the adhesive to the thread and tighten the screw again to the specified tightening torque.**

The following conditions have to be met to ensure the proper strength of the connection between the side rail and socket shells:

- The socket may be lengthened by "extending" it from the stop of the side rails by no more than 2 cm.
- The minimum overlap length between the side rail guides and side rails has to be maintained: 10 cm medial, and 8 cm lateral (sum of the length of the clamping pieces).
- Two clamping pieces have to be installed on the lateral socket shell.

**The first clamping piece must be located at the lowest possible position of the overlap. The second clamping piece must be located at the highest possible position of the overlap.**

- At least one clamping piece has to be installed on the medial socket shell.

**The clamping piece must be located at the lowest possible position of the overlap.**


The medial side rail has to extend at least 2 cm into the pocket of the medial side rail guide. The top countersunk screw has to be installed (the countersink for this screw is directly in the medial shell).

- The screws are re-tightened after 30 minutes and secured with 636K13 adhesive as described above.

### 6.3.7 Installing the lace crossover parts


To protect the laces at the crossover points, lace crossover parts must be installed. Guide grooves guide the laces at the lace crossovers on different levels and avoid direct contact between the laces.

Complete the following steps to install the lace crossover part:

- 1) First snap the lace that is lying on the socket shell (lower lace) into the groove with the larger groove depth (see fig. 29).
- 2) Then snap the upper lace into the groove with the lower groove depth (see fig. 30).

**INFORMATION: Ensure that the lace is not permanently deformed when the lace is snapped into the lace crossover parts.**

- **INFORMATION:** Ensure that the lace crossover part that is appropriate for the respective lace crossover is used (labeling on the back with “L” / “R”). The lace that is in contact with the socket shell (lower lace) has to be positioned in the groove with the larger groove depth of the lace crossover part (see fig. 32)!

### 6.3.8 Installing the socket brim pad

Complete the following steps to install the socket brim pad:

- Push the socket brim pad into the corresponding recess in the medial socket shell (see fig. 33).

### 6.3.9 Screwing the bottom plate to the connection adapter

Complete the following steps to screw the magnetic bottom shell to the connection adapter:

- 1) Insert the magnetic bottom shell into the medial side rail from the inside.
- 2) Screw the magnetic bottom shell and the connection adapter to each other from the outside with the cap screws (M6) (see fig. 34).
- 3) Tighten the cap screws with a torque of 15 Nm.
- 4) Re-tighten all screws by turns after 30 minutes and secure with 636K13 adhesive.

**INFORMATION:** When securing the screws, only loosen one screw at a time, apply the adhesive to the thread and tighten the screw to the specified torque.

The connection adapter can be installed in three positions depending on the flexion contracture.

<b>Flexion contracture</b>	<b>Recommended screw position</b>
None	Low flexion position of the socket Install the connection adapter in the anterior position.
Up to approx. 5 degrees	Moderate flexion position of the socket Install the connection adapter in the middle.
Greater than 5 degrees	High flexion position of the socket Install the connection adapter in the posterior position.

The following conditions have to be met to ensure there is an adequate connection between the magnetic bottom shell, medial side rail and connection adapter:

- All six screws always have to be used.
- The screws have to be alternately tightened several times to a torque of 15 Nm until the screws no longer turn.
- The screws are re-tightened after 30 minutes and secured with 636K13 adhesive as described above.

## **6.4 Finishing the prosthetic socket**

### **Installing the socket adapter**

The socket adapter can be selected from the list of possible combinations (see page 33).

Use the supplied countersunk head screws (M6x12 and M6x16) to install the socket adapter on the connection adapter.

A medial or lateral installation position can be used for the socket adapter depending on the residual limb position. When socket adapters with a centring collar are used, two 4Y19 pressure plates have to be installed between the socket adapter and connection adapter. When a socket adapter with pyramid receiver is used, note that a collision between the received pyramid of the connecting component and the connection adapter has to be prevented. One or more 4Y19 pressure plates may have to be installed.

The longest possible screws must be used for each screw connection. If the longest screws collide with the medial side rail, the shorter screws have to be used.

## **6.5 Verification of function and fit**

### **Inspecting the liner edge**

The liner edge should be 1-2 cm higher than the socket brim.

If pressure on the medial socket brim is excessive when checking the fit, the socket can be shortened. This shifts the pressure to the distal area of the socket.

If the pressure on the residual limb end is too high when checking the fit, the socket can be extended.

## Exchanging fabric adhesive strips for plastic adhesive strips

35

If the prosthetic suspension between the 6Y200=\* Varos liner and prosthetic socket is not sufficient (common for people with a higher mobility grade of MOBIS 2 or more and a short residual limb), the fabric adhesive strips can be replaced by plastic adhesive strips.

The highly adhesive plastic strips make it more difficult to take the socket off. Therefore, it is recommended to initially replace only the fabric adhesive strips in the lateral socket shell with plastic adhesive strips. If the prosthetic suspension is still not sufficient, the fabric adhesive strips in the medial socket shell can also be replaced.

To replace the fabric adhesive strips, remove them from the socket and remove the adhesive residue. Clean the recesses in the socket shells with isopropanol. Glue the plastic adhesive strips into the recesses in the socket shells. In doing so, glue the plastic adhesive strip into the upper and lower area of the recessed groove. Note the alignment of the plastic adhesive strips (observe the label "TOP").

### Sanding the socket brim pad

The socket brim pad can be sanded if necessary. Remove the socket brim pad from the socket for sanding. A minimum wall thickness of 1.5 mm has to be maintained to ensure sufficient hold of the socket brim pad.

#### INFORMATION

Sanding the socket brim pad creates a rough surface that needs to be sanded smooth with a fleece abrasive wheel on the socket router.

### Prosthetic alignment note

Refraining from bringing the residual limb into an adduction position has proven to be better for the prosthetic alignment of the socket. The residual limb should be brought to a neutral position. A point to the side next to the upper dial was assumed as the proximal socket reference point. This point marks the height of the hip joint and lies on the longitudinal socket axis (if no rotation via a socket adapter is installed).

We recommend that no rotation is installed between the socket and knee joint. The knee axis should be rotated slightly to the inside compared to the longitudinal axis of the foot (patella aligned in the direction of the gap between the big toe and second toe).

If the tip of the foot points too far to the outside or inside after putting on the prosthesis, even though the preceding recommendations were followed, this usually just means that the socket was applied in an incorrect position relative to the residual limb and the donning process has to be repeated.

## 7 Use

### 7.1 Putting on the liner and socket

#### Putting on the liner

- 1) Check the liner and magnetic cap for damage.
- 2) Turn the liner inside out and apply it to the end of the residual limb.
- 3) Apply the liner to the end of the residual limb and check the alignment of the liner at the same time.

**INFORMATION: The orientation thread on the inner side of the liner has to be centred on the front of the transfemoral residual limb. The liner's two seams have to be positioned on the sides of the transfemoral residual limb.**

- 4) Unroll the liner onto the residual limb so that it has no wrinkles or air pockets, without shifting the soft tissue.
- 5) Check the fit and the alignment of the liner.
- 6) In a standing position (or position with the residual limb nearly extended), unroll the liner part-way and roll it back up again to release possible skin tension.

#### Putting on the socket

- 1) To open, pull the Boa dials away from the prosthetic socket.
- 2) Open the socket by pulling the shells apart.
- 3) Sit on the front edge of a chair.
- 4) Pull the socket over the residual limb.

**INFORMATION: In doing so, note the alignment of the toes and the position of the lower leg so the desired level of exterior foot rotation is given while standing and walking. If the exterior foot rotation is not correct, take the socket off again (at least, the magnetic cap has to be released from the bottom of the socket) and put it back on in a changed position.**

- The magnetic bottom shell of the socket and the magnetic cap of the liner attract each other when they are brought close together. A clicking sound is heard, indicating the proper connection of the socket and liner.
- 5) Push the Boa dials towards the socket.
  - 6) Turn the Boa dials clockwise.

→ The Boa cables are tightened and the socket is closed.

**INFORMATION: Position the outer socket shell (lateral) under the inner socket shell (medial).**

- 7) Tighten the Boa dials as needed.

### 7.2 Taking off the socket and liner

- 1) Pull the Boa dials away from the socket.
- 2) Pull the outer socket shell (lateral) to the outside.  
If plastic adhesive strips are used in the lateral socket shell, pull until there is no longer any contact between the plastic adhesive strips and textile cover of the liner.
- 3) Push the socket off the residual limb.  
In doing so, preferably grasp the socket brim of the lateral socket shell with one hand and the socket brim of the medial socket shell with the other hand.
- 4) Unroll and remove the liner from the residual limb.

## 8 Cleaning

- 1) Clean the product with a damp cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed.
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

## 9 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Perform annual safety inspections (visual inspection for damage and wear and tear, functional test).

## 10 Legal information

### 10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 10.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

## 11 Technical data

<b>Ambient conditions</b>	
Transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage in original packaging	-25 °C/-13 °F to +50 °C/+122 °F Max. 90% relative humidity, non-condensing
Storage and transport without packaging	-25 °C/-13 °F to +50 °C/+122 °F Max. 90% relative humidity, non-condensing
Operation	-10 °C/+14 °F to +40 °C/+104 °F Max. 90% relative humidity, non-condensing
Water resistance	Fresh water, mild soap
Not waterproof	Salt water, chlorine water

<b>General</b>	
Reference number	5A60=*
Lifetime	3 years Wear and tear part that is subject to natural degeneration. The lifetime decreases or increases depending on the individual level of use. The maximum lifetime – which depends on the level of use – can only be achieved when the instructions for use are observed.

<b>General</b>	
Water resistance	Contact with fresh water and mild soap is permitted. Not designed for chlorinated water or salt water.
Magnetic bottom shell material	INOX (stainless steel)
Connection adapter material	Anodised aluminium
Medial and lateral bar material	Carbon fibre reinforced polyamide PA12
Medial and lateral socket shell material	Polyamide PA12
Clamping piece material	Die-cast zinc coated
Slot nut material	Stainless steel
Socket brim pad material	EVA blend
Lace crossover part material	Polyamide PA12

<b>Socket</b>	<b>5A60=*-S</b>	<b>5A60=*-SC</b>	<b>5A60=*-M</b>	<b>5A60=*-MC</b>	<b>5A60=*-L</b>
Adjustable socket length (in mm)	200–320	200–320	200–320	200–320	200–320
Adjustable proximal socket circumference (in mm)	410–510	460–580	460–580	520–650	520–650
Product weight (not shortened) (in g)	1260	1310	1320	1450	1490
Adjustable flexion position (socket adapter to longitudinal socket axis) (in degrees)	1.5°/4.5°/7.5°	1.5°/4.5°/7.5°	1.5°/4.5°/7.5°	1.5°/4.5°/7.5°	1.5°/4.5°/7.5°
Adjustable medial/lateral offset of the socket adapter (in mm)	10	10	10	10	10
Structural height (distance between residual limb end/connection adapter surface) (in mm)	30	30	30	30	30

## 12 Appendices

### 12.1 Symbols Used

Declaration of conformity according to the applicable European directives

Manufacturer

Lot number (PPPP YYYY WW)  
 PPPP – plant  
 YYYY – year of manufacture  
 WW – week of manufacture

Please note the instructions for use

Medical device

Article number

## 12.2 Spare parts

- Medial side rail size S/SC, art. no. 5E1=S  
 Medial side rail size M/MC, art. no. 5E1=M  
 Medial side rail size L, art. no. 5E1=L
  - Lateral side rail, straight, art. no. 5E2=1
  - Lateral side rail, curved, art. no. 5E3=1
  - Axle, art. no. 5E4=1
  - Magnetic bottom shell size S/SC/M/MC, art. no. 5E5=M  
 Magnetic bottom shell size L, art. no. 5E5=L
  - Connection adapter, art. no. 5E6=1
  - Medial socket shell size S, art. no. 5E15=\*-S  
 Medial socket shell size SC/M, art. no. 5E15=\*-M  
 Medial socket shell size MC/L, art. no. 5E15=\*-L
  - Lateral socket shell size S, art. no. 5E16=\*-S  
 Lateral socket shell size SC/M, art. no. 5E16=\*-M  
 Lateral socket shell size MC/L, art. no. 5E16=\*-L  
 Contains: screwed BOA housing, 2 pc. BOA dials, 2 pc. BOA reels
  - BOA Fit System, article number 5E11=1  
 Contains 2 BOA dials with reels, housings, assembly part, BOA lace set, 3 lace crossover parts "L" and 3 lace crossover parts "R".
  - BOA lace set, art. no. 5E12=1  
 Contains: BOA laces in these lengths: 900 mm, 1000 mm, 1300 mm and 1600 mm
  - Fabric adhesive strip set size S/SC/M, art. no. 5E7=M  
 Fabric adhesive strip set size MC/L, art. no. 5E7=L
  - Socket brim pad size S, art. no. 5E17=\*-S  
 Socket brim pad size SC/M, art. no. 5E17=\*-M  
 Socket brim pad size MC/L, art. no. 5E17=\*-L
  - Clamping piece set, art. no. 5E8=1  
 Contains: 3 pc. clamping pieces, 7 pc. slot nuts, 7 pc. M4x8 countersunk screws
  - Rubber nub set, art. no. 5E13=1  
 Contains: 14 pc. rubber nubs
  - Connection adapter screw set, art. no. 5E18=1  
 Contains: 6 pc. M6x10 cap screws
  - Socket adapter screw set, art. no. 5E14=1  
 Contains: 4 pc. M6x12 countersunk screws, 4 pc. M6x16 countersunk screws
- The "\*" placeholder in the article number stands for the "left" or "right" variant.

### **12.3 Accessories**

- Magnetic cap, threaded version size S/SC/M/MC, art. no. 5E10=M  
Magnetic cap, threaded version size L, art. no. 5E10=L
- Set of plastic adhesive strips size S/SC/M, art. no. 5E9=M  
Set of plastic adhesive strips size MC/L, art. no. 5E9=L  
Contains: 10 pc. plastic adhesive strips

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