F24=* Maverick junior

DE	Gebrauchsanweisung	7
EN	Instructions for use	16
FR	Instructions d'utilisation	25
IT	Istruzioni per l'uso	34
ES	Instrucciones de uso	44
PT	Manual de utilização	53
NL	Gebruiksaanwijzing	62
SV	Bruksanvisning	71
DA	Brugsanvisning	79
NO	Bruksanvisning	87
FI	Käyttöohje	96
PL	Instrukcja użytkowania	104
HU	Használati utasítás	113
CS	Návod k použití	122
RO	Instrucțiuni de utilizare	131
HR	Upute za uporabu	140
SL	Navodila za uporabo	149
SK	Návod na používanie	157
BG	Инструкция за употреба	166
TR	Kullanma talimatı	176
EL	Οδηγίες χρήσης	184
RU	Руководство по применению	194
JA	取扱説明書	204
ZH	使用说明书	213
KO	사용 설명서	219


ISO 22675 - P n - m kg*)

*) Body mass limit not to be exceeded!

Stiffness	Weight (kg)	Label text
1	15	ISO 22675-P3*-15-kg
2	24	ISO 22675-P3*-24-kg
3	34	ISO 22675-P3*-34-kg
4	45	ISO 22675-P3*-45-kg

1

Bench alignment, transtibial

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INFORMATION

Date of last update: 2024-07-01

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The F24=* Maverick junior prosthetic foot was designed for small foot sizes. The spring elements are made of fibreglass. Depending on the variant, the prosthetic foot is equipped either with a pyramid or a lamination anchor.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

The table that follows shows the suitable spring stiffness of the prosthetic foot, matching the patient's body weight and the size of the prosthetic foot.

Spring stiffness relative to body weight						
Body weight [kg]	Prosthetic foot size [cm]					
	16	17	18	19	20	21
9 to 15	1			-		
16 to 24	2				-	
25 to 34	-	3				

Spring stiffness relative to body weight						
Body weight [kg]	Prosthetic foot size [cm]					
	16	17	18	19	20	21
35 to 45	-		4			

2.3 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors for children. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

- **Combination with 2R37, 2R38:** Use the reduction 4X907=1 to connect the pyramid of the modular system for children to the larger pyramid receiver.
- Combination with a mechatronic prosthetic knee joint: Perform a fitting to ensure that the prosthetic knee joint works correctly for the user. Particularly make sure that the prosthetic knee joint activates the swing phase at the end of the stance phase.
- **The product cannot be combined with the following components:**
2R41=2

2.4 Environmental conditions

Allowable environmental conditions
Temperature during use: -10 °C to +45 °C
Chemicals/liquids: Fresh water, salt water, perspiration, urine, soapsuds, chlorine water
Moisture: Submersion: max. 1 h in 3 m depth, relative humidity: no restrictions
Solids: Dust, occasional contact with sand
Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and tear and damage (see page 23).
Prohibited environmental conditions
Solids: Highly liquid-binding particles (e.g. talcum), continuous contact with sand
Chemicals/liquids: Acids, continuous use in liquid media
UV light
If the product is exposed to UV light, discolouration may occur.
Storage and transport
Storage temperature: -20 °C to +60 °C, relative humidity: 20 % to 90 %, no mechanical vibrations or impacts

2.5 Lifetime

Prosthetic foot

Depending on the patient's activity level, the maximum lifetime of the product is 3 years.

Footshell, protective sock

The product is a wear part, which means it is subject to normal wear and tear.

3 General safety instructions

CAUTION! Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ Note the combination possibilities/combo exclusions in the instructions for use of the products.
- ▶ Observe the maximum lifetime of the product.
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not expose the product to excessive strain (section: "Area of application" see page 16) or to prohibited environmental conditions (see page 17).
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ The product may only be used by one person; reuse on another person is prohibited.

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover behaviour) or delamination of the spring are indications of loss of functionality. Unusual noises can indicate a loss of functionality.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	-
1	Prosthetic foot	-
1	Protective sock	2U3=1-7
1	Insert for heel wedge	2F14=*
1	Tool to remove the footshell	2C101
1	Quick reference guide	647S313=1

Only for prosthetic foot with lamination anchor		
Quantity	Designation	Reference number
1	Lamination anchor for modular system for children, with rotation	2R89=1
1	Cap screw, titanium	501T24=M4x20

Spare parts/accessories	
Designation	Reference number
Footshell	2C24=*

5 Preparing the product for use

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

NOTICE

Altering the prosthetic foot or footshell

Premature wear resulting from damage to the product

- ▶ Do not alter the prosthetic foot or footshell.

INFORMATION

A plastic part protects the connecting section of the product from scratches during the alignment and fitting.

- ▶ Remove the plastic part before the patient leaves the fitting area.

INFORMATION

A steel cap screw is pre-installed in the lamination anchor. The enclosed titanium cap screw can improve corrosion resistance.

- ▶ **Steel:** Use the steel cap screw for fitting. The screw may be used for the definitive prosthesis but is not corrosion-resistant.
- ▶ **Titanium:** Only use the titanium cap screw after the fitting if no further adjustments are required.

5.1 Applying/removing the footshell

INFORMATION

- ▶ Pull the protective sock over the prosthetic foot to prevent noises in the footshell.
- ▶ **CAUTION! Always use the prosthetic foot with the footshell. CAUTION! When handling the foot module, ensure that you do not get any splinters into your hands.**
Apply or remove the footshell as described in the footshell instructions for use.

5.2 Bench Alignment

Bench alignment, amputation through the knee or above the knee

- ▶ Observe the instructions for use for the prosthetic knee joint.

Bench alignment, transtibial prosthesis

Position the prosthetic foot and prosthetic socket in the alignment apparatus

- > **Required materials:** goniometer 662M4, heel height measuring device 743S12, 50:50 gauge 743A80, alignment apparatus (e. g. PROS.A. Assembly 743A200)
- 1) Positioning the prosthetic foot in the alignment apparatus:
 - **Heel height: Effective heel height of everyday shoe without additional elevation under the heel.**
- 2) Align the prosthetic foot with the alignment reference line:
 - **Sagittal plane:** Position the alignment reference line on the markings of the footshell.
 - **Frontal plane:** Position the alignment reference line on the markings of the footshell (**external rotation approx. 5°**).
- 3) Mark the reference lines on the prosthetic socket:
 - **Sagittal plane:** Determine the centre of the prosthetic socket proximally and distally with the 50:50 gauge and mark the centre line.
 - Mark the socket reference point on the centre line: **at the height of the centre of the patella**
- 4) Position and align the prosthetic socket in the alignment apparatus:
 - **Sagittal plane:** The alignment reference line runs vertically through the socket reference point.
 - Adjust the socket flexion by rotating around the socket reference point: **individual residual limb flexion + 5°**

- **Frontal plane:** Align the prosthetic socket with the alignment reference line, taking the abduction or adduction position into account.

Prosthetic foot with pyramid

- Install the structural components between the prosthetic foot and prosthetic socket.

Prosthetic foot with lamination anchor

NOTICE! Do not use the lamination anchor for height compensation. Screw the prosthetic foot fully into the lamination anchor and then max. 180° back.

- > **Required materials:** Bending iron 711S4*, degreasing cleaner (e. g. acetone), monofilament adhesive tape 627B2=19 (alternatively: Cellocast 699G30=*)
- 1) Measure the distance between the end of the prosthetic socket and the floor. The distance must correspond to the prosthetic foot:
 - **Required distance for size 16 to 18: 50 mm**
 - **Required distance for size 19 to 21: 55 mm**
 - 2) Set the prosthetic socket in the alignment apparatus to the socket-floor measurement, position it in the lamination anchor of the prosthetic foot and align it.
 - 3) Screw the prosthetic foot fully into the lamination anchor. In doing so, make sure the thread does not tilt.
 - 4) Unscrew the prosthetic foot by no more than half a turn to align it. Do not align the threaded clamping in the anterior or posterior direction.
 - 5) **NOTICE! Fully screw in the red pyramid receiver before bending and tighten the cap screw of the threaded clamping (3 Nm). Otherwise, the lamination anchor warps during bending.**
Adapt the arms of the lamination anchor to the prosthetic socket using the bending irons.
 - 6) Roughen the distal end of the socket.
 - 7) Clean the lamination anchor using a degreasing cleaner.
 - 8) Slightly loosen the cap screw on the threaded clamping.
 - 9) Cover the prosthetic foot to prevent it from getting dirty.
 - 10) Prepare a viscous sealing resin/talcum mixture.
 - 11) Apply the sealing resin/talcum mixture to the lamination anchor and prosthetic socket, compensating for small height deviations.
 - 12) Place the prosthetic socket back in the lamination anchor and align it.
 - 13) Allow the sealing resin/talcum mixture to cure and then remove the prosthesis from the alignment apparatus.

- 14) Apply a further sealing resin/talcum powder mixture around the arms of the lamination anchor and let it cure.
- 15) Smooth the area and wrap it with fibreglass-reinforced adhesive tape or Cellacast for reinforcement.

5.3 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.

5.4 Dynamic Trial Fitting

The optimal gait pattern is developed during dynamic fitting. The prosthesis is aligned optimally in the frontal plane and the sagittal plane for this purpose.

- ▶ **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal planes is achieved after the leg begins to bear weight following the heel strike. Avoid medial movement of the knee joint.
 - If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction.
 - If the knee joint moves in the medial direction in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.
- ▶ After completing the dynamic fitting and the walking exercises: Remove the plastic part from the pyramid.

Prosthetic foot with lamination anchor

If the alignment of the prosthesis is appropriate, the prosthetic socket is finished after the fitting. The temporary attachment of the lamination anchor is removed and the anchor is definitively laminated.

5.4.1 Optimising the heel characteristics

The prosthetic foot has a permanently installed heel wedge. An insert for the heel wedge makes it possible to adapt the behaviour of the prosthetic foot at heel strike and during heel contact in the mid stance.

- ▶ Install or remove the insert for the heel wedge on the side of the prosthetic foot (see fig. 3).

5.5 Installing the lamination anchor

- > **Required materials:** special lubricant 633G6, Loctite 636K6, torque wrench
- 1) Seal the cap screw, threaded clamping and thread of the lamination anchor with special lubricant.
- 2) Screw the red pyramid receiver into the lamination anchor to the stop. In doing so, make sure the thread does not tilt.

- 3) Tighten the cap screw of the threaded clamping to a point where the pyramid receiver is firmly seated, but can still be unscrewed and screwed in.
- 4) **INFORMATION: Do not modify the lamination dummy or remove the studs serving as spacers.**
Place the lamination dummy on the lamination anchor.
- 5) Position the lamination anchor on the prosthetic socket and reinforce it.
 - The type of reinforcement depends on the expected load.
 - Tie off the individual stockinette layers in the groove below the lamination dummy (see fig. 5).
- 6) Carry out the lamination process.
- 7) Once the resin has cured, remove the lamination dummy and sand the edges.
- 8) Unscrew the red pyramid receiver from the lamination anchor.
- 9) **If necessary:** Remove resin residue from the thread.
- 10) Screw the prosthetic foot fully into the lamination anchor. In doing so, make sure the thread does not tilt.
- 11) Unscrew the prosthetic foot by no more than half a turn to align it.
- 12) Secure the cap screw of the threaded clamping using Loctite and tighten it (3 Nm).

INFORMATION: To improve the corrosion resistance of the prosthesis after the fitting, replace the steel cap screw with the titanium screw.

6 Cleaning

- > **Allowable cleaning agent:** pH-neutral soap (e.g. 453H10 Derma Clean)
- 1) **NOTICE! To avoid product damage, only use the allowable cleaning agents.**
Clean the product with clear fresh water and a pH neutral soap.
 - 2) Remove debris from water drainage contours using a toothpick and rinse.
 - 3) Rinse the soap away with clear fresh water. In doing so, rinse the foot-shell until all debris has been removed.
 - 4) Dry the product with a soft cloth.
 - 5) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.

- ▶ During regular inspection: Inspect the prosthesis for signs of wear and check its function.
- ▶ Conduct annual safety inspections.

8 Disposal

Do not dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

10 Technical data

Sizes [cm]	16	17	18	19	20	21
Heel height [mm]	5					
System height with pyramid [mm]	26			31		
Build height with pyramid [mm]	42			47		
Required clearance for prosthetic foot with lamination anchor [mm]	50			55		
Weight with footshell [g]	249	263	285	305	328	352
Max. body weight [kg]	24	34	45			

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