




















50K306=* Agilium Vantage



DE Gebrauchsanweisung	7
EN Instructions for use	12
FR Instructions d'utilisation	17
IT Istruzioni per l'uso	22
NL Gebruiksaanwijzing	27
SV Bruksanvisning	32
DA Brugsanvisning	37
NO Bruksanvisning	42
FI Käyttöohje	47
PL Instrukcja użytkowania	52
ET Kasutusjuhhis	57


	Größe / Size	Unterschenkelumfang / Lower leg circumference		Oberschenkelumfang / Thigh circumference	
		A		B	
		[cm]	[inch]	[cm]	[inch]
	XS	28-32	11-12.5	36-41	14-16
	S	32-36	12.5-14	41-46	16-18
	M	36-41	14-16	46-53	18-21
	L	41-46	16-18	53-60	21-23.5
	XL	46-51	18-20	60-67	23.5-26.5
	XXL	51-56	20-22	67-75	26.5-29.5
	XXXL	56-61	22-24	75-81	29.5-32

Kombinierte Größen/ Combined Sizes	Mediale Kompartiment Arthrose/ Medial compartment osteoarthri- tis		Laterale Kompartiment Arthrose/ Lateral compartment osteoarthri- tis	
	Links/Left	Rechts/Right	Links/Left	Rechts/Right
	50K306*			
XS/S	*=L-XS-S-MO	*=R-XS-S-MO	*=L-XS-S-LO	*=R-XS-S-LO
M/L	*=L-M-L-MO	*=R-M-L-MO	*=L-M-L-LO	*=R-M-L-LO
XL/XXL/XXXL	*=L-XL-3XL-MO	*=R-XL-3XL-MO	*=L-XL-3XL-LO	*=R-XL-3XL-LO

Material	Aluminium/aluminium, Stahl/steel, ABS/ABS, Nylon/nylon, PU Schaumstoff/PU foam, Elasthan/spandex
-----------------	--

INFORMATION

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INFORMATION

Date of last update: 2024-06-25

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50K306 Agilium Vantage.

2 Product description

50K306 Agilium Vantage (see fig. 1)

Item	Designation	Item	Designation
1	Thigh pad	9	Positioning aid
2	Flap 4	10	Clip
3	Y-hook-and-loop	11	Release strap
4	Flap 1	12	Orthotic joint
5	Patella cutout	13	Anti-rotation strap
6	Flap 2	14	D-ring for anti-rotation strap
7	Flap 3	15	D-ring for release strap
8	Calf pad	-	Extension and flexion stops (see fig. 18)

3 Intended use

3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

3.2 Indications

- Unicompartmental osteoarthritis of the knee

- Conditions/injuries requiring unicompartamental relief (e.g. post-operative treatment following meniscus reconstruction or ligament injuries which require unilateral relief)
- Rheumatoid arthritis

Indications must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

Not known.

3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; pronounced varicose veins, especially with return flow impairment, lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support will be applied; sensory and circulatory disorders in the area of the legs, e.g. in case of diabetic neuropathy.

- Bicompartamental osteoarthritis of the knee

3.4 Lifetime

The product is designed for a maximum lifetime of **1 year**.

3.5 Effects

The orthosis provides relief for the affected compartment by means of the 3-point principle and stabilises the knee joint.

4 Safety

4.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

4.2 General safety instructions

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

Use with heat-sensitive skin

Skin irritation due to overheating

- ▶ Do not use the product in case of a known heat allergy.
- ▶ Do not continue to use the product if skin irritation occurs.

Exceeding the lifetime and reuse

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the maximum lifetime is not exceeded.
- ▶ Only use the product for one person.

⚠ CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

5 Use

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the treating physician.
- ▶ Consult a physician if any exceptional changes are noted (such as worsening of the complaint).

5.1 Size selection

- 1) Measure the circumference of the lower leg **15 cm** below the centre of the patella (circumference A).
- 2) Measure the circumference of the thigh **15 cm** above the centre of the patella (circumference B).
- 3) Determine the size of the orthosis (see size chart).

5.2 Adapting the size

Trim the product to use it in sizes XS, M, XL and XXL (see fig. 22).

- ▶ Use scissors to cut the thigh and calf pads (item 1, item 8) to the desired size at the marks.
- Trim size S to XS.
- Trim size L to M.
- Trim size XXXL to XXL or XL.

5.3 Fitting and application

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- > **Prerequisite:** The patient should be seated on the edge of a chair.
- > **Prerequisite:** The knee joint should be relieved of pressure and slightly bent.
- 1) Open all hook-and-loop closures on the orthosis.
- 2) Centre the patella cutout (item 5) over the patella (see fig. 4).
- 3) Bend the leg to a flexion position of **30°** with the foot flat on the floor (see fig. 5).
- 4) Place the orthotic joint (item 12) medially (MO) or laterally (LO) on the side being treated, positioning the centre of the orthotic joint at the height of the middle of the patella.
- 5) Use one hand to hold the orthosis on the leg (see fig. 6).
- 6) Use the other hand to tighten flap 1 (item 4) and attach it to mark 1 on the thigh pad (item 1) with the hook-and-loop closure. The hook-and-loop sections should overlap (see fig. 7).
- 7) Then, tighten flap 2 (item 6) and attach it to mark 2 on the calf pad (item 8) with the hook-and-loop closure (see fig. 8).
- 8) Tighten flap 3 (item 7) and attach it to mark 3 on the calf pad (item 8) with the hook-and-loop closure (see fig. 9).
- 9) Tighten flap 4 (item 2) and attach it to mark 4 on the thigh pad (item 1) with the hook-and-loop closure (see fig. 10).
- 10) Check the fit of the orthosis and adjust the flap positions, if necessary.
- 11) **Optional:** Use scissors to change the size or fit of the thigh and calf pads.
- 12) **Optional:** Use bending irons to change the fit of the orthotic joint bars. **INFORMATION: The joint bars are permanently integrated in the orthosis and cannot be removed from the orthosis. The orthotic joint must not be bent, as this results in loss of function of the orthotic joint.**

Adjusting the straps

- 1) Use one hand to tighten the release strap (item 11) and the other hand to slide the clip (item 10) on the strap.
- 2) Position the clip so it is centred on the side of the thigh (see fig. 12).
- 3) Guide the strap through the D-ring (item 15) of mark 5 and move the leg to a flexion position of **80°** with the foot flat on the floor (see fig. 13).
- 4) Hold the thigh pad (item 1) in place with one hand and fasten the strap with the other hand (see fig. 14).
- 5) Remove the positioning aid (item 9) from the clip; this fixes the clip to the strap (see fig. 15).
- 6) Move the leg back to a flexion position of **30°**.
- 7) Thread the anti-rotation strap (item 13) through the D-ring (item 14) (see fig. 16).
- 8) Tighten the strap and fasten the strap to itself (see fig. 17). **INFORMATION: Tighten the strap until maximum pain relief is achieved.**
- 9) **NOTICE! Once the straps have been shortened, this cannot be undone. Proceed with care!** Ensure an optimal fit by removing the Y-hook-and-loop from the respective strap and using scissors to shorten the end of the strap to the required length.
- 10) Reattach the Y-hook-and-loop to the strap end.

Optional: limiting the knee flexion and extension

The extension limitation at delivery is **0°**. The flexion is not limited.

Available extension stops (square grip):	0°, 10°, 20°, 30°, 45°
Available flexion stops (round grip):	0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

- Change the extension and flexion limits only if advised to do so by the attending physician.
- 1) Remove the cover (item 12) from the orthotic joint (see fig. 19).
- 2) Remove the placeholder from the orthosis if present.
- 3) Insert the required extension stop, alternatively the flexion stop (see fig. 18).
→ The stop will engage with slight resistance.
- 4) Place the cover on the orthotic joint, ensuring the arrows are pointing in the correct direction on the cover. The arrows should point in the posterior direction (see fig. 21).
- 5) Move the orthotic joint and check the position of the stops.

5.4 Delivery

NOTICE

Incorrect setting

Damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ The orthosis may only be adjusted by qualified personnel.
- ▶ Do not make any improper changes to the settings.

- 1) Check the product fits properly on the leg when handing the orthosis over to the patient.
- 2) Ask the patient to sit down, stand up and walk a few steps.
 - The centre of the orthotic joint should be at the height of the middle of the patella.
 - The orthosis should not be twisted on the leg.
 - The straps should be tight but not impair the patient's comfort.

Inform the patient:

Inspect the orthosis for damage

- ▶ Inspect the orthosis for damage daily. Contact the qualified personnel promptly if changes are noted.

5.5 Removal

- 1) Open both straps at the hook-and-loop closures and remove them from the D-rings.
- 2) Open flaps 1 to 4, one after the other, and remove the orthosis from the leg.

6 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the product promptly after contact with water containing salt or chlorine, or if it gets dirty.

Cleaning the orthotic joint

- ▶ The orthotic joint may only be dismantled by qualified personnel.
- 1) Rinse the orthotic joint with clean, fresh water and remove dirt.
 - 2) Dry with a cloth or allow to air dry.

Cleaning the orthosis

- 1) **Hand wash** the orthosis in warm water at **30 °C (86 °F)** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 3) Do not tumble-dry or bleach.

7 Disposal

Dispose of the product in accordance with national regulations.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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