17AD300=* Nexgear Tango ankle

DE	Gebrauchsanweisung	17
EN	Instructions for use	36
FR	Instructions d'utilisation	56
IT	Istruzioni per l'uso	77
ES	Instrucciones de uso	96
PT	Instruções de utilização	117
NL	Gebruiksaanwijzing	137
SV	Bruksanvisning	157
DA	Brugsanvisning	176
NO	Bruksanvisning	195
FI	Käyttöohje	213
PL	Instrukcja użytkowania	232
HU	Használati útmutató	252
CS	Návod k použití	272
HR	Upute za uporabu	291
SK	Návod na použitie	310
TR	Kullanım kılavuzu	330
RU	Руководство по применению	349
JA	取扱説明書	370
ZH	使用说明书	389


PDF-ONLINE


DE, EN, FR, IT, ES, PT, NL, SV,
DA, NO, FI, PL, HU, CS, HR, SK,
TR, RU, JA, ZH-CHS

<https://product-documents.ottobock.com/IFU/INT/17AD300/647G1788/06/O/S/F>

17AD300=*

		14	16	20
 Pos. 1		T20	T25	T25
		3,5 Nm	7 Nm	7 Nm


		14	16	20
 Pos. 4		T25	T30	T30
		5 Nm	8 Nm	8 Nm


3

17AD313=*
①


4

17AD310=*
⑥


17AD301=*

 Pos. 1	14	16	20
	T25	T30	T30

 Pos. 4	14	16	20
	SW12 6 Nm	SW13 8 Nm	SW13 8 Nm

17AD302=*


 Pos. 1	14	16	20
	T25	T30	T30

 Pos. 6	14	16	20
	SW12 6 Nm	SW13 8 Nm	SW13 8 Nm

17AD303=*

17AD304=*








 Pos. 2	14	16	20
	T30	T50	T50







 Pos. 3	14	16	20
	SW2,5 2 Nm	SW3 3 Nm	SW3 3 Nm

 Pos. 8	14	16	20
	T10	T15	T15

 Pos. 11	14	16	20
	SW20 6 Nm	SW24 8 Nm	SW24 8 Nm






21



22


39

INFORMATION

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INFORMATION

Date of last update: 2024-08-26

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with information on the processing of the Nexgear Tango ankle 17AD300=* orthotic ankle joint.

1.1 Notes regarding the document

The document is available in the following versions:

- 647G1788=all_INT (document reference number with all available languages), A6 format
All illustrations in this document are found at the beginning.
They are followed by the texts in all languages.
This printed document is included in the scope of delivery.
- 647G1788=XX_INT (document reference number as single-language variant), A5 format
647G1788=en_INT (example of single-language variant in en = English)
All illustrations in this document are found in the respective sections in colour.
This document is only available in digital form as a PDF file.
- XX (variable for the language codes of the languages in which the single-language version is available)

The digital instructions for use are available as a single language version via the following QR code or link (see fig. 1).

2 Product description

2.1 Available sizes

INFORMATION

- ▶ Use the next larger size if the product is being worn unilaterally and in case of flexion contractures in the knee or hip $>10^\circ$, distinct torsion or varus/valgus instabilities or increased physical activity.
- ▶ Use the product bilaterally if there is an ischial support.

The size of the product is selected according to the planned use (unilateral or bilateral) and the body weight.

Area of application	Size 14	Size 16	Size 20
	Maximum body weight in kg (lbs)		
Unilateral use (lateral or medial)	50 (110)	85 (187)	110 (242.5)
Bilateral use (lateral and medial)	80 (176)	120 (264.5)	160 (352.5)

2.2 Components/design

Scope of delivery (see fig. 2)			Size 14	Size 16	Size 20
Item	Quantity	Designation	Article number		
	1	Nexgear Tango	17AD300=14-ST 17AD300=14-T	17AD300=16-ST 17AD300=16-T	17AD300=20-ST 17AD300=20-T
		Consisting of:			
1	2	Countersunk screws			
2	2	Blind cap			
3	1	Base body			
4	1	Joint screw			
5	1	Joint nut			
6	1	Axial washers set			
7	1	Lamination dummy – upper section			

2.3 Accessories

Not in scope of delivery (see fig. 2)			Size 14	Size 16	Size 20
Item	Quantity	Designation	Article number		
8	1	Foot stirrup, straight ¹⁾	17SF300=N-14	17SF300=N-16	17SF300=N-20
9	1	Lamination bar ²⁾	17LS3=14 17LS3=14-T	17LS3=16 17LS3=16-T	17LS3=20 17LS3=20-T
¹⁾ Observe the instructions for use 647G1804=all_INT.					
²⁾ Observe the instructions for use 647G1425=all_INT.					

Lamination dummy set

(see fig. 3)			Size 14	Size 16	Size 20
Item	Quantity	Designation	Article number		
	1	Lamination dummy set	17AD313=14	17AD313=16	17AD313=20
		Consisting of:			
1	1	Lamination dummy – lower section			
2	1	Slotted oval head screw			
3	1	Insertion aid for joint nut			
The upper section of the lamination dummy is not part of the lamination dummy set but is in the scope of delivery of the orthotic ankle joint 17AD300=* (see fig. 2).					

Shaping dummy set

(see fig. 4)			Size 14	Size 16	Size 20
Item	Quantity	Designation	Article number		
	1	Shaping dummy set	17AD310=14	17AD310=16	17AD310=20
		Consisting of:			
1	1	Shaping dummy lower section			
2	4	Countersunk screws			
3	1	Countersunk screw			
4	1/1	Shaping dummy upper section for 17AD300 (L/R*)			
5	2	Shaping dummy upper section for 17AD303 (Reaction-module, size 14/16/20)			
6	2	Shaping dummy upper section for 17AD304 (Ultra-module, size 16/20)			
*) L=left, R=right					

Stop-module

(see fig. 5)			Size 14	Size 16	Size 20
Item	Quantity	Designation	Article number		
	1	Stop-module	17AD301=14	17AD301=16	17AD301=20
		Consisting of:			
1	1	Stop pin			
2	1	Spring			
3	1	Stop bolt			
4	1	Screw-in sleeve			
5	1	Cover			

Spring-module

(see fig. 6)			Size 14	Size 16	Size 20
Item	Quantity	Designation	Article number		
	1	Spring-module	17AD302=14	17AD302=16	17AD302=20
		Consisting of:			
1	1	Set screw			
2	1	Compression spring	—	x	x
3	1	Compression spring			
4	1	Spring guide housing			
5	1	Ball			
6	1	Screw-in sleeve			
7	1	Cover			

Reaction-module, Ultra-module

(see fig. 7)			Size 14	Size 16	Size 20
Item	Quantity	Designation	Article number		
	1	Reaction-module	17AD303=14	17AD303=16	17AD303=20
	1	Ultra-module	—	17AD304=16	17AD304=20
		Consisting of:			
1	1	Cover cap for cover			
2	1	Spring bolt			
3	1	Counter screw			
4	1	Reaction-Spring, blue (accessory)			
5	1	Reaction-Spring, yellow (accessory)			
6	1	Spring dummy for Reaction-module and Ultra-module (accessory)			
7	1	Ultra-Spring, grey (accessory)			
8	1	Spring holder			
9	1	Spring			
10	1	Stop bolt			
11	1	Screw-in with O-ring			
12	1	Cover			

Spring dummy set

(see fig. 8)			Size 14	Size 16	Size 20
Item	Quantity	Designation	Article number		
	1	Spring dummy set	17AD308=14	17AD308=20	17AD308=20
		Consisting of:			
1	1	Reaction-module mounting adapter	x	x	x
	1	Ultra-module mounting adapter	—	x	x
2	2	Reaction-module spring dummy	x	x	x
	1	Ultra-module spring dummy	—	x	x
3	2	Guide screw			

Springs for Reaction- and Ultra-module

Springs for Reaction-/Ultra-module (see fig. 9, see fig. 10)			Size 14	Size 16	Size 20
Item	Quantity	Designation	Article number		
	1	Reaction-Spring, blue (strong)	17AD305=14	17AD305=20	17AD305=20
	1	Reaction-Spring, yellow (extra strong)	17AD306=14	17AD306=20	17AD306=20
	1	Ultra-Spring, grey (ultra-strong)	-	17AD307=20	17AD307=20

2.3.1 Optional accessories

- Fit-Kit 17AT3=1 (see fig. 11)
- Dummy kit 17AT4=1 (see fig. 11)

3 Intended use

3.1 Indications for use

The orthotic joint is intended **exclusively** for orthotic treatments of the lower limbs on a user.

For unilateral use, the manufacturer recommends using the lamination or prepreg technique to process the orthotic joint.

3.2 Indications

- Partial or total paralysis of the leg muscles
- Orthopaedic diseases of the lower limbs

Indications must be determined by the physician.

3.3 Combination possibilities

Subject to the permissible weight limits, the following orthotic knee joints can be combined with the 17AD300=*:

- 17KF110=*
- 17KF210=*
- 17KL100=*
- 17LK3=*

Observe the information in the instructions for use for all products when combining the orthotic joints.

Combination with 17KO1=* C-Brace joint unit

Orthotic ankle joint	Max. body weight [kg] (lbs)	Orthotic knee joint	Area of application
17AD300=16-T	85 (187.4)	17KO1=*	bilateral
17AD300=20-T	110 (242.5)		bilateral

Further information about the processing of the 17KO1=* C-Brace joint unit is found in the 647G1337 instructions for use.

Combination with 17B203=* E-MAG Active

Orthotic ankle joint	Max. body weight [kg] (lbs)	Orthotic knee joint	Area of application
17AD300=16*, =20*	85 (187.4)	17B203=*	unilateral
17AD300=16*, =20*	100 (220.5)	17B203=* with 17B206 medial support	bilateral

INFORMATION

Do not make any further static or dynamic changes to the 17AD300=* after automatic calibration of the 17B203=* E-MAG Active.

If user-specific adjustments to the 17AD300=* are required, repeat the process steps for automatic calibration of the 17B203=* E-MAG Active according to the 647G1165 instructions for use.

Further information about the processing of the 17B203=* E-MAG Active is found in the 647G1165 instructions for use.

3.4 Lifetime

The product is designed for a lifetime of **3 years** when used as intended and assembled professionally.

3.5 Qualification

Users may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

4 Safety

4.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

4.2 Safety instructions

CAUTION

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.

- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

⚠ CAUTION

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

NOTICE

Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

NOTICE

Thermal overloading of product components

Loss of function due to improper thermal treatment

- ▶ Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below **-10 °C (14 °F)** or above **+60 °C (140 °F)** (e.g. sauna, excessive sunlight, drying on a radiator).

5 Preparing the product for use

⚠ CAUTION

Incorrect alignment, assembly or adjustment

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

5.1 Fabricating the orthosis

5.1.1 Using the lamination dummy

Required materials: 17AD313=* lamination dummy set, 30Y594=* lamination dummy upper section, 17LS3=* lamination bar, 17SF300=N-* foot stirrup, 501S137=* countersunk screws, 743Y56=4 alignment insert, 636K8 Plastaband

> **Prerequisite: The orthosis model has been fabricated specifically for the user.**

- 1) Mount the rail (see fig. 2, item 9) with the screws (item 1) on the lamination dummy lower section (see fig. 3, item 1).
- 2) Place the foot stirrup (item 8) without the bearing bushing in the upper section of the lamination dummy (item 7).
- 3) Place the upper section on the lower section and secure it with the slotted oval head screw (see fig. 3, item 2). (Depending on the design, the foot stirrup with the upper section is brought to the neutral position.)
- 4) Attach the lamination dummy to the orthosis model and screw it to the alignment insert.

5.1.2 Optional use of the shaping dummy to determine the space requirement

When using the Reaction- and/or Ultra-module, Ottobock recommends using the shaping dummy to determine the subsequent space requirement of the modules and for shaping on the orthosis model.

- 1) Use the screws (item 2) to mount the upper sections of the shaping dummy (see fig. 4, item 5 and/or item 6) of the corresponding module to the shaping dummy lower section (item 1).
- 2) Loosen and remove the upper screw of the bar on the joint bar connection of the lower section of the lamination dummy (see fig. 2, item 1).
- 3) Mount the shaping dummy to the upper hole of the joint bar connection of the lamination dummy with the screw (see fig. 4, item 3) (see fig. 12).
- 4) Determine the subsequent course of the modules on the orthosis model and mark them with a pen.
- 5) Unscrew the shaping dummy.

5.1.3 Laminating the orthosis

- 1) Reattach the upper screw to the bar.
- 2) Seal the screws of the bar and lamination dummy with Plastaband.

- 3) Laminate the orthosis according to the patient's needs.

5.1.4 Optional use of the shaping dummy for shaping

- 1) To prevent the modules from subsequently colliding with the cured laminate, position the shaping dummy on the bar joint bar connection and press it onto the **uncured laminate** (see fig. 13).
- 2) Secure the shaping dummy to the orthosis model with adhesive tape (see fig. 14).

5.1.5 Finishing the orthosis

- 1) After the laminate has cured, remove the laminated blank from the orthosis model. To do so, cut the lamination dummy free with a knife and remove it.
- 2) Remove the screws from the Plastaband.
- 3) Remove the two screws from the bar.
- 4) Loosen the slotted oval head screw and remove the lamination dummy.
- 5) Use the oscillating saw to expose the edges for removal.
- 6) Remove the foot component and lower leg shell from the model.
- 7) Sand off sharp edges.

5.2 Mounting the orthotic joint

Axial washers

Make sure to use the same axial washer material thickness when selecting the axial washer, e.g. 0.90 and 0.90 mm. Optionally, you can use axial washers of varying thicknesses. Use the next larger or smaller thickness of the axial washer, e.g. 0.90 and 0.95 or 0.90 and 0.85. A combination of, for example, 0.85 and 1.10 mm is **not permitted** and may cause damage to the orthotic joint.

Size of the orthotic joint	Thickness of axial washers [mm]
17AD300=14*, 17AD300=16*, 17AD300=20*	0.85; 0.90; 0.95; 1.0; 1.05; 1.10

- 1) Select the axial washers and set them on the foot stirrup (thickness of the axial washers – see table). Use one axial washer per side for this purpose (see fig. 15).
- 2) Optimise the seat by replacing the axial washers step-by-step: Thicker washers if there is play in the medial-lateral direction, thinner washers if there is resistance in plantar flexion/dorsal extension.

- 3) Set the base body of the orthotic joint onto the foot stirrup with the axial washers (see fig. 16).
- 4) Use the insertion aid for the joint nut from the lamination dummy set 17AD313=* , to facilitate installation of the joint nut (see fig. 17).
- 5) ① Screw the insertion aid onto the joint nut.
- 6) ② Slide the joint nut through the opening on the orthotic joint to the stop.
- 7) ③ Unscrew the insertion aid from the joint nut.
- 8) ④ Remove the insertion aid.
- 9) Insert the joint screw into the orthotic joint (see fig. 18).
- 10) Tighten the joint screw with a torque wrench (see fig. 19).
- 11) Check whether the foot component can be moved smoothly and without play.
- 12) Before handing the orthosis over to the user, clean the screw connection of the joint screw using a degreasing cleaner, secure it with Loctite 241 and tighten it to the corresponding tightening torque (see following table).

17AD300=*		14	16	20
	Joint screw	T25 5 Nm	T30 8 Nm	T30 8 Nm

5.3 Installing the modules

INFORMATION

Bilateral design of the orthosis

Ensure that the forces (e.g. compression spring force) are evenly distributed for the patient-specific setting.

INFORMATION

Ottobock recommends checking the alignment of the orthosis using the 3D L.A.S.A.R. Posture and adapting it according to the patient's needs.

- We recommend positioning the sagittal load line approx. **15 mm** in front of the knee compromise pivot point according to Nietert (see fig. 20).
- Thanks to the individual possibilities for combining the specific modules with one another, the orthotic ankle joint can be adjusted to the needs of the user. The adjustment of the alignment can be controlled via the angle adjustment of the orthotic ankle joint. The angle adjustment can be set via the screw depth of the 17AD301=* , 17AD302=* set screws or the 17AD303=* , 17AD304=* stop pins.

5.3.1 Installing the Stop-module

Required materials: hexagon socket spanner, torque wrench, Torx key, Stop-module 17AD301=*

Prerequisite: The orthotic ankle joint is mounted to the foot stirrup.

- 1) Insert the screw-in sleeve through the cover from above (see fig. 21).
- 2) Screw the cover with the screw-in sleeve into the base body by hand. To do this, set the cover with the flattened side facing the insertion area onto the base body of the orthotic ankle joint.
- 3) Tighten the screw-in sleeve using a torque wrench and observe the specified tightening torque (see table).
- 4) Insert the stop bolt into the screw-in sleeve with the curve pointing towards the foot stirrup.
- 5) Insert the spring into the stop bolt and secure it in the screw-in sleeve.
- 6) Screw the stop pin into the screw-in sleeve using a Torx key and adjust the stop for the user.

17AD301=*		14	16	20
	Set screw	T25	T30	T30
	Screw-in sleeve	AF12 6 Nm	AF13 8 Nm	AF13 8 Nm

5.3.2 Installing the Spring-module

Required materials: hexagon socket spanner, torque wrench, Torx key, degreasing cleaner, 636K13 Loctite 241, 17AD302=* Spring-module

Prerequisite: The orthotic ankle joint is mounted to the foot stirrup.

- 1) Insert the screw-in sleeve through the cover from above (see fig. 22).
- 2) Screw the cover with the screw-in sleeve into the base body by hand. To do this, set the cover with the flattened side facing the insertion area onto the base body of the orthotic ankle joint.
- 3) Tighten the screw-in sleeve using a torque wrench and observe the specified tightening torque (see table).
- 4) Insert the ball into the screw-in sleeve.
- 5) **INFORMATION: If a lower dorsiflexion force is required, the thinner compression spring (see fig. 22) can be removed.** Insert the compression spring(s) into the spring guide housing and secure it in the screw-in sleeve.

- 6) Screw the set screw into the screw-in sleeve with a Torx key and adjust the spring force to the user.
- 7) Before handing the orthosis over to the user, clean the set screw using a degreasing cleaner and then secure it with Loctite 241; while doing so, observe the angle adjustment of the ankle joint.

17AD302=*		14	16	20
	Set screw	T25	T30	T30
	Screw-in sleeve	AF12 6 Nm	AF13 8 Nm	AF13 8 Nm

5.3.3 Installing the Reaction-/Ultra-module

The Reaction-/Ultra-module is delivered installed in the cover.

- ▶ Disassemble the Reaction-/Ultra-module into individual components.

Required materials: hexagon socket spanner, hexagon key, Torx key, torque wrench, 17AD303=*/Ultra-Modul 17AD304=* Reaction-module, 17AD308=* spring dummy set

Prerequisite: The orthotic ankle joint is mounted to the foot stirrup.

- 1) Insert the screw-in sleeve with O-ring through the cover from above and then screw it into the base body by hand. Set the cover with the flattened side facing the insertion zone onto the base body of the orthotic ankle joint (see fig. 24).
- 2) Tighten the screw-in sleeve using a torque wrench and observe the specified tightening torque (see the following table).
- 3) Screw the guide screw into one of the 3 bores of the screw-in sleeve so that it protrudes approx. **1 mm** out of the sleeve (see fig. 25). Our recommendation is to use the bore facing away from the body.
- 4) Position the stop bolt in the screw-in sleeve with the curve facing down (see fig. 26).
- 5) Slide the spring dummy onto the spring holder, noting the slots on the lower edge of the spring dummy (see fig. 27).
- 6) Insert the mounting adapter into the spring dummy to the stop until the magnetic contact is felt.
- 7) Hold the components horizontally and slide the spring onto the stop pin of the spring holder from below.

- 8) Hold the screw-in sleeve horizontally so the components do not fall out of the sleeve. Then position the mounting adapter with the spring dummy, spring holder and spring in the screw-in sleeve in the stop bolt, aligning the groove of the spring dummy with the guide screw on the screw-in sleeve (see fig. 28).
- 9) Hold the spring dummy and remove the mounting adapter (see fig. 29).
- 10) Put on the spring bolt and tighten it with a Torx key, turning in the spring bolt until the spring dummy has no more play in the module (see fig. 30).
- 11) Continue tightening the guide screw in until you feel resistance.

17AD303=* / 17AD304=*		14	16	20
	Spring bolt	T30	T50	T50
	Counter screw	AF2.5 2 Nm	AF3 3 Nm	AF3 3 Nm
	Stop pin	T10	T15	T15
	Screw-in sleeve	AF20 6 Nm	AF24 8 Nm	AF24 8 Nm

5.4 Assembling the lower leg shell

- ▶ Attach the lower leg shell to the base body of the orthotic joint using the countersunk screws, clean the screw connections using a degreasing cleaner and secure them with Loctite 241 and the corresponding tightening torque.

17AD300=*		14	16	20
	Countersunk screws	T20 3.5 Nm	T25 7 Nm	T25 7 Nm

5.5 Reaction-/Ultra-module – static fitting

Required materials: 743A8 knee pivot point gauge, Torx key

Prerequisite: The Reaction-/Ultra-module is mounted to the spring dummy in the orthotic ankle joint.

- 1) Put the orthosis on the patient.
- 2) Use the knee pivot point gauge to mark the compromise pivot point according to Nietert on the knee joint (see fig. 31).

- 3) Relieve the lower leg by hand so that the stop pin moves easily in the thread (see fig. 32).
- 4) Guide the Torx key through the spring bolt and adjust the stop to change the angle of the upper ankle joint (see fig. 33).
- 5) If 2 Reaction-/Ultra-modules are used, readjust the stop so the orthosis does not have any play in the anterior/posterior direction.
- 6) After the static fitting, use an Allen wrench to screw the counter screw through the hole of the spring bolt into the spring holder to secure the stop pin with the counter screw (see fig. 34).
- 7) Check the static alignment of the orthosis using the 3D L.A.S.A.R. Posture (see fig. 35).

5.6 SUB Reaction-/Ultra-module – dynamic fitting

- > **Required materials:** Reaction-module: 17AD305=* spring (strong, blue) or 17AD306=* spring (extra strong, yellow)
Ultra-module: 17AD307=* spring (ultra strong, grey)
 - > **Prerequisite:** The Reaction-/Ultra-module is mounted to the spring dummy in the orthotic ankle joint.
- 1) Use a Torx key to loosen the spring bolt and remove it from the module (see fig. 36).
 - 2) Remove the spring dummy and replace it with the Reaction-/Ultra-spring (see fig. 37), and remove the guide screw from the screw-in sleeve.
 - 3) Then screw in the spring bolt flush with the screw-in sleeve.
 - 4) Put the orthosis on the user and perform a fitting.
 - 5) Check the sagittal alignment of the user on the 3D L.A.S.A.R. Posture.
 - 6) Set the desired pretension of the spring.
 - 7) Have the user walk wearing the orthosis.
 - 8) After the final adjustment, attach the cover cap to the Reaction-/Ultra-module (see fig. 7, item 1).

5.7 Adjusting the Reaction-/Ultra-modules in a bilateral orthosis

When using multiple Reaction-/Ultra-modules in a bilateral orthosis, note the adjustment order. The stops and spring tensions of the Reaction-/Ultra-modules that face each other medially and laterally must be equally adjusted. The identical adjustment of Reaction-/Ultra-modules 1 and 3 as well as 2 and 4 prevents deformation and damage to the orthosis.

Adjusting dorsiflexion and plantar flexion

The stops of the Reaction-/Ultra-modules are adjusted according to the following adjustment sequence. If a Reaction-/Ultra-module is not required, the corresponding instructions can be ignored.

- > **Prerequisite:** The orthosis has been applied to the user.
- > **Prerequisite:** All Reaction-/Ultra-modules are mounted with the spring dummy.
 - 1) Adjust the stop of Reaction-/Ultra-module 1 (see fig. 38).
 - 2) Adjust the stop of Reaction-/Ultra-module 2 until the orthotic joint has no more anterior-posterior play.
 - 3) If Reaction-/Ultra-modules 1 and 2 have been adjusted and Reaction-/Ultra-module 3 is to follow: Loosen the spring bolt of Reaction-/Ultra-module 1.
 - 4) Adjust the stop of Reaction-/Ultra-module 3 (if Reaction-/Ultra-module 1 is present, then identically to Reaction-/Ultra-module 1).
 - 5) Now retighten the previously loosened spring bolt of Reaction-/Ultra-module 1.
 - 6) If Reaction-/Ultra-module 4 is to be adjusted: Loosen the spring bolt of Reaction-/Ultra-module 2.
 - 7) Adjust the stop of Reaction-/Ultra-module 4 until the orthotic joint no longer exhibits any anterior-posterior play (if Reaction-/Ultra-module 2 is present, then identically to Reaction-/Ultra-module 2).
 - 8) Now retighten the previously loosened spring bolt of Reaction-/Ultra-module 2.
 - 9) If a screw-in hole is left unused, secure the hole with a dummy cap.
 - 10) Optionally adjust the settings in the static fitting and repeat the steps.

Adjusting the spring tension

The spring tensions of the Reaction-/Ultra-modules are set according to the following adjustment sequence. If a Reaction-/Ultra-module is not required, the corresponding instructions can be ignored.

- > **Prerequisite:** The orthosis has been applied to the user.
- > **Prerequisite:** All Reaction-modules are mounted with the Reaction-Springs or the Ultra-modules are mounted with the Ultra-Spring.
 - 1) Adjust the spring tension of the spring bolt of Reaction-/Ultra-module 1 (see fig. 39).
 - 2) Adjust the spring tension of the spring bolt of Reaction-/Ultra-module 2.
 - 3) Using callipers, measure and note down the insertion depth of the spring bolt in the screw-in sleeve of Reaction-/Ultra-modules 1 and 2 (see previous figure).

- 4) Adjust the spring tension of the spring bolt of Reaction-/Ultra-module 3. If Reaction-/Ultra-module 1 is mounted, transfer the measurement value from Reaction-/Ultra-module 1 to Reaction-/Ultra-module 3.
- 5) Adjust the spring tension of the spring bolt of Reaction-/Ultra-module 4. If Reaction-/Ultra-module 2 is mounted, transfer the measurement value from Reaction-/Ultra-module 2 to Reaction-/Ultra-module 4.
- 6) If a screw-in sleeve is unused, secure it with a dummy cap.
- 7) Have the user walk wearing the orthosis. Readjust the springs if necessary, observing the adjustment sequence.

6 Cleaning

Clean the orthosis joint promptly after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Rinse the orthosis joint with clean, fresh water.
- 2) Dry the orthosis joint with a cloth or allow it to air-dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

7 Maintenance

INFORMATION

The product may be exposed to increased loads on a user-specific basis.

- Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

Corresponding service sets can be ordered for servicing the product.

The spare parts / wear and tear parts are listed in the section "Maintenance".

The product is designed to be free of lubricants.

Do not use lubricants.

7.1 Maintaining the orthotic joint

- 1) Disassemble the orthotic joint.
- 2) In case of wear and tear on individual components or on the orthotic joint, replace them with new components.

Spare parts / wear and tear parts

(see fig. 2)			Size 14	Size 16	Size 20
Item	Quantity	Designation	Article number		
6	1	Axial washers set	17AD315=14	17AD315=16	17AD315=20
7	1	Lamination dummy – upper section	17AD314=14	17AD314=16	17AD314=20

7.2 Maintaining the Stop-module

- 1) Remove the module and disassemble it into individual components.
- 2) If there is wear and tear on the Stop-module, use the corresponding 17AD318=* service set or replace the module with a new one (see page 48, see fig. 21).

(see fig. 5)			Size 14	Size 16	Size 20
Item	Quantity	Designation	Article number		
not illustrated	1	Service set for Stop-module	17AD318=14	17AD318=20	17AD318=20
		consisting of:			
		Spring			
		Sleeve for stop pin			
		Stop bolt			

7.3 Maintaining the Spring-module

- 1) Remove the module and disassemble it into individual components.
- 2) In case of wear and tear, replace the module with a new one.

7.4 Maintaining the Reaction-/Ultra-module

- ▶ **Replace the individual components of the Reaction/Ultra-module every six months, and use the corresponding 17AD317=* service set for this purpose.**

Required materials: hexagon socket spanner, hexagon key, Torx key

- 1) Disassemble the Reaction-/Ultra-module into individual components (see fig. 7).
- 2) Remove the cover cap from the spring bolt.
- 3) Use a Torx key to loosen and remove the spring bolt.
- 4) Remove the Reaction- or Ultra-Spring from the screw-in sleeve.
- 5) Remove the spring holder with the spring and stop bolt.

Checking

- 1) Inspect the stop bolt and stop pin of the spring holder and replace them in case of wear and tear.
- 2) Replace the counter screw when the static alignment is readjusted or adjusted. In all other cases, it is generally not necessary to replace the counter screw.
- 3) Carefully check and clean all components prior to assembling them. Replace the components depending on the degree of wear.

For assembly of the Reaction-/Ultra-module, see the section "Installing the Reaction-/Ultra-module".

(see fig. 7)			Size 14	Size 16	Size 20
Item	Quantity	Designation	Article number		
not illustrated	1	Service set for Reaction-/Ultra-module	17AD317=14	17AD317=16	17AD317=20
		consisting of:			
		Spring			
		Sleeve for stop pin			
		Stop bolt			
		Counter screw			

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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