







50K14 Genu Arexa PCL

DE Gebrauchsanweisung	4
EN Instructions for use	9
FR Instructions d'utilisation	14
IT Istruzioni per l'uso	19
ES Instrucciones de uso	24
PT Instruções de utilização	29
NL Gebruiksaanwijzing	34
SV Bruksanvisning	39
DA Brugsanvisning	44
NO Bruksanvisning	49
FI Käyttöohje	54
PL Instrukcja użytkowania	59
CS Návod k použití	64
SK Návod na použitie	69
RU Руководство по применению	74
ZH 使用说明书	79

	Seite/Side	Größe/Size	Umfang/Circumference			
			A (cm)	A (inch)	B (cm)	B (inch)
	L/R	S	32-37	12.6-14.6	39-46	15.4-18.1
	L/R	M	37-41	14.6-16.1	46-53	18.1-20.9
	L/R	L	41-46	16.1-18.1	53-61	20.9-24.0
	L/R	XL	46-50	18.1-19.7	61-69	24.0-27.2
	L/R	XXL	50-56	19.7-22.1	69-79	27.2-31.1

Material	AL, SS, klettbarer Velours aus PA/velastic velour of PA, PU-Schaum/PU-foam, EVA-Schaum/EVA-foam, PA, POM, TPU, PP, PK
-----------------	---

INFORMATION

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INFORMATION

Date of last update: 2025-01-06

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50K14 Genu Arexa PCL knee joint orthosis.

2 Product description

Scope of delivery (see fig. 1)

Item	Component	Item	Component
1	Pad button	6	Pad retainer
2	Strap guide loop	7	Extension stop (angular hand grip)
3	Tibia pad	8	Flexion stop (round hand grip)
4	Clip rivet	-	Closures 1 to 6
5	Condyle pad		

Scope of delivery (see fig. 2)

Item	Component	Item	Component
1	Calf pad	5	Fastening strap
2	Dial	6	Quick-release closure
3	Closure holder	7	Posterior lower leg pad
4	Screw connection	8	Additional posterior lower leg pad (optional)

3 Intended use

3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

3.2 Indications

- Strain to the posterior cruciate ligament
- Partial rupture and rupture of the posterior cruciate ligament (acute treatment, conservative therapy, postoperative treatment)

Indications must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

Not known.

3.3.2 Relative Contraindications

- Sensory disorders of the lower limb (e.g. polyneuropathy)
- Circulatory disorders
- Skin changes in the treatment area (e.g. skin diseases, skin injuries, scars)
- Inflammation, reddening and hyperthermia in the treated body region
- Circulatory disorders, lymphatic and venous insufficiencies, severe varicose veins in the lower limbs

3.4 Effects

The orthosis stabilises and guides the knee joint.

3.5 Lifetime

The product is designed for a maximum lifetime of **2 years**.

3.6 Qualification

Users may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

4 Safety

4.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

4.2 General safety instructions

Use of the product on damaged skin

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product on healthy skin.

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

Improper moulding or application

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

⚠ CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

NOTICE

Incorrect setting

Damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ The orthosis may only be adjusted by qualified personnel.
- ▶ Do not make any improper changes to the settings.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

5 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

5.1 Size Selection

- 1) Measure the circumference of the knee joint (circumference A).
- 2) Measure the circumference of the thigh (circumference B) 15 cm above the centre of the patella.
- 3) Determine the size of the orthosis (see page 3, Fig. 7).

5.2 Adaptation

Adjusting the range of motion

INFORMATION

The joint stops (medial and lateral) must have the same degree settings.

Available extension stops (square grip):	0°, 10°, 20°, 30°, 45°
Available flexion stops (round grip):	0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

- 1) Remove the pad retainer (and the spacer rings if applicable) (see fig. 4).
- 2) Remove the pinch protection stop (see fig. 1, item 8).
- 3) Insert the desired extension and flexion stops into the joint as specified by the doctor (see fig. 5).
→ The stops engage with slight resistance.
- 4) **Optional:** Snap the required number of spacer rings onto the joint (see fig. 6).
- 5) Attach the pad retainer.

Adapting the orthosis

- > **Prerequisite:** The patient is sitting on the front edge of a chair and bends the knee at **70°**. All straps are opened.
- 1) Set the maximum width of the top lower leg closure on dial 5 (see fig. 2, item 2).
 - 2) Place the orthosis on the leg.
 - 3) Align the centre of the orthotic joints according to the compromise pivot point according to Nietert (see fig. 3).

Adjusting the strap length

- > The orthosis has been taken off.
- 1) Remove the pad.
 - 2) Adjust the strap length on the hook-and-loop strap.
 - 3) Attach the pad.

Adjusting the knee width

The knee width of the orthosis can be adjusted with the included spacer rings and the condyle pads.

- 1) Remove the condyle pads and the pad retainers (see fig. 4).
- 2) Attach the required number of spacer rings to the orthotic joint (see fig. 6).
- 3) Attach the pad retainers.
- 4) Attach the condyle pads to the pad retainers.

Adjusting the top lower leg closure

- 1) Fasten closure 1 of the top lower leg closure (see fig. 2, item 6).
INFORMATION: The position of the calf pad is specified by the mounting hole in the orthosis frame.
- 2) **Optional:** Move the top lower leg closure down by one hole.
- 3) Pre-tension the calf pad using the dial.
NOTICE! To avoid a risk of breakage, do not over-tighten the dial.
- 4) Check the pre-tension with the patient standing up and correct as needed.

Repositioning the top lower leg closure

- 1) Remove the screw connections of the closure holder and the closure on both sides.
- 2) Reposition the closure holder and then insert the screws and tighten them.

Adjusting the posterior lower leg pad

- > **Prerequisite:** The orthosis has been applied to the patient.
- 1) Check the length of the posterior lower leg pad medially and laterally (see fig. 2, item 7).
 - 2) **Optional:** Shorten the posterior lower leg pad (item 7) equally on both sides at the pre-stamped cut edges using scissors. To do so, take the orthosis off of the patient to avoid cutting them.
 - 3) **Optional:** If the adjustment range of the calf pad is not sufficient to generate sufficient contact pressure, use the additionally supplied posterior lower leg pad (item 8). To do so, attach the additional posterior lower leg pad (item 8) to the posterior lower leg pad (item 7).
 - 4) Put the orthosis on the patient and check whether the calf padding fits. Repeat the steps until the desired result is obtained.

5.3 Application

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

> **Prerequisite:** All straps have been opened, and the orthosis has been applied.

- 1) Fasten the top lower leg closure 1 (see fig. 3).
- 2) Fasten closures 2, 3 and 4 one after the other.
- 3) Turn dial 5 until the calf pad is sitting firmly and comfortably against the calf.
- 4) Fasten closure 6.
- 5) Tighten the closures so the orthosis is sitting firmly and comfortably.
- 6) **Optional:**After **15 min of** active movement, retighten the closures in the specified sequence.

5.4 Removal

- 1) Loosen the pre-tension of the calf pad with dial 5 (see fig. 3).
- 2) Remove quick-release closure 1 from the pad button.
- 3) Open closures 2, 3 and 4.
- 4) Remove the orthosis from the leg.

6 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the orthosis regularly:

Cleaning the straps and pads

- 1) Remove the straps and pads from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Hand-wash the straps and pads in warm water at **30 °C (86 °F)** using neutral detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 5) Attach the straps and pads to the orthosis.

Cleaning the orthosis

If the orthosis comes into contact with salt water or dirt:

- ▶ Rinse the orthosis with clear water and allow to air dry.

7 Disposal

Dispose of the product in accordance with national regulations.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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