LP2-W2 Freestyle Swim

DE	Gebrauchsanweisung	3
EN	Instructions for use	6
FR	Instructions d'utilisation.....	9
NL	Gebruiksaanwijzing	12

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Freestyle Swim Foot Instructions for Use

Product Number: LP2-W2


Assembly

The Freestyle Swim Foot module is pre-assembled consisting of a graphite component (keel), Spectra™ sock, and foot shell. Foam *stiffening bumpers* for increasing heel stiffness are provided. After dynamic alignment, torque pyramid adjustment screws to the manufacturer's specifications. Secure pyramid adjustment screws with a thread locking adhesive (i.e., Loctite 242).

Bench Alignment

Prior to donning the prosthesis:

- Plantarflex/Dorsiflexion foot to match the shoe heel height.
- Adduct/Abduct socket to provide appropriate frontal plane angle.
- Flex/Extend socket to provide appropriate sagittal plane angle.
- Move the socket linearly to ensure the weight line falls along the anterior edge of the pylon (see illustration).

Dynamic Alignment

During *loading response*, the heel lever stores energy and releases it during *midstance*. This action provides momentum for the keel to store energy and release it during *terminal stance*. To optimize the heel to toe rollover motion, adjust the following variables:

- Anterior/posterior foot placement
- Dorsiflexion/plantarflexion
- Heel stiffness

Troubleshooting

Heel too soft

Symptoms

- Foot flat occurs too rapidly
- Toe feels excessively stiff
- Knee hyperextension

Solutions

- Shift socket anteriorly in relation to the foot
- Attach foam *stiffening bumpers*

Heel too hard

Symptoms

- Rapid knee flexion, instability
- Heel to toe progression to rapid
- Lack of energy return sensation

Solutions

- Shift socket posteriorly in relation to the foot
- Verify appropriate foot module category

	Foot module too stiff	<p><i>Symptoms</i></p> <ul style="list-style-type: none"> • Flat spot in rollover motion at slow cadences <p><i>Solutions</i></p> <ul style="list-style-type: none"> • Consider a lower category foot module
	Foot module too soft	<p><i>Symptoms</i></p> <ul style="list-style-type: none"> • Clicking noise at <i>initial contact</i> • Excessive toe deflection with high impact activity <p><i>Solutions</i></p> <ul style="list-style-type: none"> • Consider a higher category foot module
Stiffening Bumpers	Foam <i>stiffening bumpers</i> are included to adjust the heel stiffness during <i>loading response</i> . The bumpers may be temporarily attached with tape between the heel lever and keel. If the bumpers provide excessive heel stiffness, exclude or trim as appropriate. For permanent placement, adhere bumpers using Barge or similar contact cement.	
Spectra™ Sock	A Spectra™ sock is provided to protect the foot shell and minimize noise. Spectra™ socks must be replaced at intervals appropriate to the user's activity level. Failure to inspect and replace the Spectra™ socks may prematurely wear the foot module, and will void the warranty.	
Foot Shell	When removing or installing the foot shell, use the Foot Shell Removal Tool (ACC-00-10300-00) to prevent damage to the foot module.	
Freestyle Swim Foot	Minimum clearance: 102mm-117mm Maximum user weight: 116 kg (255 lbs) Available sizes: 22cm-31cm Heel height: 10mm	
Operation	The Freestyle Swim Foot articulates from a neutral position (0°) to plantarflexed position (70°) for aquatic use. The Freestyle Swim Foot may be used in both fresh and salt water environments. To articulate, follow these steps: <ul style="list-style-type: none"> • Depress spring-loaded button and plantarflex foot module. • The spring-loaded button will automatically return to its original position once the foot module is fully plantarflexed to 70°. 	
Maintenance	The foot module requires periodic maintenance. <ul style="list-style-type: none"> • Inspect the foot module every six months. If the user is more active, more frequent inspection may necessary. Service as necessary. Replace Spectra™ sock and/or foot shell if worn to prevent damage to the graphite components. • The foot module may be cleaned and/or disinfected with soap and warm water. If the foot is exposed to water, remove the foot shell to clean and disinfect its interior. 	

Warnings

Failure to adhere to the guidelines of the *Instructions for Use* will void the warranty.

- Never use the foot module without a foot shell. Failure to comply may cause premature wear, loss of function, and/or product failure.
- Always use the foot module with a sock and shoe for non-aquatic activities. Failure to comply may cause premature wear, loss of function, and/or product failure.
- Never allow aggregates such as sand to remain in the foot shell. Upon exposure to aggregates, immediately disassemble foot module and rinse with water. The abrasive properties of aggregates will quickly wear the graphite components of the foot module.
- The Freestyle Swim Foot is intended to only be used in the plantarflexed position for aquatic sports.
- Do not ambulate with the Freestyle Swim Foot in the plantarflexed position.
- Do not use without the spring-loaded button in its original position to prevent bodily injury.

- Freedom Innovations foot modules are manufactured to fit industry standard pyramids and receivers. It is the prosthetist's responsibility to select and/or fabricate properly fitting attachment components.
- Never attempt to loosen the bolt affixing the pyramid connector.
- Discontinue use and consult your prosthetist if any part of the prosthesis starts to make noise.
- Inform your prosthetist if you lose or gain a significant amount of weight.
- Freedom Innovations foot products are manufactured and tested for a particular weight and activity impact level. Use by another user for whom it was not originally manufactured may cause injury and shall void any written or implied warranty.

Instructions for use Heel Bumper

To temporarily place foam bumpers on the LP Symes or Freestyle Swim Foot, use two-sided tape (included), and secure foam to the underside of the keel. For permanent placement, we recommend the use of Barge Cement.

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Advena Limited. Tower Business Centre, 2nd Flr., Tower